
## 1S49, 1S66, 1S67, 1S90 SACH

<b>DE</b>	Gebrauchsanweisung .....	3
<b>EN</b>	Instructions for use .....	11
<b>FR</b>	Instructions d'utilisation .....	19
<b>IT</b>	Istruzioni per l'uso .....	28
<b>ES</b>	Instrucciones de uso .....	37
<b>PT</b>	Manual de utilização .....	45
<b>NL</b>	Gebruiksaanwijzing .....	54
<b>SV</b>	Bruksanvisning .....	62
<b>DA</b>	Brugsanvisning .....	70
<b>NO</b>	Bruksanvisning .....	79
<b>FI</b>	Käyttöohje .....	87
<b>PL</b>	Instrukcja użytkowania .....	95
<b>HU</b>	Használati utasítás .....	103
<b>CS</b>	Návod k použití .....	111
<b>RO</b>	Instrucțiuni de utilizare .....	119
<b>HR</b>	Upute za uporabu .....	128
<b>SL</b>	Navodila za uporabo .....	136
<b>SK</b>	Návod na používanie .....	144
<b>BG</b>	Инструкция за употреба .....	152
<b>TR</b>	Kullanma talimatı .....	160
<b>EL</b>	Οδηγίες χρήσης .....	168
<b>RU</b>	Руководство по применению .....	178
<b>JA</b>	取扱説明書 .....	187
<b>ZH</b>	使用说明书 .....	195
<b>KO</b>	사용 설명서 .....	201

## INFORMATION

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### INFORMATION

Last update: 2017-06-22

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

### 1.1 Construction and Function

The 1S49, 1S66, 1S67 and 1S90 SACH (Solid Ankle Cushion Heel) prosthetic feet are suitable for use in both modular prostheses and exoskeletal prostheses. The functional properties are achieved through the combination of a contoured core and functional foam.

## 1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

Permissible ankle blocks and screw connections	
Foot size [cm]	Reference number (ankle block and screw connection)
24 to 25	2K34=25 and 2Z22=M10
26 to 30	2K34=30 and 2Z22=M10

Permissible foot adapters			
Prosthetic foot	Size [cm]	Body weight [kg]	Permissible foot adapter
1S67	22 to 25	≤100	2R8=M8, 2R31=M8, 2R54=M8
	26		2R8=M10, 2R31=M10, 2R54=M10
	26	101 to 125	2R8=M10, 2R31=M10
1S49, 1S66, 1S90	21 to 30	≤100	2R8=M10, 2R31=M10, 2R54=M10
	26 to 30	101 to 125	2R8=M10, 2R31=M10

## 2 Intended use

### 2.1 Indications for use

The product is to be used solely for lower limb prosthetic fittings.

### 2.2 Area of application

Our components perform optimally when paired with appropriate components based upon weight and mobility grades identifiable by our MOBIS classification information and which have appropriate modular connectors.

The product is recommended for mobility grade 1 (indoor walker) and mobility grade 2 (restricted outdoor walker).

- The maximum approved body weight is specified in the technical data (see Page 19).

### 2.3 Environmental conditions

Allowable environmental conditions
Temperature range for use: -10 °C to +60°C
Allowable relative humidity 0 % to 90 %, non-condensing

<b>Unallowable environmental conditions</b>	
Mechanical vibrations or impacts	
Perspiration, urine, fresh water, salt water, acids	
Dust, sand, highly hygroscopic particles (e.g. talcum)	
<b>Other environmental conditions</b>	
<b>UV light</b>	Not resistant
<b>Hydrolysis</b>	Not resistant

## 2.4 Service life

This product has been tested by the manufacturer with 2 million load cycles. Depending on the patient's activity level, this corresponds to a service life of 2 to 3 years.

## 3 Safety

### 3.1 Explanation of warning symbols

 <b>CAUTION</b>	Warning regarding possible risks of accident or injury.
 <b>NOTICE</b>	Warning regarding possible technical damage.

### 3.2 General safety instructions

<p> <b>CAUTION</b></p> <p><b>Exceeding the service life and reuse on another patient</b></p> <p>Risk of injury due to loss of functionality as well as damage to the product</p> <ul style="list-style-type: none"> <li>▶ Ensure that the approved service life is not exceeded.</li> <li>▶ Only use the product on a single patient.</li> </ul>
<p> <b>CAUTION</b></p> <p><b>Excessive strain on the product</b></p> <p>Risk of injury due to breakage of load-bearing components</p> <ul style="list-style-type: none"> <li>▶ Use the product according to the specified area of application (see Page 12).</li> </ul>
<p> <b>CAUTION</b></p> <p><b>Unallowable combination of prosthetic components</b></p> <p>Risk of injury due to breakage or deformation of the product</p> <ul style="list-style-type: none"> <li>▶ Only combine the product with prosthetic components that are approved for that purpose.</li> </ul>

- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

#### NOTICE

#### Mechanical damage to the product

Change in or loss of functionality due to damage

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

#### NOTICE

#### Use under unallowable environmental conditions

Damage to product due to unallowable environmental conditions

- ▶ Do not expose the product to unallowable environmental conditions.
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop, etc.).

#### Signs of changes in or loss of functionality during use

Decreased forefoot resistance or changes in roll-over behaviour are noticeable indications of loss of functionality.

## 4 Scope of delivery

Quantity	Designation
1	Instructions for use
1	Prosthetic foot

#### Spare parts/accessories (not included in the scope of delivery)

Fig.	Item	Designation	Reference number
2	1	Foot adapter with screw connection	2R8*, 2R31*, 2R54*
2	2	Connection plate	2R14
3	-	Ankle block (exoskeletal design)	2K34=*

Spare parts/accessories (not included in the scope of delivery)			
Fig.	Item	Designation	Reference number
4	–	Screw connection for ankle block (exoskeletal design)	2Z22=M10
–	–	Ottobock special adhesive with hardener	636W28 or 636W18 and 636W19

## 5 Preparation for use

### CAUTION

#### Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- Observe the alignment and assembly instructions.

### INFORMATION

- Not all of the materials may be available in your country. In this case, please contact the local branch of the manufacturer to obtain information on alternative materials.

### 5.1 Installing the foot adapter

- Install the foot adapter on the prosthetic foot according to its instructions for use.

### 5.2 Assembling the ankle block (exoskeletal design)

- > **Required materials :** Screw connection (see Page 14), 710D4 torque wrench, 40 grit sandpaper, 617H21 sealing resin, 617P37 hardener powder, suitable solvent (e.g. 634A28), only for 1WR95: 636W17 plastic adhesive with 636W26 cross-linking agent
- 1) Insert the threaded socket of the screw connection in the proximal side of the hole in the ankle block.
  - 2) Insert the screw on the foot sole side in the prosthetic foot and screw the ankle block to the prosthetic foot (mounting torque: **30 Nm**).
  - 3) Conduct the fitting of the ankle block.
  - 4) Loosen the screw connection between the prosthetic foot and the ankle block.
  - 5) Sand the ankle block according to the contour of the prosthetic foot.
  - 6) Roughen the connection surface of the prosthetic foot in an a – p direction using the sandpaper.
  - 7) Remove the dust from the sanded surfaces.

- 8) Mix the sealing resin with **2%** to **3%** hardener powder for the connection between the ankle block and the prosthetic foot.
- 9) **INFORMATION: Do not put the resin in the hole for the screw connection.**  
Apply the resin evenly to the connection surfaces of the ankle block and prosthetic foot. At least **2 applications** are required due to the absorbency of the wood.
- 10) When the resin is slightly cured (threads on the spatula), align the ankle block and prosthetic foot with each other.
- 11) Insert the screw on the foot sole side in the prosthetic foot and screw the ankle block to the prosthetic foot (mounting torque: **30 Nm**).
- 12) Wipe off resin that comes out of the joint using a cloth.
- 13) **1WR95 only:** Apply plastic adhesive to the plug for the screw connection and press it fully into the prosthetic foot.
- 14) Remove any remnants of the adhesive using a suitable solvent.
- 15) Allow the adhesive to cure.

### 5.3 Alignment

#### NOTICE

#### Grinding the prosthetic foot

Premature wear resulting from damage to the prosthetic foot

- ▶ Do not grind the prosthetic foot.

#### 5.3.1 Bench Alignment

##### TT bench alignment

Bench alignment process	
<b>Required materials:</b> 662M4 goniometer, 743S12 heel height measuring apparatus, 743A80 50:50 gauge, alignment apparatus (e.g. 743L200 L.A.S.A.R. Assembly or 743A200 PROS.A. Assembly)	
Perform the assembly and alignment of the prosthetic components in the alignment apparatus according to the following specifications:	
<b>Sagittal plane</b>	
<b>①</b>	Heel height: <b>Effective heel height</b> (shoe heel height - sole thickness of fore-foot) + <b>5 mm</b>
<b>②</b>	Exterior foot rotation: <b>approx. 5°</b>
<b>③</b>	Anterior placement, middle of the prosthetic foot to the alignment reference line: <b>30 mm</b>
<b>④</b>	Connect the prosthetic foot and prosthetic socket using the chosen adapters. The instructions for use of the adapters must be observed.

<b>Bench alignment process</b>	
<b>5</b>	Determine the centre of the prosthetic socket with the 50:50 gauge. Align the prosthetic socket centrally to the alignment reference line. Socket flexion: <b>Individual residual limb flexion + 5°</b>
<b>Frontal plane</b>	
<b>6</b>	Alignment reference line of prosthetic foot: <b>between big toe and second toe</b> Alignment reference line of prosthetic socket: <b>along the lateral patella edge</b>
<b>7</b>	Observe the abduction or adduction position.

### **TF bench alignment**

- ▶ Observe the information in the prosthetic knee joint instructions for use.

### **5.3.2 Static Alignment**

- Ottobock recommends checking the alignment of the prosthesis using the L.A.S.A.R. Posture and adapting it as needed.
- If necessary, the alignment recommendations (TF modular leg prostheses: **646F219\***, TT modular leg prostheses: **646F336\***) may be requested from Ottobock.

### **5.3.3 Dynamic Trial Fitting**

- Adapt the setup of the prosthesis in the frontal plane and the sagittal plane (e.g. by making angle or slide adjustments) so that correct heel contact and optimal roll-over are ensured.
- **TT fittings:** Make sure that physiological knee movement is achieved after the heel makes impact and the leg begins to bear weight.

### **5.4 Optional: Installing the foam cover**

A connecting element (such as a connection plate, connection cap or foam connection cap) provides a detachable connection between the foam cover and the prosthetic foot.

- > **Required materials:** Degreasing cleaner (e.g. 634A58 isopropyl alcohol), 636N9 contact adhesive or 636W17 plastic adhesive with 636W26 cross-linking agent
- 1) Measure the distance from the knee rotation point to the desired end of the foam cover and add the compression allowance: For PE foam **10 mm**, for PUR soft foam **30 mm**. For TF prostheses, add twice the compression allowance.
  - 2) Cut the pre-shaped foam block to length.
  - 3) Pull the foam cover over the prosthesis.

- 4) Set the connecting element onto the footshell or prosthetic foot. Depending on the version, the connecting element engages in the edge or sits on the foot adapter.
- 5) Mount the prosthetic foot on the prosthesis.
- 6) Mark the outer contour of the connecting element on the distal face of the foam cover.
- 7) Dismount the prosthetic foot and remove the connecting element.
- 8) Clean the connecting element using a degreasing cleaner.
- 9) Glue the connecting element onto the distal face of the foam cover according to the marked outer contour.
- 10) Let the glue dry (approx. **10 minutes**).
- 11) Mount the prosthetic foot and adapt the exterior cosmetic cover shape. Take into account compression due to cosmetic stockings or SuperSkin.

## **6 Cleaning**

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

## **7 Maintenance**

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

## **8 Disposal**

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in the country of use regarding return, collection and disposal procedures.

## **9 Legal information**

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### **9.1 Liability**

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The

manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

## 9.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH