Instructions for Use

4-Hole Socket Adapter AL
4-Hole Socket Adapter TI
4-Hole Socket Adapter Junior

EN Instructions for Use	3	JA 取扱説明書	71
DE Gebrauchsanweisung	7	ZH 中文说明书	76
FR Notice d'utilisation	12	KO 사용 설명서	79
ES Instrucciones para el uso	16	CS Návod k použití	83
IT Istruzioni per l'uso	21	SK Indikácie použitia	88
NO Bruksanvisning	26	RO Indicații de utilizare	92
DA Brugsanvisning	30	HR Upute za uporabu	97
SV Bruksanvisning	34	HU Használati javallat	101
EL Οδηγίες Χρήσης	39	BG Указания за употреба	106
FI Käyttöohjeet	44	SL Indikácie pre použitie	110
NL Gebruiksaanwijzing	48	UA Інструкція з використання	115
PT Instruções de Utilização	53	AR تعليمات الاستخدام	120
PL Instrukcja użytkowania	57		
TR Kullanım Talimatları	62		
RU Инструкция по использованию	66		

Medical Device

INTENDED USE

The device is intended as a connection component of a prosthetic system that replaces a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

For more detailed device information, see **Table 1**.

Indications For Use and Target Patient Population

- Lower limb loss, amputation, or deficiency
- No known contraindications

For weight limit and impact level, see **Table 2**.

GENERAL SAFETY INSTRUCTIONS

Warning: Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

Warning: Only use Össur replacement parts for the device for ensured durability.

The device is for single patient use.

ASSEMBLY INSTRUCTIONS

Instructions for integration into laminated sockets:

Note: The socket must be laminated in two layers.

1. Isolate cast with PVA bag.
2. Apply appropriate layup material for first lamination.
3. Apply single 50 mm carbon tapes crossed over distal end of cast, in AP/ML direction.
4. Apply outer PVA bag and laminate.
5. Once cured, remove PVA bag and roughen outer surface of lamination.
6. Attach the Adapter in correct alignment and rotation with putty.
7. Position Lamination Cap on Adapter and secure it by locking the screw.
8. Tie 50 mm carbon tapes into grooves on Adapter, position in AP/ML direction.
9. Seal putty and Lamination Cap with sealing resin.
10. Apply layup material and combinations of material appropriate to the individual's weight and impact level.
11. Layers of material may be trimmed to different lengths to provide a thinner, lighter proximal socket.
12. Tie 50 mm carbon tapes into grooves on Adapter, position in AP/ML direction.
13. Continue with appropriate layup material.
14. Laminate with 80:20 acrylic lamination resin.

Warning: Tie 50mm carbon reinforcements into both grooves, between lamination cap and adapter, and groove in adapter.

Note: Remove the Lamination Cap after lamination by unscrewing from Adapter.

USAGE

Cleaning and Care

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

Environmental Conditions

See which environmental conditions apply to the device in **Table 1**.

A Weatherproof device can be used in a wet or humid environment and can tolerate being splashed by fresh water (e. g., rain), no submersion is allowed.

No contact with salt water or chlorinated water is allowed.

Dry with a cloth after contact with fresh water or humidity.

Clean with fresh water in case of accidental exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

A Waterproof device can be used in a wet or humid environment and submerged in up to 3-meter-deep water for a maximum of 1 hour.

It can tolerate contact with: Salt water, chlorinated water, perspiration, urine, and mild soaps.

It can also tolerate occasional exposure to sand, dust, and dirt. Continuous exposure is not allowed.

Dry with a cloth after contact with fresh water or humidity.

Clean with fresh water after exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

SPECIFICATIONS

Table 1 - Device Specifications

Part #	Description	Material	Environmental Conditions
A-122100	4-Hole Socket Adapter AL	Aluminum	Weatherproof
AJ-122100	4-Hole Socket Adapter Junior		Waterproof
A-125100	4-Hole Socket Adapter TI	Titanium	

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance

This device has been tested according to ISO 10328 standard to three million load cycles. Depending on patient activity, this may correspond to 3-5 years of use.

Table 2 - Usage Specifications

Part #	Impact Level	Weight Limit	ISO 10328 - Label Text
AJ-122100	The device is intended for children and adolescent patients.	60 kg	ISO 10328 - P4 - 60 kg
A-122100	The device is for low to moderate impact use, e.g., walking.	100 kg	ISO 10328 - P5 - 100 kg
A-125100	The device is for low to high impact use, e.g., walking and occasional running.	166 kg	ISO 10328 - P7 - 166 kg

ISO 10328 - "P" - "m"kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use see manufacturer's written instructions on intended use!

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