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AliMed® Therapy Carrot™ Hand Contracture Orthosis Kit

Instructions for Use and Care

MD

Intended Use

The AliMed Therapy Carrot Hand Contracture Orthosis Kit treats patients by increasing extension at the MP, PIP, and IP joints of the fingers where flexion contracture exists. The device helps protect from excessive moisture, pressure, infection, and risk of nail puncture.

Instructions for Use

For use under clinician guidance. See Sizing Chart below.

 Single-patient use only. Reusable.

1. Thread the loop at the end of the device through the eye of the placement wand.
2. Pull the loop over the tip of the placement wand and pull the wand to hold the loop taut (*Figure 1*).
3. Gently slide the placement wand and the attached Therapy Carrot through the patient's fist either at the index finger or little finger, determined by clinician based on contracture severity.
4. Pull through the fist until resistance is felt and the patient is settled on a comfortable diameter. If the loop is within reach at the other side of the fist, remove the placement wand as instructed in Step 5, and pull by the loop.
 - 4.1. May pull the device back through original entrance to decrease diameter.
5. Remove the placement wand by lifting the loop over the placement wand tip and pulling away.
6. Track progress by taking note on the number visible at the little finger or index finger. May mark progress on the device. An increase in number indicates a larger diameter.
7. (Optional) Write patient's name on device for identification.


Figure 1

Sizing Chart				
REF	Color	Length (Tip to Maximum Circumference)	Minimum (Documentable) Circumference	Maximum Circumference
#2970001888	Orange	8½"	2¼"	5"
#2970001889	Green	9"	2¼"	5¾"
#2970002637, #2970002644	Antimicrobial Blue	9"	2½"	5½"

Contraindications

There are no known contraindications for these devices.

Warnings

- Do not pull too hard when inserting the device into the fist.
- Use caution when removing the device from the fist to reduce pain, abrasion, and skin damage.
- Stop if the patient is experiencing pain. Consider using a smaller sized Therapy Carrot.
- Device can be safely used for several hours.
- Monitor the patient for skin irritation. If found, discontinue use of the device.
- Establish a schedule for use of the device for the patient.
- Regularly clean the device.
- Do not use the device if it is damaged (i.e., torn, punctured, inside is exposed).
- Consult your institution's guidelines for assessing the risk of each individual patient.
- Report any serious incident to AliMed and the appropriate regulatory body.

AliMed®, LLC

297 High St., Dedham, MA 02026 USA
PHONE 800.225.2610 • FAX 800.437.2966

Care and Cleaning

Clean with mild detergent and cold water. Air dry. Do not wring. Using improper cleaning solutions or methods can result in damage to the product or irritation to the patient.

Materials

Cotton/polyester or Aglon-coated polyester fleece casing, natural wool or LLDPE bead filling, polyester ribbon, and nylon placement wand.
Not made with natural rubber latex.

Storage/Disposal

Store in a clean, dry environment protected from direct sunlight. Storage at ambient room temperature is advised. Dispose of in accordance with existing federal, state, and local environmental control laws.

Symbol Glossary

 Medical Device	 Caution	 Batch Code	 Device is non-sterile	 Does not contain latex	 Manufacturer's catalog number
 Single-Patient, Multiple Use	 Manufacturer	 Date of Manufacture	 UK Conformity Assessment	 Conformité Européenne	 European Authorized Representative

Item Numbers

To Reorder:

	2970001888	2970001889	2970002637	2970002644	2970018882	2970002645	2970018886
(Former Item #)	(51024)	(51025)	(52274)	(52308)	(D51024)	(52310)	(D51323)

Replacement Items:

	2970001900	2970001915	2970002638	2970018883	2970001934
(Former Item #)	(51026)	(51027)	(52275)	(D51026)	(51028)

Additional Accessory:

	2970001962
(Former Item #)	(510381)

Warranty

ALIMED warrants to Purchaser that the Product(s) sold hereunder, is at the time of shipment to Purchaser from ALIMED, free from defects in material and workmanship. As Purchaser's sole exclusive remedy under this warranty ALIMED agrees either to repair or replace, at ALIMED's sole option and free of charge to Purchaser, any part or parts of such Product(s) which under proper and normal conditions of use prove to be defective within six months from the date of shipment. Visit [AliMed.com/terms](https://www.AliMed.com/terms) for full warranty info.

California Proposition 65 Statement:

65  **WARNING:** This product can expose you to chemicals including Di(2-Ethylhexyl) Phthalate (DEHP), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to <https://www.p65warnings.ca.gov>.

 <i>UK Responsible Person</i> MDSS-UK RP Ltd. Parkway House, Palatine Rd, Northenden, Wythenshawe, Manchester M22 4DB United Kingdom	  MDSS GmbH Schiffgraben 41 Hannover 30175 Germany
---	--

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©2026 AliMed®, LLC. Wszelkie prawa zastrzeżone.
AliMed® LLC