# ottobock.

**4R76**

**4R78**


## 4R76, 7R78

<b>DE</b>	Gebrauchsanweisung .....	3
<b>EN</b>	Instructions for use .....	5
<b>FR</b>	Instructions d'utilisation .....	7
<b>IT</b>	Istruzioni per l'uso .....	9
<b>ES</b>	Instrucciones de uso .....	12
<b>PT</b>	Manual de utilização.....	14
<b>NL</b>	Gebruiksaanwijzing.....	16
<b>SV</b>	Bruksanvisning .....	18
<b>DA</b>	Brugsanvisning.....	20
<b>NO</b>	Bruksanvisning .....	22
<b>PL</b>	Instrukcja użytkowania.....	24
<b>HU</b>	Használati utasítás .....	27
<b>CS</b>	Návod k použití.....	29
<b>HR</b>	Upute za uporabu .....	31
<b>TR</b>	Kullanma talimatı.....	34
<b>EL</b>	Οδηγίες χρήσης .....	36
<b>JA</b>	取扱説明書.....	38
<b>ZH</b>	使用说明书.....	40

**4R76**

①

**4R78**

②

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English

Last update: 2015-12-03

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

## 1 Indications for use

The 4R76 and 4R78 Double Adapters are to be **exclusively** used for the prosthetic fitting of amputations of the lower limb.

Field of application according to the **Ottobock Mobility System MOBIS**:

Approved for a patient weight of up to **150 kg/330 lbs**

## 2 Description

The 4R76 and 4R78 Ottobock Double Adapters provide the connection between two prosthetic components with pyramid receiver and add to the adjustment possibilities of the Ottobock Modular Lower Limb Prosthesis.

The **4R76 Double Adapter** (Fig. 1) allows for horizontal adjustment on the flattened side and for the proven angle adjustment on the convex side.

The **4R78 Double Adapter** (Fig. 2) allows for angle adjustment on both sides.

## 2.1 Technical data

Article no.	4R76	4R78
Weight	95 g	115 g
System height	-32 mm	-30 mm
Material	Stainless steel	Stainless steel
Max. patient weight	<b>150 kg / 330 lbs</b>	<b>150 kg / 330 lbs</b>
Mobility grade	–	–

### Attention!

Please avoid exposing prosthetic components to corrosive elements such as fresh water, salt water, and acids. Using this medical product in such environmental conditions will render all claims against Otto Bock HealthCare null and void.

**Please also inform your patients!**

## 3 Handling

### 3.1 Maintenance instructions

#### Note:

As a basic principle, all Ottobock modular adapters are subjected to tests involving three million load cycles. Depending on the amputee's activity this corresponds to a service life of three to five years.

We recommend carrying out regular safety checks once a year.

## 4 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 4.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

## 4.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

## 4.3 Warranty

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction and that are reported to the manufacturer within the warranty period.

Further information on the warranty terms and conditions can be obtained from the competent manufacturer distribution company.

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Otto Bock HealthCare GmbH