
## 5R1=1, 5R1=2, 5R1=6, 5R1=6-H

<b>DE</b> Gebrauchsanweisung .....	3	<b>CS</b> Návod k použití.....	26
<b>EN</b> Instructions for use .....	4	<b>RO</b> Instrucțiuni de utilizare.....	28
<b>FR</b> Instructions d'utilisation.....	6	<b>HR</b> Upute za uporabu .....	29
<b>IT</b> Istruzioni per l'uso .....	8	<b>SL</b> Navodila za uporabo .....	31
<b>ES</b> Instrucciones de uso .....	10	<b>SK</b> Návod na používanie .....	33
<b>PT</b> Manual de utilização.....	12	<b>BG</b> Инструкция за употреба .....	34
<b>NL</b> Gebruiksaanwijzing.....	14	<b>TR</b> Kullanma talimatı .....	36
<b>SV</b> Bruksanvisning .....	15	<b>EL</b> Οδηγίες χρήσης .....	38
<b>DA</b> Brugsanvisning .....	17	<b>RU</b> Руководство по применению .....	40
<b>NO</b> Bruksanvisning .....	19	<b>JA</b> 取扱説明書.....	42
<b>FI</b> Käyttöohje.....	21	<b>ZH</b> 使用说明书.....	44
<b>PL</b> Instrukcja użytkowania.....	22	<b>KO</b> 사용 설명서.....	45
<b>HU</b> Használati utasítás .....	24		

# 1 Produktbeschreibung

Deutsch

## INFORMATION

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##### **INFORMATION**

Last update: 2016-07-06

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

#### 1.1 Construction and Function

The 5R1=1, 5R1=2, 5R1=6 and 5R1=6-H socket attachment blocks are bonded to the distal end of the prosthetic socket and then laminated together with it. They serve to connect the prosthetic socket to a socket adapter.

#### 1.2 Combination possibilities

This prosthetic component is part of the Ottobock modular system and can be combined with other products of this modular system.

#### 2 Intended use

##### 2.1 Indications for use

The product is to be used solely for lower limb prosthetic fittings.

##### 2.2 Area of application

- **5R1=1, 5R1=2, 5R1=6:** Approved for a body weight of up to **125 kg**.
- **5R1=6-H:** Approved for a body weight of up to **150 kg**.

## 2.3 Environmental conditions

### Allowable environmental conditions

Temperature range for use: -10 °C to +60°C  
Allowable relative humidity 0 % to 90 %, non-condensing

### Unallowable environmental conditions

Mechanical vibrations or impacts  
Perspiration, urine, fresh water, salt water, acids  
Dust, sand, highly hygroscopic particles (e.g. talcum)

## 2.4 Service Life

This prosthetic component has been tested by the manufacturer for 3 million load cycles according to ISO 10328. Depending on the patient's activity level, this corresponds to a service life of 3 to 5 years.

## 3 Safety

### 3.1 Explanation of warning symbols

**⚠ CAUTION** Warning regarding possible risks of accident or injury.

### 3.2 General safety instructions

#### ⚠ CAUTION

##### Excessive strain on the product

Risk of injury due to breakage of load-bearing components  
► Use the product according to the specified area of application (see Page 4).

#### ⚠ CAUTION

##### Unallowable combination of prosthetic components

Risk of injury due to breakage or deformation of the product  
► Only combine the product with prosthetic components that are approved for that purpose.  
► Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

#### ⚠ CAUTION

##### Use under unallowable environmental conditions

Risk of injury due to damage to the product  
► Do not expose the product to unallowable environmental conditions.  
► If the product has been exposed to unallowable environmental conditions, check it for damage.  
► If damage is apparent or in case of doubt, do not continue using the product.  
► Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

#### ⚠ CAUTION

##### Exceeding the service life

Risk of injury due to change in or loss of functionality and damage to the product  
► Ensure that the approved service life is not exceeded.

#### ⚠ CAUTION

##### Mechanical damage to the product

Risk of injury due to change in or loss of functionality  
► Use caution when working with the product.  
► If the product is damaged, check it for proper function and readiness for use.  
► In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).  
► Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

## Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

## 4 Scope of delivery

see fig. 1, item	Quantity	Designation	Reference number
–	1	Instructions for use	–
1	1	Socket attachment block	5R1=1
2			5R1=2
3			5R1=6-H
6			5R1=6
4	1	Lamination dummy	4X6
5	4	Countersunk head screw	501S41=M6x-25

## 5 Preparation for use

### ⚠ CAUTION

#### Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components  
► Observe the alignment and assembly instructions.

### INFORMATION

The reinforcement described in this document was approved for the maximum product user body weight. The prosthetist assumes full responsibility for any change to the reinforcement.

## 5.1 Fabricating the Prosthetic Socket

### First lamination process

- > **Required materials:** 99B81 PVA bag, 616G6 Dacron® felt, 623T3 Perlon stockinette, 616G12 carbon fibre cloth, 617H119 Orthocryl lamination resin 80:20 PRO
- 1) Pull a soaked PVA bag over the model.
  - 2) Place a layer of Dacron® felt over the entire model.
  - 3) Place a layer of Dacron® felt around the proximal edge of the model.
  - 4) Pull 2 layers of Perlon stockinette over the model.
  - 5) Place 3 layers of carbon fibre cloth (e.g. 15 cm x 15 cm) on the distal end of the model, offsetting the fibre alignment.
  - 6) Pull 2 layers of Perlon stockinette over the model.
  - 7) Pull a soaked PVA bag over the model.
  - 8) Complete the lamination process with lamination resin.
  - 9) Once the lamination resin has cured, remove the PVA bag.

### Foaming and fastening the socket attachment block

- > **Required materials:** 617H32 Pedilen rigid foam 300, 617H21 Orthocryl sealing resin
- 1) Sand down the surplus lamination resin on the distal end of the model.
  - 2) Foam the distal end of the model with rigid foam.
  - 3) Trim the rigid foam to the maximum possible extent and sand it smooth in the desired position.
  - 4) Fasten the socket attachment block to the rigid foam with sealing resin.
  - 5) Mill a transition between the socket attachment block and the rigid foam (see fig. 2). While doing this, do not sand the shape or size of the adapter connection surface, which is determined by the socket attachment block.
  - 6) Insert the lamination dummy in the threaded holes of the socket attachment block (see fig. 3).
  - 7) Apply sealing resin to the outside of the prosthetic socket and to the lamination dummy (see fig. 4).

### Second lamination process

- > **Required materials:** 99B81 PVA bag, 623T3 Perlon stockinette, 616B1=50 carbon fibre cloth strap, 616G12 carbon fibre cloth, 616G15 woven carbon fibre stockinette, 617H119 Orthocryl lamination resin 80:20 PRO
- 1) Pull a Perlon stockinette over the model.
  - 2) Apply 2 strips of carbon fibre cloth strap in a cross over the socket attachment block up to the middle of the model.
  - 3) Place 2 layers of carbon fibre cloth (e.g. 15 cm x 15 cm) on the distal end of the model, offsetting the fibre alignment.
  - 4) Pull a Perlon stockinette over the model.
  - 5) Pull 2 layers of woven carbon fibre stockinette over the model. Do not allow the stockinette to protrude beyond the socket edge.

- 6) Pull 2 layers of Perlon stockinette over the model.
- 7) Pull a soaked PVA bag over the model.
- 8) Complete the lamination process with lamination resin.
- 9) Once the lamination resin has cured, remove the PVA bag.

## 5.2 Final assembly

### ⚠ CAUTION

#### Improper assembly of the screw connections

Risk of injury due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified torque values.
- ▶ Follow the instructions regarding the length of the screws and about how to secure the screws.

> **Required materials:** sanding tool, 710D1 torque wrench

- 1) **NOTICE! Do not sand down the laminate around the edge of the product.**

Sand down the laminate on the distal end of the prosthetic socket to the surface of the lamination dummy and remove the lamination dummy.

- 2) **CAUTION! Do not use thread lock (e.g. Loctite®) to mount the socket adapter. Chemical thread locks damage the laminate.**

Fasten the socket adapter to the socket attachment block using the countersunk head screws (see fig. 5):

Screw in 2 countersunk head screws on the posterior side (12 Nm).

Screw in 2 countersunk head screws on the anterior side (12 Nm).

## 6 Cleaning

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

## 7 Maintenance

- ▶ The prosthetic components should be inspected after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

## 8 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in the country of use regarding return, collection and disposal procedures.

## 9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

### 9.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

## 10 Technical data

Reference number	5R1=1	5R1=2	5R1=6	5R1=6--H
Material	Wood			
Wood connection diameter [mm]	147	120		

Reference number	5R1=1	5R1=2	5R1=6	5R1=6--H
Weight [g]	445	355	155	240
System height [mm]	Min. 46, max. 64		30	33
Max. body weight [kg]	125			150

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Otto Bock HealthCare GmbH