
4R98

DE	Gebrauchsanweisung	3
EN	Instructions for use	8
FR	Instructions d'utilisation	14
IT	Istruzioni per l'uso	20
ES	Instrucciones de uso	26
PT	Manual de utilização	32
NL	Gebruiksaanwijzing	37
SV	Bruksanvisning	43
DA	Brugsanvisning	49
NO	Bruksanvisning	54
FI	Käyttöohje	60
PL	Instrukcja użytkowania	65
HU	Használati utasítás	71
CS	Návod k použití	76
RO	Instrucțiuni de utilizare	82
HR	Upute za uporabu	88
SL	Navodila za uporabo	93
SK	Návod na používanie	98
BG	Инструкция за употреба	104
TR	Kullanma talimatı	110
EL	Οδηγίες χρήσης	115
RU	Руководство по применению	122
JA	取扱説明書	128
ZH	使用说明书	133
KO	사용 설명서	138

1

INFORMATION

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INFORMATION

Last update: 2015-10-15

► Please read this document carefully before using the product.

- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

1.1 Construction and Function

The 4R98 Sliding Tube Clamp Adapter permits translation adjustments at the proximal end of a tube adapter. It therefore permits the distal components of the modular prosthesis to be shifted in parallel, regardless of the angle adjustment by the pyramid adapters. The adjustments can be made either in the frontal plane (medial or lateral) or in the sagittal plane (anterior or posterior). This is also possible while wearing the prosthesis.

1.2 Combination possibilities

This prosthetic component is part of the Ottobock modular system. As a result of its design, it may be combined with other prosthetic components of the modular system. You can find possible restrictions in this section. Contact the manufacturer with questions.

2 Application

2.1 Indications for use

The product is to be used solely for lower limb prosthetic fittings.

2.2 Area of application

Approved for a body weight of up to **75 kg**.

Ottobock recommends the use of structural components with a diameter of 34 mm when elevated loads are expected (e.g. for a transtibial prosthesis).

2.3 Environmental conditions

Allowable environmental conditions
Temperature range for use: -10 °C to +60°C
Allowable relative humidity 0 % to 90 %, non-condensing
Unallowable environmental conditions
Mechanical vibrations or impacts
Perspiration, urine, fresh water, salt water, acids
Dust, sand, highly hygroscopic particles (e.g. talcum)

2.4 Service Life

In general, all modular adapters are tested by the manufacturer for 3 million load cycles. Depending on the patient's level of activity, this corresponds to a service life of 3 to 5 years.

3 Safety

3.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

3.2 General safety instructions

 **CAUTION**

Excessive strain on the product

Risk of injury due to breakage of load-bearing components

- ▶ Use the prosthetic components according to the classification (see section "Area of application").

 **CAUTION**

Unallowable combination of prosthetic components

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved according to the section "Combination possibilities".
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

 **CAUTION**

Use under unallowable environmental conditions

Risk of injury due to damage to the product

- ▶ Do not expose the product to unallowable environmental conditions (see the section "Environmental conditions").
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

 **CAUTION**

Exceeding the service life

Fall due to change in or loss of functionality and damage to the product

- Ensure that the approved service life is not exceeded (see section "Service life").

⚠ CAUTION

Mechanical damage to the product

Risk of injury due to change in or loss of functionality

- Use caution when working with the product.
- If the product is damaged, check it for proper function and readiness for use.
- In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

4 Scope of delivery

The following single components and accessories are included in the scope of delivery in the specified quantities and are available for reordering as single components (■), single components with minimum order quantities (▲) or single-component packs (●):

4R98 Sliding Tube Clamp Adapter					
Fig.	Item no.		Quantity	Designation	Reference number
–	–	■	1	Instructions for Use	647H66
2	1	▲	1	Cap Screw	501Z2=M6x25
2	2	▲	1	Cap Screw	501Z2=M6x35
2	3	▲	4	Set Screw	506G3=M8x16

5 Preparation for use

⚠ CAUTION

Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- Observe the alignment and assembly instructions.

CAUTION

Incorrect mounting of the tube

Fall due to breakage of load-bearing components

- ▶ Slide the tube all the way to the stop in the intended prosthetic component when mounting.

5.1 Mounting the adapter

> **Recommended tools and materials:**

710D4 Torque Wrench, degreasing cleaner (e.g. 634A3 Acetone)

- 1) Loosen the tube clamp cap screw by **2 turns**.
- 2) **For final assembly:** clean the insertion range of the tube clamp adapter and tube adapter with degreasing cleaner.
- 3) Fully insert the tube adapter into the tube clamp adapter up to the stop.
Align the clamping slot:
Tube clamp adapter: **anterior**
Tube clamp adapter, sliding: **anterior** or **medial**
- 4) Tighten the tube clamp cap screw (installation torque: M4: **5 Nm**, M5: **10 Nm**, M6: **13 Nm**).

5.2 Installation in the modular prosthesis

> **Required tools and materials:**

710D4 Torque Wrench, 636K13 Loctite®

- ▶ Perform the installation of the modular prosthesis and replacement of the set screws as described in the instructions for use of the connecting component.

5.3 Adjusting the adapter

The proximal section of the tube clamp adapter is movable in a dovetail guide. It is secured with a cap screw.

> **Required tools:**

710D4 Torque Wrench

- > No load on the prosthesis, or the prosthesis is not being worn.
- 1) Loosen the cap screw of the dovetail guide (see fig. 1).
 - 2) Move the proximal section of the sliding tube clamp adapter in the dovetail guide. Note the maximum displacement.
 - 3) Tighten the cap screw (assembly torque: **15 Nm**).

6 Cleaning

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.

- 3) Allow to air dry in order to remove residual moisture.

7 Maintenance

- ▶ The prosthetic components should be inspected after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

8 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in the country of use regarding return, collection and disposal procedures.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

9.3 Warranty

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction and that are reported to the manufacturer within the warranty period.

Further information on the warranty terms and conditions can be obtained from the competent manufacturer distribution company.

10 Technical data

Reference number	4R98
Weight [g]	150
System height [mm]	57
Material	Aluminium
Diameter [mm]	30
Max. displacement [mm]	±9
Max. body weight [kg]	75

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Otto Bock HealthCare GmbH