# GENOU 1P650

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
# 1P650 KNEE

## *Assembling and adjustment instructions for the CPO*

**The Easy Ride knee has been designed by handisport athletes to meet the specific needs of many sports activities**

### I. FEATURES AND ADVANTAGES

- Flexion and extension assist are adjustable thanks to the specific preparation of the shock absorber.
- Rotating adjustment of the upper connector.
- Can be used in every environment.
- Lightweight (1200g) and strong.
- Full height : 257 mm.
- Build height : 227 mm.
- Knee structure in anodised aeronautics aluminium alloy.
- Axis and screw fittings in stainless steel.

### II. PACK CONTENTS

(x1) Air pump

(x2) Screw M6x14 WMTX6014Z

(x2) T30 Torx wrench

(x1) Axis change tool

(x1) 1P650 EASY RIDE knee with its new generation air/oil shock absorber and its 2 connectors, the upper one being rotating

### III. USE

This knee joint has been designed to make possible sports activities for patients up to 125 kg (carried load included) with L4 activity level. Its design combining an air/oil shock absorber and a rigid structure enables the patient to practice sports activities both in summer (road biking, skateboard, surf) and in winter (snowboard).


- As the shock absorber has been lubricated during preparation, it may show slight traces of oil. This is normal and they should be wiped with a soft cloth.
- As the knee includes several moving items, please keep your hands away while functioning to avoid any risk of getting stuck.
- Wearing flotation equipment is highly recommended during any water activity. We also advise you to wear a light foot (in composite or EVA).
- We remind you that there is a risk of falling when cycling.

### IV. ASSEMBLING AND BUILD HEIGHTS

The upper and lower linking parts of this knee are two male pyramids. They receive a female pyramid of 1K185 or 1D40 type.

The male pyramids provided can be replaced by female pyramids (1K185 and 1K66) that will be ordered separately.

## VI. RECOMMENDED ALIGNMENTS


### Preparation :

Trace the build line and the knee joint center line (with patient standing, knee in complete extension).

### Sliding position :

Sagittal plane : Assemble the socket in a flexed position, so that the load line goes through the hip joint center, at about 3 cm (according to the activity) behind the front-upper axis of the knee and at the level of the joint between the metatarsus and the foot phalanges. The heel can get off the ground.


Frontal plane : Assemble the knee in line with the socket. Assemble the leg segment with a adduction angle of about 10° from the knee (according to the activity). The load line goes through the joint center of the socket and the medial side of the foot (according to the activity).

### Biking position :

Sagittal plane : Assemble the knee so that it has the same axis as the socket. The build line should go through the hip joint center, the socket center and behind the rotation axis of the knee. It goes approximately through the third of the foot (for further information refer to the foot alignment recommendations).

Frontal plane : Assemble the knee so that it has the same axis as the socket. The build line goes through the knee center and the foot center. A medial rotation of the leg segment may be necessary.


Nota : These recommendations are based on laboratory tests and trials. When the patient gets into the real situation of its sport activities, these theoretical alignments can be modified provided that the safety rules and the rules of good use for the assembly of this type of material be respected. The patient's safety and comfort depend on them.

## VII. CHANGING POSITION

To switch from SLIDING mode to BIKING mode, you have to change the upper axis of the cylinder from position 1 to position 2 as follows :

- Undo the two screws on either side of the axis using the two T30 Torx wrenches provided.
- Remove the axis from its place using the provided axis change tool.
- Insert the axis in the required place.
- Re-position the two screws with their washers and with a drop of XC048 glue, then tighten them up to a 6 Nm torque.

- To tighten and to loose the screws, please use a **torque wrench** fitted with a T30 Torx tip.
- In case of deformation of the impression on the screw head, replace this screw by one of the two screws provided.
- The maximal flexion stops in BIKING position aim only at **protecting the cylinder** and can in no case be subject to regular and repeated stresses.

## VIII. SELECTION OF THE POSITION AND AJUSTMENTS ACCORDING TO THE ACTIVITY

The cylinder has three controls for adjustment:

- ⇒ **Blue control** (FIRM / MEDIUM /OPEN) for respectively low, medium and fast compression speed.
- ⇒ **Black control** (ADJUSTING TRAIL MODE) for fine adjustment of compression speed. It can be used in OPEN Mode only and has 3 positions : 1, 2 and 3 from soft to firm.
- ⇒ **Red control** (REBOUND) for the speed of extension assist.

## IX. ADJUSTING PRESSURE IN THE SHOCK ABSORBER

Use the pump provided to set the pressure required for the shock absorber. If pressure is too high, it can be decreased with the purge button located on the pump.

Use only the pump provided or a pump that is equipped with a 20 bar manometer and adapted to FOX shock absorbers.

The combinations of positions and adjustments shown in the table below are given as a guide only and may vary according to the patient's activities and weight.

<u>Use</u>	<u>Position</u>	<u>Shock absorber</u>		
		<u>Extension adjustment (REBOUND)</u>	<u>Flexion adjustment (OPEN MEDIUM FIRM)</u>	<u>Air pressure</u>
Road biking / Cross-country MTB	BIKING	MEDIUM / MAX	OPEN	1 to 3 bars
Downhill MTB	BIKING	MAX	MEDIUM / FIRM	4 to 6 bars
Snowboard / Track ski	SLIDING	MAX	MEDIUM / FIRM	4 to 6 bars
Snowboard / Off track ski	SLIDING	MEDIUM / MAX	FIRM	5 to 7 bars
Water ski / Wakeboard	SLIDING	MEDIUM	OPEN	1 to 3 bars
Surfing / Paddling	SLIDING	MEDIUM	MEDIUM	1 to 4 bars
Fitness/Rackets sports/Fightening	SLIDING	MEDIUM / MAX	FIRM	5 to 12 bars
Motocross	BIKING	MAX	FIRM	5 to 9 bars


For further information about shock absorber, please visit <http://www.ridexox.com/OwnersManuals>.

- It is possible that traces of oil appear during the first uses of the knee. Wipe them with a soft cloth.
- Air pressure in the shock absorber should never exceed **20,6 bars**.
- The **air in the shock absorber should necessarily be cleared** at high altitude (travel by plane, ...).
- After adjusting the air pressure, do not forget to **put the cap on the valve again**.
- As the shock absorber has undergone a specific preparation, **never dismantle it** and never use it for another purpose.
- Let the shock absorber cool after use. **Risk of burns**.
- If the temperature of the shock absorber varies within a wide range, its behaviour can differ.
- The shock absorber is delivered deflated to avoid any problem during transport.

## X. MAINTENANCE


- Regularly check the state of the flexion stop screws.
- Regularly check the tightening of the M6 screws (6Nm tightening torque).
- Clean the knee after each "polluting" activity (water, mud, underwater use).
- After immersion in sea water, rinse thoroughly with fresh water, then dry completely. Do not use any solvent or other aggressive cleaning products.
- In case of malfunction, premature wear or rust stains, contact your CPO as soon as possible.
- For the maintenance of the shock absorber, refer to the instructions sheets of the manufacturer supplied in the knee packaging.
- An inspection of the knee must be made by Proteor every 100 hours of use (ie each year for a moderate activity and every 6 months for a high level activity).


## XI. SPARE PARTS

During the normal life cycle of the EASY RIDE sports knee, ageing may be faster for some parts than for others according to the activity performed. Therefore PROTEOR supplies replacement kits for the most concerned parts.

### 1P652 (KIT OF AXIS, BUSHINGS AND CYLINDER)

### 1P653 (COVER KIT)

## XII. ASSEMBLING THE KITS

### 1P652 (KIT OF AXIS + SHOCK ABSORBER)

- Take away the shock absorber by removing the M6 screws on one side of the knee using the two T30 Torx wrenches provided. If a screw is hard to remove, heat it cautiously with a heat gun.

**Remark :** The T30 Torx wrench provided is compatible only with the M6 screws that must be removed. The other screws are Allen screws.

- Carrefully remove the flange that is now released from its screws.

If the flange cannot be removed easily, tap with a plastic mallet alternately on both ends.

- Remove the axis (cylinder-plate) (on the left on the picture) by tapping on the provided axis change tool with a plastic mallet.
- Remove with your hand the upper plate and the cylinder, then remove the axis (flanches) (on the right on the picture) by undoing the last screw that holds it.

Do not use a tool because this might damage the plate or the cylinder. Pull with your hand as close as possible to the axis while making small circular motions.

- Replace successively:
  - ◇ The provided axis (flanges).
  - ◇ The assembled cylinder (Caution : for an easy rotation of the cylinder a slight adjustment of the bushing length with sandpaper may be necessary).
  - ◇ The upper plate that was previously removed.
  - ◇ The provided axis (cylinder-plate), using the axis change tool.
  - ◇ The flange that was previously removed.
- Replace the previously removed screw fittings by the one provided for use in the event of damage. Do not forget to coat them with XC047 glue and to tighten them up to a 6 Nm torque.

The screws should be tightened with a torque wrench fitted with a T30 Torx tip.

## 1P653 (COVER KIT)

- Remove the 8 screws M3x8\* using an Allen wrench of size 2, then unstick the three shells (it may be necessary to use a solvent such as acetone).

- Clean glue residues on the surfaces of the flanges and the upper plate.
- Use a dry solvent to degrease the bonding surfaces (aluminium surfaces and covers), then coat one of the two surfaces with cyanoacrylate glue XC053 (a thread along the border of the surface is generally sufficient).
- Glue the positioned covers and apply an even pressure during at least 3 minutes

Once the two surfaces are in contact, do not carry out micro-motions and do not re-position the parts.

- Use XC047 glue to position the 8 washers and the 8 screws provided. Tighten up to a 0.25Nm torque.

- The screws should be tightened with a torque wrench fitted with an hexagonal tip of size 2.
- After gluing, the knee should dry for at least 24 hours before use.

### XIII. RECYCLING

This product includes main parts in aluminium, axis and screws in stainless steel, connectors in titanium and bushes in plastics. It also contains oil. These items must be recycled according to the laws in force.

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