ottobock.


CE

3R106=P

| | | |
|----|------------------------------------|-----|
| DE | Gebrauchsanweisung | 10 |
| EN | Instructions for use | 18 |
| FR | Instructions d'utilisation | 26 |
| IT | Istruzioni per l'uso | 35 |
| ES | Instrucciones de uso | 44 |
| PT | Instruções de utilização | 53 |
| NL | Gebruiksaanwijzing | 61 |
| SV | Bruksanvisning | 69 |
| DA | Brugsanvisning | 77 |
| NO | Bruksanvisning | 85 |
| FI | Käyttöohje | 93 |
| PL | Instrukcja użytkowania | 101 |
| HU | Használati útmutató | 110 |
| CS | Návod k použití | 118 |
| TR | Kullanım kılavuzu | 126 |
| EL | Οδηγίες χρήσης | 135 |
| RU | Руководство по применению | 144 |
| UK | Інструкція щодо використання | 153 |
| JA | 取扱説明書 | 162 |
| ZH | 使用说明书 | 169 |

1

1 3R106=P

2 2R49

3 2Z11=KIT

4 4G1244=1


5 647G1865

2


3


743Y580=*

4


704Y14*


719R5

2


718R1

3


5

1

X

2


6


7


8


X

9



10



11




3D-Mode

12


2D-Mode

15

1

2

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INFORMATION

Date of last update: 2025-11-07

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

The product has the following key features:

- Polycentric prosthetic knee joint
- Product components for stance phase stability
 - Polycentric kinematics
- Large flexion angle
 - Up to a maximum of 172°
- Auto-adaptive swing phase control
 - Enables walking at different speeds without having to adapt swing phase control
- Product components for swing phase control
 - Flexion and extension damping can be adjusted separately and individually

2 Safety-related information

CAUTION! Risk of injury and risk of product damage

- ▶ To prevent mechanical damage, treat the product with care.
- ▶ Observe the information on product combinations in the instructions for use for the products.
- ▶ Observe the maximum lifetime of the product (see page 20).
- ▶ Observe the product's operating conditions (see page 19).
- ▶ Do not expose the product to prohibited environmental conditions (see page 19).
- ▶ Check the product for damage and readiness for use prior to each use.
- ▶ Do not use the product if it is damaged or if you are not sure whether it is working properly. Have the product inspected, repaired or replaced by the manufacturer or a specialist workshop.
- ▶ The product is intended for use by one person only and may not be reused by other persons.
- ▶ To avoid the risk of pinching, do not reach into the joint mechanism.
- ▶ Observe the instructions for cleaning and care (see page 25).
- ▶ Follow the maintenance instructions (see page 25).

3 Intended use

3.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

3.2 Operating conditions

- Maximum permissible body weight: **100 kg**

The product is recommended for:

- Mobility grade 2: restricted outdoor walker
- Mobility grade 3: unrestricted outdoor walker

3.3 Target group

This document is intended for both the qualified personnel that ordered the product and the user and their attendant.

All assembly tasks and settings in the section **Preparing for use** (see page 20) may only be carried out by authorised, qualified personnel. All manufacturer specifications and all applicable legal provisions must be complied with. Please contact the manufacturer's service department for further information.

3.4 Environmental conditions

| Storage and transport |
|--|
| Storage temperature: -20 °C to +60 °C, relative humidity: 20 % to 90 %, no vibrations or impacts |

| Allowable environmental conditions |
|--|
| Temperature for use: -10 °C to +45 °C |
| Moisture: relative humidity: 20 % to 90 %, non-condensing |

| Prohibited environmental conditions |
|---|
| Chemicals/liquids: fresh water, salt water, perspiration, urine, acids, soapsuds, chlorine water |
| Solids: highly liquid-binding particles (e.g. talcum), dust, sand |

3.5 Product combinations

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

Permissible and impermissible product combinations with the 3R106=P prosthetic component are listed in the tables (see page 7-8, see page 9).

3.6 Lifetime

This prosthetic component has been load tested by the manufacturer according to ISO 10328. The maximum lifetime is 5 years.

4 Scope of delivery

| 3R106=P (see fig. 1) | | |
|----------------------|------------------|-----------------------|
| Item | Reference number | Designation |
| ① | 3R106=P | Prosthetic knee joint |
| ② | 2R49 | Tube adapter |
| ③ | 2Z11=KIT | Adapter protector |
| ④ | 4G1244=1 | Stop bumper |
| ⑤ | 647G1865 | Instructions for use |

5 Accessories

The following components can be ordered separately as needed:

| 3R106=P | |
|------------------|-------------|
| Reference number | Designation |
| 4G1244=1 | Stop bumper |

6 Preparing the product for use

CAUTION

Incorrect alignment, assembly or adjustment

Injury due to incorrectly installed or adjusted as well as damaged prosthetic components

- Observe the alignment, assembly and adjustment instructions.

6.1 Optional: Using the 2Z11=KIT adapter protector

The adapter protector protects the connecting section of the prosthetic joint from scratches during alignment and fitting.

> The following instructions refer to the illustration: see fig. 2

- 1) Follow the instructions in the accompanying document for the adapter protector.
- 2) Remove the adapter protector before delivery of the product.
- 3) Readjust the prosthetic alignment.

6.2 Optional: Using the 743A220 PROS.A.* Assembly alignment apparatus

An alignment apparatus such as the 743A220 PROS.A. Assembly can be used for the alignment of modular leg prostheses.

> The following instructions refer to the illustration: see fig. 3

- 1) Follow the instructions in the accompanying document for the alignment device.
- 2) Mount the 743Y580=* adapter inserts.
- 3) Carry out the alignment of the prosthesis.

6.3 Bench alignment of the prosthesis

The prosthesis must be aligned using the alignment images (see fig. 9; see fig. 10) as well as the instructions for use of all prosthetic components used. The prosthesis should be aligned so as to enable the user to stand safely.

- > **Required materials and tools:** alignment apparatus (e.g. 743A220 PROS.A. Assembly), 743A80 50:50 gauge, 743S12 heel height measuring device
- 1) Position the prosthetic foot in the alignment apparatus.
- 2) **Heel height**
 - The effective heel height of everyday shoe (x) **(+ 5 mm)**.
- 3) **Frontal plane**
 - Align the prosthetic foot with the alignment reference line: **Position the middle of the prosthetic foot on the alignment reference line (external rotation approx. 5 – 7°)**.
- 4) **Sagittal plane**
 - Align the prosthetic foot with the alignment reference line: **Place the alignment reference line behind the middle of the prosthetic foot (+ 30 mm)**.
- 5) Position the prosthetic knee joint and prosthetic socket in the alignment apparatus.
- 6) **Frontal plane**
 - Align the prosthetic socket with the alignment reference line.
 - Take the abduction or adduction of the residual limb into account.
- 7) **Sagittal plane**
 - Align the prosthetic socket with the alignment reference line.
 - Determine the centre of the prosthetic socket proximally and distally with the 50:50 gauge in the direction of movement and mark the centre line.
 - Mark the socket reference point on the centre line proximal to the tuberosity reference point **(30 mm distance)**.
 - Align the prosthetic socket so the alignment reference line runs vertically through the socket reference point.
 - Adjust the socket flexion by rotating around the socket reference point: **individual residual limb flexion/flexion contracture (+ 5°)**.

6.4 Shortening the tube adapter

CAUTION

Incorrect processing of the tube

Risk of injury due to damage to the tube

- ▶ Do not clamp the tube in a vice.
- ▶ To shorten the tube, use only a tube cutter or a cutting device.

- > The following instructions refer to the illustration: see fig. 4
- > **Recommended tools:** 704Y14=* cutting device or 719R5 tube cutter, 718R1 tube deburring knife
- 1) Shorten the tube adapter with a cutting device or a tube cutter ①.
- 2) Deburr the cut edge with a tube deburring knife ②.
- 3) Make sure that the tube is evenly deburred on the inside and outside ③.

6.5 Installing the tube adapter

CAUTION

Improper assembly of the screw connections

Risk of injury due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified torque values.
- ▶ Follow the instructions for thread lock.

- > The following instructions refer to the illustration: see fig. 5
- > **Required materials and tools:** 634A58 isopropyl alcohol, 710D20 torque wrench, hexagon socket screw adapter 4 mm
- 1) Clean the contact surfaces of the tube and of the tube receiver using a degreasing cleaning agent (such as isopropyl alcohol).
- 2) **CAUTION! Insert the tube adapter into the tube receiver up to the stop. Otherwise, components may become overloaded and break.**
Insert the tube adapter into the tube receiver up to the stop (0 mm)①.
- 3) Tighten the cap screw **without** additional Loctite to **10 Nm**②.

6.6 Static Alignment

Ottobock recommends optimising the alignment with the 743L500 3D L.A.S.A.R. Posture measuring device. The 743L100 L.A.S.A.R. Posture can be used as an option.

The following optimisation methods are available:

- 3D mode
- 2D mode

Optimising the prosthetic alignment:

Make sure that the user is standing comfortably and that the reference points in the illustrations are observed.

Noting the lines:

Dotted lines = auxiliary lines

Solid lines = load lines

Maintaining distances:

The distances between the auxiliary and load lines shown in the illustrations have to be maintained.

3D mode (see fig. 11)

- ① **Prerequisite:** The prosthetic knee joint is fully extended.
- ② The auxiliary line lies on the sagittal reference point of the prosthetic knee joint (pivot point).
- ③ The auxiliary line lies on the frontal reference point (anterior superior iliac spine).
- ④ The auxiliary line lies on the frontal reference point of the prosthetic knee joint (pyramid).

2D mode (see fig. 12)

- ① **Prerequisite:** The prosthetic knee joint is fully extended.
- ② The auxiliary line lies on the sagittal reference point of the prosthetic knee joint (pivot point).
- ③ The auxiliary line lies on the frontal reference point (anterior superior iliac spine).
- ④ The auxiliary line lies on the frontal reference point of the prosthetic knee joint (pyramid).
- ⑤ The auxiliary line lies on the frontal reference point of the prosthetic foot (pyramid).

6.7 Dynamic Trial Fitting

The optimal gait pattern is developed during dynamic fitting. The prosthesis is aligned optimally in the frontal plane and the sagittal plane for this purpose.

6.7.1 Restoring factory settings

| 3R106=P factory defaults | | | |
|--------------------------------|---|---|---|
| Function | Adjustment position | Factory settings | Meaning |
| Swing phase flexion resistance | Adjustment screw F (see fig. 6) | Turned clockwise from left-hand stop (number of turns: 2) | Adjustment screw closed slightly – low resistance |

| 3R106=P factory defaults | | | |
|----------------------------------|---|------------------|--------------------|
| Function | Adjustment position | Factory settings | Meaning |
| Swing phase extension resistance | Adjustment screw E (see fig. 7) | Left-hand stop | Minimum resistance |

6.7.2 Checking the settings

> **Prerequisite:** The product is set to factory settings on delivery.

- 1) Check the settings of the product prior to the dynamic fitting.
- 2) In case of deviations, reset the product to the factory settings.

6.7.3 Adjusting the swing phase flexion resistance

> The following instructions refer to the illustrations: see fig. 6; see fig. 8

> **Required materials and tools:** hexagon socket screw 2 mm

- 1) Use the hexagon socket screw to adjust the swing phase flexion resistance to the user's gait pattern.
 - Make sure that the prosthetic foot swings through sufficiently but not too far.
- 2) **For more resistance**, turn adjustment screw **F** clockwise .
- 3) **For less resistance**, turn adjustment screw **F** anticlockwise.

6.7.4 Adjusting the swing phase extension resistance

> The following instructions refer to the figures: see fig. 7; see fig. 8

> **Required materials and tools:** hexagon socket screw 2 mm

- 1) Use the hexagon socket screw to adjust the swing phase extension resistance to the user's gait pattern.
 - **CAUTION! The prosthetic knee joint has to reach full extension, even at slow walking speeds and with small steps.**
Make sure that the prosthetic knee joint does not swing too gently against the extension stop.
- 2) **For more resistance**, turn adjustment screw **E** clockwise .
- 3) **For less resistance**, turn adjustment screw **E** anticlockwise.

6.8 Notes for socket stop

> **Checking the socket design (see fig. 13)**

▶ Check whether the socket touches the product at maximum flexion (e.g. when kneeling) in the area of the swing phase adjustment screws.

→ If the swing phase adjustment screws are not touched, the socket design is OK.

> **Optional: adaptation of the socket design (see fig. 13)**

▶ If contact between the socket and the product cannot be avoided, the socket is designed so that the product lies flat in the area of the tube clamp.

→ The socket design is optimised.

6.9 Finishing the prosthesis

Finish working on the prosthesis as described in **Bench alignment of the prosthesis** (see page 20) and **Static alignment** (see page 22) as well as the reference image (see fig. 14) and the following instructions.

- 1) Replace set screws that are too short with longer ones to prevent components from coming loose accidentally.
- 2) Replace set screws that are too long with shorter ones to prevent collisions and product damage.
- 3) Clean the threads of the screw connections on the prosthetic socket and foot using a degreasing cleaner.
- 4) Only use Loctite 241 – 636K13 on the screw connections of the prosthetic socket and foot.
- 5) Finish all screw connections. Observe the assembly tightening torques and thread locks for all prosthetic components.

→ The prosthesis is completed.

6.9.1 Optional: Installing the foam cover

CAUTION

Use of liquid-binding substances (e.g. talcum)

Risk of injury, damage to the product due to lack of lubrication

- ▶ Prevent contact of the product with liquid-binding substances.

> **Prerequisite:** Use the **3S107** foam cover for the product.

> **Required materials and tools:** 519L5 silicone spray

- 1) Spray the silicone spray directly onto the contact surfaces in the foam cover to improve the gliding characteristics and reduce noise.
- 2) Attach the foam cover to the prosthesis.
 - After attaching the foam cover, perform another **dynamic fitting** (see page 22) to check the settings.
- 3) Finish the prosthesis.

7 Use

7.1 Information for use

▶ **Prior to use**, check that **full flexion**, **extension** and **damping resistance** are working.

→ In case of insufficient functionality, have the prosthesis adjusted by qualified personnel.

Information on flexion angle

The prosthetic knee joint automatically pushes against the flexion stop during kneeling.

The prosthetic foot can lift itself off the floor. This is normal and due to the design.

7.2 Walking on stairs or ramps

CAUTION

Walking step-over-step on stairs or ramps not possible

Risk of falling due to unintentional buckling of the prosthetic knee joint

- ▶ Only walk on stairs or ramps with the prosthetic knee joint extended.

Walking step-over-step is not possible because the prosthetic knee joint flexes in an uncontrolled manner due to low flexion and the simultaneous load.

The prosthetic knee joint provides optimum stability only in full extension.

Walking down stairs or ramps

> **Prerequisite:** Practise walking correctly on stairs or ramps.

- 1) Extend the prosthetic knee joint in front of the stairs or ramp.
- 2) Hold on to the railing for safety.
- 3) First set the prosthetic leg down in the extended state.
- 4) Repeat with the sound leg.
- 5) Maintain the correct sequence until the end.

→ The user should be able to walk safely on stairs and ramps.

Walking up stairs or ramps

> **Prerequisite:** Practise walking correctly on stairs or ramps.

- 1) Extend the prosthetic knee joint in front of the stairs or ramp.
- 2) Hold on to the railing for safety.
- 3) Take the first step with the sound leg.
- 4) Then continue with the prosthetic leg in the extended state.
- 5) Maintain the correct sequence until the end.

→ The user should be able to walk safely on stairs and ramps.

7.3 Cleaning

| Permissible cleaning agents and disinfectants | |
|---|---|
| Requirements: | <ul style="list-style-type: none">• Solvent-free• Chlorine-free• Phosphate-free |
| Permissible cleaning agents: | <ul style="list-style-type: none">• pH-neutral soap |
| Permissible disinfectants: | <ul style="list-style-type: none">• Water-based disinfectant |

Cleaned using soap and water

> **Required materials and tools:** 2 microfibre cloths, pH-neutral soap (e.g. 453H10 Derma Clean), tap water

- 1) **In case of light soiling (e.g. small deposits):** Moisten the microfibre cloth with water.
- 2) **In case of heavy soiling (e.g. stubborn deposits):** Moisten the microfibre cloth with pH-neutral soap and water.
- 3) Clean the product with the moistened cloth.
- 4) Dry the product with a dry microfibre cloth.
- 5) Allow to air dry in order to remove residual moisture.

Cleaning with Disinfectants

> **Required materials and tools:** 2 microfibre cloths, water-based disinfectant

- 1) Apply disinfectant to the microfibre cloth.
- 2) **NOTICE! Verify material compatibility when using disinfectants.**
- 3) Disinfect the product with the moistened cloth.
- 4) Dry the product with a dry microfibre cloth.
- 5) Allow to air dry in order to remove residual moisture.

8 Maintenance

- 1) **Every 12 months:** Have maintenance performed by qualified personnel. Reduce the maintenance interval based on the use of the prosthesis.
- 2) Check the prosthesis for **wear and tear, functionality and unusual noises**.
→ In case of insufficient functionality, have the prosthesis adjusted by qualified personnel.
- 3) Do not lubricate or grease the prosthetic joint.
→ Repair work must be performed exclusively by manufacturer service.

8.1 Replacing the stop bumper

Replace the stop bumper in case of visible, pronounced degeneration and/or loss of comfort (see fig. 15).

> **Required materials and tools:** flat screwdriver, compressed air, 4G1244=1 stop bumper

- 1) Move the product to maximum flexion ❶.
- 2) Carefully lever out the stop bumper using the flat screwdriver ❷.
- 3) Remove dirt and residues with compressed air.
- 4) Insert the new stop bumper.
- 5) Bring the product to maximum extension.
→ The product has a damped stop at full extension.

9 Disposal

Do not dispose of the product in household waste. Improper disposal is harmful to health and the environment. Observe the instructions of the local authority regarding return, collection and disposal.

10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

The full text of the EU Declaration of Conformity is available at the following Internet address: <http://www.ottobock.com/conformity>

11 Technical data

| 3R106=P | |
|--|---------|
| Proximal connection | Pyramid |
| Distal connection [mm] | 30 |
| Knee flexion angle [°] | 172 |
| Total weight [g] | 730 |
| Overall height [mm] | 222 |
| Build height [mm] | 207 |
| Mobility grade | 2, 3 |
| Permissible body weight, maximum [kg] | 100 |

| 2R49 | |
|-----------------------------------|------------------|
| Distal connection | Pyramid receiver |
| Weight [g] | 255 |
| Maximum body weight [kg] | 125 |
| Diameter [mm] | 30 |
| Minimum system height [mm] | 97 |
| Maximum system height [mm] | 472 |
| Material | Aluminium |

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