## 1E91 Runner, 1E93 Runner junior 1E97 Sprinter junior

<b>DE</b>	Gebrauchsanweisung .....	5
<b>EN</b>	Instructions for use .....	13
<b>FR</b>	Instructions d'utilisation .....	21
<b>IT</b>	Istruzioni per l'uso .....	29
<b>ES</b>	Instrucciones de uso .....	38
<b>PT</b>	Instruções de utilização .....	46
<b>NL</b>	Gebruiksaanwijzing .....	55
<b>SV</b>	Bruksanvisning .....	63
<b>DA</b>	Brugsanvisning .....	71
<b>NO</b>	Bruksanvisning .....	79
<b>FI</b>	Käyttöohje .....	86
<b>PL</b>	Instrukcja użytkowania .....	95
<b>HU</b>	Használati útmutató .....	103
<b>CS</b>	Návod k použití .....	111
<b>RO</b>	Instrucțiuni de utilizare .....	119
<b>HR</b>	Upute za uporabu .....	127
<b>SL</b>	Navodila za uporabo .....	136
<b>SK</b>	Návod na použitie .....	143
<b>BG</b>	Инструкция за употреба .....	151
<b>TR</b>	Kullanım kılavuzu .....	160
<b>EL</b>	Οδηγίες χρήσης .....	168
<b>RU</b>	Руководство по применению .....	177
<b>JA</b>	取扱説明書 .....	185
<b>ZH</b>	使用说明书 .....	193
<b>KO</b>	사용 설명서 .....	200

## 1E91, 1E93 Bench alignment, transtibial


### Socket

Sport socket: Flexion of everyday socket

### 1E91, 1E93

0 mm

Alignment reference line

## 1E91, 1E93 Bench alignment, transfemoral

### Socket

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### INFORMATION

Date of last update: 2024-11-25

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

The 1E91 Runner, 1E93 Runner junior and 1E97 Sprinter junior prosthetic feet are intended for use in sports prostheses.

The carbon spring contour provides high propulsion and low resistance. The prosthetic foot is also distinguished by its low weight.

This document also contains information about individual accessory products (foot adapters, outer soles, protectors).

## 2 Intended use

### 2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

## 2.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

- **The prosthetic socket must be able to withstand the higher loading when practising sports.**

Permitted combinations			
		Designation	Reference number
Outer soles and protectors	1E91	Runner sole	2Z540=6
		Runner spike sole	2Z541=6
	1E93, 1E97	Runner junior sole	2Z543=6
Foot adapter	1E91	Runner adapter with four-hole connector	4R216=6
		Runner adapter with pyramid, rotatable	4R218=6
	1E93	Runner junior adapter with pyramid, rotatable	4R224=6
	1E97	Posterior connection plate (set)	4R420
Prosthetic knee joints	1E91, 1E93	Sports knee joint	3S80
		Prosthetic sport knee joint with low viscosity oil	3S80=1

### 1E91

- **Only use prosthetic components approved for the respective sport or a body weight of 150 kg.**

### 1E93, 1E97

- **Only use prosthetic components approved for the respective sport or a body weight of 100 kg.**

## 2.3 Area of application

- The maximum approved body weight is specified in the technical data (see page 20).

This product is **not** suitable as an everyday prosthesis.

The product is **not** suitable for the long jump or sports with a similar load on the prosthesis.

## 1E93, 1E97

Approved for a body height up to **max. 145 cm.**

### 2.4 Environmental conditions

<b>Storage and transport</b>
Storage temperature: -20 °C to +60 °C, relative humidity: 20 % to 90 %, no mechanical vibrations or impacts
<b>Allowable environmental conditions</b>
<b>Temperature for use:</b> -10 °C to +45 °C
<b>Moisture:</b> relative humidity: 20 % to 90 %, non-condensing
<b>Chemicals/liquids:</b> fresh water as dripping water, occasional contact with salty air (e.g. near the ocean)
<b>Solids:</b> dust
<b>Prohibited environmental conditions</b>
<b>Chemicals/moisture:</b> salt water, perspiration, urine, acids, soapsuds, chlorine water
<b>Solids:</b> highly liquid-binding particles (e.g. talcum), dust in high concentrations (e.g. construction site), sand

### 2.5 Lifetime

#### Prosthetic foot, foot adapter

**Recreational sports:** The maximum lifetime is 2 years.

**Competitive sports:** The lifetime is 1 year.

#### Outer sole, protector

The product is a wear part, which means it is subject to normal wear and tear.

### 3 General safety instructions

#### CAUTION! Risk of injury and risk of product damage

- ▶ To prevent mechanical damage, use caution when working with the product.
- ▶ Note the combination possibilities/combo exclusions in the instructions for use of the products.
- ▶ Observe the maximum lifetime of the product.
- ▶ Check the product for damage and readiness for use prior to each use.
- ▶ Do not expose the product to excessive strain (section: "Area of application" see page 14) or to prohibited environmental conditions (see page 15).

- ▶ Do not use the product if it is damaged or in a questionable condition. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).
- ▶ The product may only be used by one person; reuse on another person is prohibited.

### Signs of changes in or loss of functionality during use

Decreased forefoot resistance or changes in roll-over behaviour are noticeable indications of loss of functionality.

## 4 Scope of delivery

Quantity	Designation
1	Instructions for use
1	Prosthetic foot

### Spare parts/accessories

Designation	Reference number
Hexagon socket screw (for foot adapter)	501S89=M8X25

## 5 Preparing the product for use

### CAUTION

#### Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- ▶ Observe the alignment and assembly instructions.

### 5.1 Selecting the stiffness

#### 1E91

Otto Bock recommends that the stiffness version is selected according to body weight and running style.

Body weight	Long-distance running	Sprint
	Stiffness version	
<b>40 kg to 50 kg</b> (90 lbs to 110 lbs)	<b>SPR-1</b>	<b>SPR-2</b>
<b>50 kg to 60 kg</b> (110 lbs to 130 lbs)	<b>SPR-2</b>	<b>SPR-3</b>
<b>60 kg to 72 kg</b> (130 lbs to 160 lbs)	<b>SPR-3</b>	<b>SPR-4</b>
<b>72 kg to 86 kg</b> (160 lbs to 190 lbs)	<b>SPR-4</b>	<b>SPR-5</b>

Body weight	Long-distance running	Sprint
	Stiffness version	
86 kg to 104 kg (190 lbs to 230 lbs)	SPR-5	SPR-6

### 1E93, 1E97

The stiffness of the prosthetic foot is selected depending on the body weight.

Body weight	Stiffness version
15 kg to 20 kg (35 lbs to 44 lbs)	SPR-1
20 kg to 25 kg (44 lbs to 55 lbs)	SPR-2
25 kg to 30 kg (55 lbs to 66 lbs)	SPR-3
30 kg to 37 kg (66 lbs to 81 lbs)	SPR-4
37 kg to 45 kg (81 lbs to 100 lbs)	SPR-5

### 5.2 For Runner only: Installing the adapter

The prosthetic foot is connected to the proximal prosthetic components with a connection adapter. The connection adapter is constructed in two parts. The prosthetic foot is secured between the upper and lower adapter sections. The connection adapter has either a rotating pyramid adapter or a four-hole connector for mounting of a socket adapter.

> **Recommended materials:** Torque wrench 710D20, Loctite 241 636K13

- 1) Set the upper section of the adapter onto the prosthetic foot. Align the adapter on the prosthetic foot in the middle of the scale.
- 2) Apply the lower section of the adapter on the opposite side.
- 3) Apply Loctite to the two hexagon socket screws and screw them in (**12 Nm**).

### 5.3 Installing the outer sole

#### CAUTION

#### **Incorrect use of prosthetic foot and outer sole**

Risk of injury due to lack of surface grip and damage to the product

- ▶ Only use the prosthetic foot with an outer sole.
- ▶ Only use outer soles with spikes on sports floors suitable for this purpose.

> **Recommended materials:** Loctite 241 636K13

- 1) Apply Loctite to the thread.
- 2) Set the sole against the prosthetic foot.
- 3) Secure the sole with the screws provided.

## **Spike sole only: Installing spikes**

> **Recommended materials:** Spike key 710Z3, spikes

- 1) Using a spike key, screw the spikes completely into the threads on the outer sole and hand-tighten them.
- 2) Check the tightness of the spikes before each use.

## **5.4 Alignment**

### **Bench alignment**

- ▶ **Sagittal plane:** Perform the bench alignment of the prosthesis (see illustrations at the beginning of the document).
- ▶ **Frontal plane:** The alignment reference line runs centrally through the prosthetic foot, centrally through the entrance plane of the prosthetic socket for transfemoral prostheses, and at the lateral edge of the patella for transtibial prostheses. Observe the abduction/adduction position. Align the prosthetic foot without external rotation.

**INFORMATION: To compensate for the deflection of the prosthetic sports foot, the sports prosthesis may have to be aligned longer than the everyday prosthesis.**

- ▶ Observe the instructions for use for the selected prosthetic components.

**The mounting instructions for the 4R420, 2R176=T, 2R177=5, 2R177=18 adapters are included with the products.**

### **Static alignment and dynamic fitting**

During the static alignment and the dynamic fitting, the settings of the prosthesis (length, angle, displacements) are checked and adapted on the user. Adjusting the prosthetic foot has a major influence on the behaviour of the prosthesis:

#### **Sagittal plane**

- **Moving in the anterior direction:** Increases the knee extending moment and makes the knee joint more secure (while standing and while walking). The prosthesis has diminished forward thrust. The user may feel like walking uphill.
- **Moving in the posterior direction:** Reduces the knee extending moment. The prosthesis has increased forward thrust. The user needs to use more muscle strength to keep the knee joint extended.
- **Adjusting the angle to be steeper:** Shortens the forefoot lever. This increases the forward thrust of the prosthesis and makes the prosthetic foot feel stiffer.
- **Adjusting the angle to be flatter:** Extends the forefoot lever. The prosthesis' forward thrust decreases and the prosthetic foot feels softer.

## Frontal plane

- Adjust the prosthesis so the user steps down flat and rolls over the virtual big toe.

## 5.5 For Runner only: Aligning the connection adapter

### Adjusting the angle

The connection adapter can be shifted on the prosthetic foot. It is shifted along a radius and only affects the angle at ground contact. The position in the a-p direction remains unchanged.

<b>Shift the prosthetic foot to anterior</b>	Adjusts the angle to be steeper
<b>Move the prosthetic foot in the posterior direction</b>	Adjusts the angle to be flatter

> **Required materials:** Loctite 241 636K13, torque wrench 710D20

- 1) Loosen the two screws on the underside of the adapter.
- 2) Adjust the angle of the prosthetic foot.
- 3) Secure the two screws with Loctite and tighten them (**12 Nm**).

### Turning

**Connection adapter with pyramid adapter only:** The pyramid adapter can be rotated to any position. The pyramid adapter is fixed when the set screws of the adapter that is attached proximally are tightened.

> **Required materials:** Loctite 241 636K13, torque wrench 710D20

- 1) Unscrew the two set screws that are screwed in the furthest.
- 2) **If the pyramid cannot be turned: Carefully tap the pyramid from the top to loosen it.**  
Set the desired rotation angle.
- 3) Secure the set screws with Loctite and tighten (for the required tightening torque, see instructions for use of the respective adapter).

## 6 Cleaning

- 1) Rinse the product with clear fresh water.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

## 7 Maintenance

- ▶ Schedule regular maintenance appointments with the user according to use.
- ▶ During regular inspection: Inspect the prosthesis for signs of wear and check its function.
- ▶ Conduct annual safety inspections.

## 8 Disposal

Do not dispose of the product in household waste. Improper disposal is harmful to health and the environment. Observe the instructions of the local authority regarding return, collection and disposal.

## 9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

## 10 Technical data

1E91					
Stiffness version	SPR-1	SPR-2	SPR-3	SPR-4	SPR-5
Build height [cm]	Approx. 30				
Build height, under load [cm]	Approx. 27				
Weight [g]	460	490	515	550	585
Max. body weight [kg]	50	60	72	86	104

1E93					
Stiffness version	SPR-1	SPR-2	SPR-3	SPR-4	SPR-5
Structural height [cm]	approx. 18.5				
Structural height, under load [cm]	approx. 16				
Weight [g]	170	180	190	200	210
Max. body weight [kg]	20	25	30	37	45

1E97					
Stiffness version	SPR-1	SPR-2	SPR-3	SPR-4	SPR-5
Min. build height [cm]	approx. 11				
Min. build height, under load [cm]	approx. 8.5				

1E97					
Stiffness version	SPR-1	SPR-2	SPR-3	SPR-4	SPR-5
Weight [g]	260	275	295	315	335
Max. body weight [kg]	20	25	30	37	45

Reference number	4R216=6	4R218=6	4R224=6
Weight [g]	175	180	175
System height [mm]	13	5	3
Max. body weight [kg]	125		45

Reference number	2Z540=6	2Z541=6	2Z543=6
Weight [g]	160	135	90
Structural height [mm]	15	12	13

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