dynion explore 3R84=P

DE	Gebrauchsanweisung (Fachpersonal)	10
EN	Instructions for use (qualified personnel)	24
FR	Instructions d'utilisation (personnel spécialisé)	38
IT	Istruzioni per l'uso (personale tecnico specializzato)	53
ES	Instrucciones de uso (personal técnico especializado)	68
PT	Instruções de utilização (pessoal técnico)	83
NL	Gebruiksaanwijzing (deskundig personeel)	98
SV	Bruksanvisning (fackpersonal)	112
DA	Brugsanvisning (faguddannet personale)	125
NO	Bruksanvisning (fagpersonell)	139
FI	Käyttöohje (ammattihenkilöstö)	152
PL	Instrukcja używania (personel fachowy)	166
HU	Használati útmutató (szakszemélyzet)	181
CS	Návod k použití (odborný personál)	195
RO	Instrucțiuni de utilizare (personal de specialitate)	208
HR	Upute za uporabu (stručno osoblje)	223
SL	Navodila za uporabo (strokovno osebje)	237
SK	Návod na použitie (odborný personál)	250
BG	Инструкция за употреба (за специалисти)	264
TR	Kullanım kılavuzu (uzman personel)	280
EL	Οδηγίες χρήσης (Τεχνικό προσωπικό)	294
RU	Руководство по применению (для специалистов)	310
JA	取扱説明書(有資格担当者)	326
ZH	使用说明书(专业人员)	339
KO	사용 설명서(전문가)	351

1

1 3R84=P

2 4G1127=1.9

3 2R37/2R38

4 2Z11=KIT

5 647G1866

6 647H1866

2

1 4G1127=1.9

2 4G1127=12

3 4G1127=5.9

4 4G1127=5.12

5 4G1127


6 4D34

3









4


743Y580=*
(3R80, 3R106=HD, =KD, =ST)


743A220
PRO.S.A. Assembly

8



9



10



3D-Mode



11


2D-Mode

20


1

Isopropyl alcohol
634A58

1 Min


2

Loctite 3090
636K44

3

5 Min

4

21


4 mm
15 Nm
Loctite 241 – 636K13

4 mm
10 Nm

4 mm
15 Nm
Loctite 241 – 636K13

506G3=M8X6

506G3=M8X10


506G3=M8X14

506G3=M8X16

506G3=M8X20


22

2


< 5 mm

3


≥ 5 mm

23


1

2

3

INFORMATION

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INFORMATION

Date of last update: 2026-04-23

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

The 3R84=P dynion explore prosthetic knee joint has the following key features:

- Monocentric prosthetic knee joint with rotation hydraulics
- Stance phase control
 - Adjustable stance phase flexion resistance (hydraulic damping)

- Adjustable switching threshold to deactivate stance phase flexion resistance at the end of the stance phase
- Lock (activated and deactivated by the user)
- Swing phase control
 - Adjustable swing phase flexion resistance (hydraulic damping)
 - Adjustable swing phase extension resistance (hydraulic damping)
 - Hydraulic unit with extension assist function (spring force)
- Adjustable modes (by the user)
 - Standard mode: Stance phase flexion resistance activated
 - Cycling mode: Stance phase flexion resistance deactivated

2 Safety-related information

CAUTION! Risk of injury and risk of product damage

- ▶ To prevent mechanical damage, treat the product with care.
- ▶ Observe the information on product combinations in the instructions for use for the products.
- ▶ Observe the maximum lifetime of the product (see page 27).
- ▶ Observe the product's operating conditions (see page 25).
- ▶ Do not expose the product to prohibited environmental conditions (see page 26).
- ▶ Check the product for damage and readiness for use prior to each use.
- ▶ Do not use the product if it is damaged or if you are not sure whether it is working properly. Have the product inspected, repaired or replaced by the manufacturer or a specialist workshop.
- ▶ Keep the product away from open flames, embers and other sources of heat.
- ▶ The product may become warm during use. Do not touch hot product components.
- ▶ If there are changes in functionality, reduce all activities so the product can cool down.
- ▶ If there are changes in functionality, have the product checked by the manufacturer or a specialist workshop.
- ▶ To avoid the risk of pinching, do not reach into the joint mechanism.
- ▶ In the event of oil leaks, stop using the product and have it repaired.
- ▶ The product is intended for use by one person only and may not be reused by other persons.
- ▶ Follow the maintenance instructions (see page 37).
- ▶ Observe the instructions for cleaning and care (see page 36).

3 Intended use

3.1 Indications for use

The product is to be used **solely** for lower limb exoprosthetic fittings.

3.2 Operating conditions

- Permissible body weight: **25 – 60 kg**

The product is recommended for:

- Mobility grade 3: unrestricted outdoor walker
- Mobility grade 4: unrestricted outdoor walker with particularly high demands

3.3 Target group

This document is intended for qualified personnel with manufacturing knowledge of lower limb prostheses.

3.4 Environmental conditions

Storage and transport

- Storage temperature: -20 °C to +60 °C
- Relative humidity: 20 % to 90 %
- No mechanical vibrations or impacts

INFORMATION

Clean the product after contact with humidity/chemicals/solids, in order to avoid increased wear and tear and damage (see page 36).

Allowable environmental conditions

- Temperature range: -10 °C to +45 °C
- Relative humidity: 20 % to 90 %
- Contact with soap solution (e.g. when taking a shower)
- Contact with perspiration
- Submersion in fresh water/chlorinated water (e.g. swimming pool):
 - Maximum immersion depth: 2 m
- Submersion in salt water with a low salt content (14 days/year, 1 hour/day):
 - Permissible salt content: max. 3.5 % (e.g. ocean)
 - Maximum immersion depth swimming pool: 1 m
 - Maximum immersion depth sea: 0.5 m
- Contact with dust, drifting sand (e.g. walking on the beach)
- Contact with salty, condensing air (e.g. in a 40 °C steam sauna)
- UV resistant

Unallowable environmental conditions

- Contact with liquid-binding particles (e.g. talcum)
- Contact with large amounts of sand and dust (e.g. digging, kneeling in sand, construction site work)
- Contact with acids
- Contact with urine
- Submersion in salt water with excessive salt content
- Salt content: > 3.5 % (e.g. brine bath, Dead Sea)
- Cleaning agents and disinfectants with solvents, chlorine and phosphate
- High water pressure (e.g. diving, jumping into water)

3.5 Indications

- Hip disarticulation, unilateral
- Transfemoral amputation, unilateral/bilateral
- Knee disarticulation, unilateral/bilateral

3.6 Product combinations

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

Permissible and impermissible product combinations with the 3R84=P prosthetic component are listed in the following tables.

Permissible foot combinations

1A30	1C51	1C61	F22
------	------	------	-----

Permissible foot combinations			
1C10	1C52	1C62	F23
1D35	1C53	1C63	F24
1C30-1	1C56	1C64	1K10
1C40	1C58	1C68	1K30
1E56	1C59	1C70	R11
1C50	1C60	F21	R16

Reference number	Designation	Product combination
7E7	Hip joint	✓
7E9	Hip joint	✓
7E10	Hip joint	✗

3.7 Lifetime

This prosthetic component has been load tested by the manufacturer according to ISO 10328. The maximum lifetime is 5 years.

4 Scope of delivery

3R84=P (see fig. 1)		
Item	Reference number	Designation
①	3R84=P	Prosthetic knee joint
②	4G1127=1.9	Front cover (grey – pre-assembled)
③	2R37/2R38	Tube adapter
④	2Z11=KIT	Adapter protector
⑤	647G1866	Instructions for use (qualified personnel)
⑥	647H1866	Instructions for use (user)

5 Accessories

INFORMATION

As a supplement to this document, a quickstart guide for qualified personnel is also available via the QR code and link provided.

Quickstart guide – 3R84=P

Adjustments and settings

<https://product-documents.ottobock.com/IFU/INT/3R84/647G1867/01/O/S/F>

3R84=P (see fig. 2)		
Item	Reference number	Designation
①	4G1127=1.9	Front cover (grey)
②	4G1127=12	Front cover (pink)
③	4G1127=5.9	Front cover (light blue)

3R84=P (see fig. 2)

Item	Reference number	Designation
④	4G1127=5.12	Front cover (azure blue)
⑤	4G1127	Front cover (dark grey)
⑥	4D34	Single component pack for foam cover (push buttons)

6 Preparing the product for use**⚠ CAUTION****Incorrect alignment, assembly or adjustment**

Injury due to incorrectly installed or adjusted as well as damaged prosthetic components

- Observe the alignment, assembly and adjustment instructions.

⚠ CAUTION**Initial use of the prosthesis by the user**

Falling due to lack of experience of the user, incorrect alignment or adjustment of the prosthesis

- For the safety of the user when standing and walking for the first time, use a suitable aid (e.g. parallel bars, handrail and anterior walker).

6.1 Optional: Using the adapter protector

The adapter protector protects the connecting section of the prosthetic joint from scratches during alignment and fitting.

- > The instructions refer to the following illustration: see fig. 3
- 1) Follow the instructions in the accompanying document for the adapter protector.
- 2) Remove the adapter protector from the product before handing it over to the user.
- 3) Readjust the prosthetic alignment after removing the adapter protector.

6.2 Optional: Using the 743A220 PROS.A.* Assembly alignment apparatus

An alignment apparatus such as the 743A220 PROS.A. Assembly can be used for the alignment of modular leg prostheses.

- > The instructions refer to the following illustration: see fig. 4
- 1) Follow the instructions in the accompanying document for the alignment device.
- 2) Mount the 743Y580=* adapter inserts with the labels, **3R80, 3R106=HD, =KD, =SD**.
- 3) Carry out the alignment of the prosthesis.

6.3 Shortening the tube adapter**⚠ CAUTION****Incorrect processing of the tube**

Risk of injury due to damage to the tube

- Do not clamp the tube in a vice.
- To shorten the tube, use only a tube cutter or a cutting device.

- > The instructions refer to the following illustration: see fig. 5
- > **Required materials and tools:** 704Y14=* cutting device or 719R5 tube cutter, 718R1 tube deburrer

- 1) Shorten the tube adapter with a cutting device or a tube cutter ①.
- 2) Deburr ③ the inside and outside of the cut edge with a tube deburrer ②.

6.4 Installing the tube adapter

- > The instructions refer to the following illustration: see fig. 6
- > **Required materials and tools:** 634A58 isopropyl alcohol, 710D20 torque wrench, hexagon socket screw adapter 4 mm
- 1) Clean the contact surfaces of the tube and of the tube receiver using a degreasing cleaning agent (such as isopropyl alcohol).
- 2) Insert the tube adapter into the tube receiver up to the stop **(0 mm) ①** and make a mark on the tube adapter (red line in the illustration) **②**.
- 3) **CAUTION! Pulling the tube adapter out further than permitted may result in the breakage of components.**
Pull the tube adapter out a **maximum of 13 mm** during adjustment. Use the marking on the tube adapter as a reference **②**.
- 4) Tighten the cap screw **without** additional Loctite to **10 Nm ③**.

6.5 Bench Alignment

- Align the prosthesis using the alignment images (see fig. 8; see fig. 9) as well as the instructions for use of all prosthetic components used.
- The prosthesis should be aligned so as to enable the user to stand safely.

NOTICE! Design the prosthetic socket in such a way that no metal parts press against the posterior area of the prosthetic knee joint at maximum flexion. This prevents damage to the product.

- > **Required materials and tools:** alignment apparatus (e.g. 743A200 PROS.A. Assembly), 743A80 50:50 gauge, 743S12 heel height measuring device
- 1) **Lock**
 - Perform the bench alignment of the prosthesis with the prosthetic knee joint locked.
 - During activation and deactivation of the lock, push the prosthetic knee joint against the extension stop.
- 2) Position the prosthetic foot in the alignment apparatus.
- 3) **Heel height**
 - The effective heel height of everyday shoe (x) **(+ 5 mm)**.
- 4) **Frontal plane (see fig. 8)**
 - Align the prosthetic foot with the alignment reference line: **Position the middle of the prosthetic foot on the alignment reference line (external rotation approx. 5 – 7°).**
- 5) **Sagittal plane (see fig. 9)**
 - Align the prosthetic foot with the alignment reference line: **Place the alignment reference line behind the middle of the prosthetic foot (+ 30 mm).**
- 6) Position the prosthetic knee joint and prosthetic socket in the alignment apparatus.
- 7) **Frontal plane (see fig. 8)**
 - Align the prosthetic socket with the alignment reference line.
 - Take the individual adduction position into account.
- 8) **Sagittal plane (see fig. 9)**
 - Align the prosthetic socket with the alignment reference line.
 - Determine the centre of the prosthetic socket proximally and distally with the 50:50 gauge in the direction of movement and mark the centre line.
 - Mark the socket reference point on the centre line proximal to the tuberosity reference point **(30 mm)**.
 - Align the prosthetic socket so the alignment reference line runs through the socket reference point.
 - Adjust the socket flexion by rotating around the socket reference point: **individual residual limb flexion/flexion contracture (x) (+ 5°)**

6.6 Static Alignment

Ottobock recommends optimising the alignment with the 743L500 3D L.A.S.A.R. Posture measuring device. The 743L100 L.A.S.A.R. Posture can be used as an option.

Optimising the prosthetic alignment:

Make sure that the user is standing comfortably and that the reference points in the illustrations are observed.

Noting the lines:

Dotted lines = auxiliary lines

Solid lines = load lines

Maintaining distances:

Use the distances between the auxiliary line and load line shown in the illustrations for static prosthetic alignment.

3D mode (see fig. 10)

- ❶ **Prerequisite:** The prosthetic knee joint is fully extended. The lock is activated.
- ❷ The auxiliary line lies on the sagittal reference point of the prosthetic knee joint (pivot point).
- ❸ The auxiliary line lies on the frontal reference point (anterior superior iliac spine).
- ❹ The auxiliary line lies on the frontal reference point of the prosthetic knee joint (pyramid).

2D mode (see fig. 11)

- ❶ **Prerequisite:** The prosthetic knee joint is fully extended. The lock is activated.
- ❷ The auxiliary line lies on the sagittal reference point of the prosthetic knee joint (pivot point).
- ❸ The auxiliary line lies on the frontal reference point (anterior superior iliac spine).
- ❹ The auxiliary line lies on the frontal reference point of the prosthetic knee joint (pyramid).
- ❺ The auxiliary line lies on the frontal reference point of the prosthetic foot (pyramid).

6.7 Dynamic Trial Fitting

CAUTION

Adjusting the settings

Falling due to improper or unfamiliar settings

- ▶ Only adapt the settings to the user gradually.
- ▶ Explain to the user the effects of the adjustments on the use of the prosthesis.

The optimal gait pattern is developed during dynamic fitting. The prosthesis is aligned optimally in the frontal plane and the sagittal plane for this purpose.

6.7.1 Front cover

The front cover has to be removed to adjust the prosthetic knee joint or to replace it.

> The instructions refer to the following illustration: see fig. 7

- 1) Open the cover at the lower notch using a blunt object ❶.
Alternatively, the cover can be lifted at the side edge using a fingernail ❷.
- 2) Pull the opened cover off the prosthetic knee joint ❸.
→ The cover has been removed and can be replaced.
→ The prosthetic knee joint can be adjusted.
- 3) After adjusting the prosthetic knee joint or replacing the cover, attach it. Make sure that the studs are properly fitted into the stud openings ❹.
→ The cover is attached correctly.

6.7.2 Checking the settings

- 1) Manually check the settings of the product prior to the dynamic fitting.
- 2) In case of deviations, reset the product to the factory settings.

6.7.3 Factory settings

NOTICE! Do not turn the switching threshold past its stops, or you will damage the thread on the switching threshold.

3R84=P factory settings (see Fig. 12–15)			
Function	Adjustment position	Factory settings	Meaning
Lock	Lock symbols on the sides of the joint (see fig. 16)	Open lock symbol has been pressed	The lock is deactivated
Cycling mode	Bicycle symbol and walking symbol on the sides of the joint (see fig. 17)	Walking symbol has been pressed	The cycling mode is deactivated
Switching threshold	Lowest adjustment screw (see fig. 15)	3 o'clock (from left-hand stop, 1/2 turn clockwise)	Increased force required to transition from the stance phase to the swing phase
Stance phase flexion resistance	Adjustment screw S (see fig. 12)	Middle setting	Medium resistance
Swing phase flexion resistance	Adjustment screw F (see fig. 13)	Middle setting	Medium resistance
Swing phase extension resistance	Adjustment screw E (see fig. 14)	Middle setting	Medium resistance

6.7.4 Adjusting and setting

INFORMATION

As a supplement to this section, the following videos are available for qualified personnel via the QR codes and links provided.

Dynion and Dynion Explore – adjustments and settings

Adjustments and settings – available languages: English

<https://youtu.be/ukZ1Q-dgm5A>

Dynion and Dynion Explore – User training

User training – available languages: English

<https://youtu.be/zMZZBA0-h0>

The recommended sequence for adapting and adjusting the prosthetic knee joint is as follows:

1. Sitting down, 2. Walking, 3. Walking on ramps and stairs, 4. Sitting down, 5. Walking

Sitting down

- > The settings for sitting down should provide the user with adequate safety while not generating excessive resistance. To do so, the stance phase flexion resistance (see page 33) must be adjusted properly.

- 1) Adjust the stance phase flexion resistance.
- 2) Have the user practise sitting down safely several times and adjust the settings as necessary.

- 3) Check the stance phase flexion resistance again after making the adjustments for walking on ramps and stairs.

Walking

CAUTION

Risk of injury due to incorrect adjustment of switching threshold

Falling due to flexing the prosthetic knee joint too early.

- Allow the user to walk safely only with safeguards in place once the switch threshold has been adjusted correctly.

INFORMATION

When setting the switching threshold, it must first be turned anticlockwise to the stop. Then turn the switching threshold clockwise two full turns.

- > Adjust the settings first at normal walking speed, then with short, fast steps and finally with long, fast steps. Check the gait pattern and the effects after each adjustment. Repeat the walking exercises several times.
- > The procedure described ensures that the switching threshold is not enabled too early or not at all. The user gains confidence in how it works.
- 1) At the start of the walking exercises: Set the switching threshold (see page 33) high enough that the swing phase is **not** enabled. This requires several turns.
- 2) Decrease the switching threshold in small increments (**max. 30°** per step) until the swing phase is initiated.
 - Check the gait pattern after each change.
- 3) Adjust the swing phase flexion resistance (see page 33) so the lower leg of the prosthesis does not swing too far in dorsal position and reaches full extension in time for the next heel strike.
- 4) Adjust the swing phase extension resistance (see page 34) so the prosthetic knee joint does not swing too hard and reaches full extension in time for the next heel strike.
 - Check the switching threshold after each change. If the stance phase flexion resistance (see page 33) is deactivated under heel load, the switching threshold must be increased.
 - After making adjustments (especially for settings in the range marked with!), test the settings while walking with support (e.g. between parallel bars) at different walking speeds.
 - The switching threshold buffer exhibits settlement (especially during the first 10 minutes). Therefore, while walking on a level surface, check whether short, quick steps are still possible or whether the settings need to be readjusted.

Switching between the stance phase and swing phase

CAUTION! Fast, hard movements may deactivate stance phase flexion resistance. Therefore, follow the instructions below and avoid the movements mentioned.

- > The instructions refer to the following illustration: see fig. 23
- 1) Avoid letting the forefoot strike the ground hard. This is particularly critical in combination with strong extension (e.g. at a kerb) followed immediately by knee flexion ❶.
- 2) Avoid quick, forceful steps and forward jumps, as well as strong hip extension moment at heel strike followed immediately by knee flexion ❷.
- 3) Avoid walking backward with a load on the prosthetic forefoot, with simultaneous knee flexion ❸.

Walking on ramps and stairs

- > The settings for walking on ramps and stairs should provide the user with adequate safety while not generating excessive resistance. To do so, the stance phase flexion resistance (see page 33) and swing phase extension resistance (see page 34) must be adjusted properly.
- 1) Adjust the stance phase flexion resistance.

- 2) Set the swing phase extension resistance so that full extension is reached at heel strike if possible.
- 3) Have the user practise walking on ramps and stairs several times and at different walking speeds, adjusting the settings as necessary.
 - Advise the user that when going down stairs, they should place the prosthetic foot no further forward than the middle of foot.

Using the lock

- 1) Explain the function of the lock and the position of the switch to the user (see page 35).
- 2) Practise using the lock with the user.

Using cycling mode

- 1) Explain the cycling mode and the position of the switch to the user (see page 36).
- 2) Practise using cycling mode with the user.

6.7.4.1 Stance phase flexion resistance

INFORMATION

Noises during flexion in the stance phase

There may be noise during flexion in the stance phase up to 5° and subsequent extension. This is normal and harmless.

The stance phase flexion resistance is activated at heel strike. It prevents the prosthetic knee joint from flexing too early during the stance phase.

- > The instructions refer to the following illustration: see fig. 12
- > **Required materials and tools:** hexagon socket screw 2 mm
- ▶ Adjust the stance phase flexion resistance so it provides the user with adequate safety while not generating excessive resistance.
 - **For more resistance**, turn adjustment screw **S** clockwise.
 - **For less resistance**, turn adjustment screw **S** anticlockwise.

6.7.4.2 Switching threshold

The switching threshold deactivates the stance phase flexion resistance at the end of the stance phase.

- > The instructions refer to the following illustration: see fig. 15
- > **Required materials and tools:** hexagon socket screw 4 mm
- 1) Adjust the switching threshold so that it provides the user with adequate safety while not generating excessive resistance.
 - Turn the adjustment screw ❶ clockwise to **secure stance phase damping for a longer period**.
 - Turn the adjustment screw ❶ anticlockwise to **secure stance phase damping for a shorter period**.
- 2) The spring cartridge ❷ cannot be adjusted!

6.7.4.3 Swing phase flexion resistance

- > The instructions refer to the following illustration: see fig. 13
- > **Required materials and tools:** hexagon socket screw 2 mm
- ▶ Adjust the swing phase flexion resistance so that the prosthetic foot swings through sufficiently but not too far.
 - **For more resistance**, turn adjustment screw **F** clockwise.
 - **For less resistance**, turn adjustment screw **F** anticlockwise.

6.7.4.4 Swing phase extension resistance

CAUTION! Set the swing phase extension resistance so that the prosthetic knee joint reaches full extension, even when taking slow and small steps.

- > The instructions refer to the following illustration: see fig. 14
- > Take extra care when using the adjustment range marked with an !. Adjustments in this range may cause stance phase flexion resistance to be deactivated under heel load.
- > **Required materials and tools:** hexagon socket screw 2 mm
- ▶ Adjust the swing phase extension resistance so the prosthetic knee joint reaches full extension but does not swing too hard against the extension stop.
 - **For more resistance**, turn adjustment screw **E** clockwise.
 - **For less resistance**, turn adjustment screw **E** anticlockwise.

6.8 Optional: Installing the foam cover

CAUTION

Use of liquid-binding substances (e.g. talcum)

Risk of injury, damage to the product due to lack of lubrication

- ▶ Prevent contact of the product with liquid-binding substances.

- > The instructions refer to the following illustration: see fig. 20
- > **Prerequisite:** Use the foam cover (**3S107**) for the product.
- > **Required materials and tools:** foam cover (3S107), silicone spray (519L5), isopropyl alcohol (634A58), Loctite 3090 (636K44), microfibre cloth
- > **Optional:** push buttons for foam cover (4D34)
- 1) Clean the push buttons on the back and the prosthetic knee joint at the switches with isopropyl alcohol and a microfibre cloth **①**. Wait **1 minute**.
- 2) **NOTICE! Do not use too much adhesive so it does not run into the slits of the product.** Apply Loctite to the back of the push buttons only at specific points **②**.
- 3) Attach the push buttons to the product **③** and allow the adhesive to cure for **5 minutes** **④**.
- 4) Spray the silicone spray directly onto the contact surfaces in the foam cover to improve the gliding characteristics and reduce noise.
- 5) Attach the foam cover to the prosthesis.
 - Ensure that the product can be operated safely. One possible measure is to make holes in the foam cover on the controls.
 - After attaching the foam cover, perform another **dynamic fitting** (see page 30) to check the settings.
- 6) Finish the prosthesis.

6.9 Finishing the prosthesis

Finish working on the prosthesis using the reference image (see fig. 21) and the following instructions.

- 1) Remove the adapter protector from the product before handing it over to the user.
- 2) Readjust the prosthetic alignment after removing the adapter protector.
- 3) Reattach the cover on the prosthetic knee joint.
- 4) Replace set screws that are too short with longer ones to prevent components from coming loose accidentally.
- 5) Replace set screws that are too long with shorter ones to prevent collisions and product damage.
- 6) Clean the threads of the screw connections on the prosthetic socket and foot using a degreasing cleaner.
- 7) Only use Loctite 241 (636K13) on the screw connections of the prosthetic socket and foot.
- 8) Finish all screw connections. Observe the assembly tightening torques and thread locks for all prosthetic components.

→ The prosthesis is completed.

7 Use

INFORMATION

If the product is stored at a low temperature (-10 °C to +10 °C), this can lead to changes in functionality due to the cold hydraulic oil. The hydraulic oil heats up again during use of the product.

► Flex and extend the joint several times to warm the oil.

7.1 Switching between the stance phase and swing phase

- The prosthetic knee joint is pushed into hyperextension when the forefoot is loaded at the end of the stance phase. If you then initiate a knee flexion moment, stance phase flexion resistance is deactivated and swing phase flexion resistance is activated.
- Stance phase flexion resistance is activated when switching from knee flexion to knee extension in the swing phase.
- The prosthetic knee joint is in extension at heel strike. Stance phase flexion resistance is activated, securing the knee.

7.2 Using the lock

⚠ CAUTION

Locking function

Risk of injury due to unintentional activation or deactivation of the locking function

- Check that the joint is functioning correctly each time you activate or deactivate the lock.
- Make sure that the foam cover and clothing cannot activate the lock.

INFORMATION

Noises may occur due to the play between the extension stop and lock pin when walking with the lock activated. The lock can only be released by actively pressing the deactivation button.

Activating the lock

- > Activate the lock for use in wet areas (e.g. showering and swimming) and relaxed standing.
- > The instructions refer to the following illustration: see fig. 16
- 1) Put weight on the prosthetic knee joint and extend it to the stop.
- 2) Press the switch to activate the lock ❶.
 - **TIP:** Press the switch in the marked area ❷; this requires the least amount of operational force.
- The lock is activated.

Deactivating the lock

- > Deactivate the lock for normal walking.
- > The instructions refer to the following illustration: see fig. 16
- 1) Put weight on the prosthetic knee joint and extend it to the stop.
- 2) Press the switch to deactivate the lock ❸.
 - **TIP:** Press the switch in the marked area ❹; this requires the least amount of operational force.
- The lock is deactivated.

7.3 Using cycling mode

CAUTION

Walking in cycling mode

Risk of injury due to lack of stance phase flexion resistance in cycling mode

- ▶ Switch back to standard mode each time after using cycling mode.
- ▶ Carefully verify that standard mode has been activated after switching.

Activating cycling mode

- > Activate cycling mode for cycling and similar movement patterns.
- > The instructions refer to the following illustration: see fig. 17
- ▶ **INFORMATION: The standard mode and the stance phase flexion resistance are deactivated by activating the cycling mode.**
Activate cycling mode by pressing the switch all the way in ①.
 - **TIP:** Press the switch in the marked area ②; this requires the least amount of operational force.
 - When the switch for cycling mode is pressed, the switch for normal walking protrudes slightly.
- Cycling mode is activated.

Deactivating cycling mode

- > Deactivate cycling mode for normal walking.
- > The instructions refer to the following illustration: see fig. 17
- ▶ **INFORMATION: Activating standard mode deactivates cycling mode and activates the stance phase extension resistance.**
Deactivate cycling mode by pressing the switch all the way in ③.
 - **TIP:** Press the switch in the marked area ④; this requires the least amount of operational force.
 - When the switch for normal walking is pressed, the switch for cycling mode protrudes slightly.
- The cycling mode is deactivated.

7.4 Use in water

Prior to use in water

NOTICE! Activate the lock and cycling mode before using the product in water to avoid product damage.

- > The instructions refer to the following illustration: see fig. 18
- 1) Activate the lock ①.
- 2) Activate cycling mode as well ②.
- The product can be used in water ③.

After use in water

- > The instructions refer to the following illustration: see fig. 19
- 1) Deactivate cycling mode ①.
- 2) Then deactivate the lock ②.
- The product can be used normally again ③.

7.5 Cleaning

Requirements for cleaning agents and disinfectants

- Solvent-free
- Chlorine-free
- Phosphate-free

Cleaned using soap and water

- > **Required materials and tools:** 2 microfibre cloths, pH-neutral soap (e.g. 453H10 Derma Clean), tap water
- 1) **In case of light soiling (e.g. small deposits):** Moisten a microfibre cloth with water.
- 2) **In case of heavy soiling (e.g. stubborn deposits):** Moisten a microfibre cloth with pH-neutral soap and water.
- 3) Clean the product with the moistened cloth.
- 4) Dry the product with a dry microfibre cloth.
- 5) Allow the product to air dry.

Cleaning after contact with salt water

NOTICE! Remove salt crystals from the prosthetic knee joint to prevent corrosion.

- > **Required materials and tools:** 2 microfibre cloths, warm tap water
- 1) Activate cycling mode.
- 2) Fully immerse the joint in warm tap water and flex and extend it several times.
 - If crystals stick, carefully wipe them off with a microfibre cloth. Immerse the joint in water again and rinse it off.
- 3) Remove the joint from the water and allow it to drain.
- 4) Dry it thoroughly with a microfibre cloth.
- 5) Check the functionality of the joint (flexing, extending, locking).
- 6) Deactivate cycling mode.
- After maximum use in salt water (**1 hour/day for a maximum of 14 days/year**), have it checked by qualified personnel.

Disinfecting the product

- > **Required materials and tools:** 2 microfibre cloths, water-based disinfectant
- 1) Apply disinfectant to the microfibre cloth.
- 2) Disinfect the product with the moistened cloth.
- 3) Dry the product with a dry microfibre cloth.
- 4) Allow the product to air dry.

8 Maintenance

- ▶ **Every 12 months:** Have maintenance performed by qualified personnel. Reduce the maintenance interval based on the use of the prosthesis.
 - In case of insufficient functionality, adjust the prosthesis.
- ▶ Check the prosthesis for **wear and tear, functionality and unusual noises**.
- ▶ Do not lubricate or grease the prosthetic joint.
- Repair work must be performed exclusively by manufacturer service.
- ▶ If the product is used in salt water: check it closely for salt residues and signs of corrosion (e.g. set screws).
- ▶ Check all product functions for safety. In particular, check the locking function, activated and deactivated cycling mode, flexion and extension and the switching threshold.
- ▶ In case of functional limitations and damage, send the product to the manufacturer's service department.

8.1 Checking the sacrificial anode

The sacrificial anode ❶ protects the product against electrochemical corrosion by corroding itself first, thereby protecting the product from damage (see fig. 22).

- ▶ Have the prosthesis checked as soon as possible if corrosion is visible. Example: In case of rust on the set screws that attach the proximal adapter to the prosthetic knee joint.
- ▶ Check the sacrificial anode **every 12 months**:
 - 1) Remove dirt and crystals from the area around the sacrificial anode.
 - 2) Check the wear and tear of the sacrificial anode using callipers:
 - Result < **5 mm** ❷: The product can continue to be used.

- Result $\geq 5 \text{ mm}$ ③: The sacrificial anode is too worn. Have a new sacrificial anode installed by the manufacturer's service department.

9 Disposal

Do not dispose of the product in household waste. Improper disposal is harmful to health and the environment. Observe the instructions of the local authority regarding return, collection and disposal.

10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

The full text of the EU Declaration of Conformity is available at the following Internet address: <http://www.ottobock.com/conformity>

11 Technical data

3R84=P	
Proximal connection	Pyramid
Distal connection [mm]	30
Knee flexion angle [°]	145
Total weight [g]	1217
Overall height [mm]	231
Build height [mm]	216.5
Mobility grade	3, 4
Permissible body weight [kg]	25 — 60

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This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

A series of horizontal lines for writing, consisting of 20 lines.

A series of horizontal lines for writing, spanning the width of the page. The lines are evenly spaced and extend from the left margin to the right margin.

Ottobock SE & Co. KGaA
Max-Näder-Straße 15 · 37115 Duderstadt · Germany
T +49 5527 848-0 · F +49 5527 848-3360
healthcare@ottobock.de · www.ottobock.com