
Kneuro Microprocessor Knee

Instructions For Use

V2.3.3

REF BR4B REF BR4C

www.brainrobotics.com/resource/downloads/ifu

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EN

These Instructions for Use are intended for the patient (user) and their certified prosthetist.

All images, screenshots, and figures are for reference only and are subject to change.

Important Information for the Certified Prosthetist

To ensure the safe and effective use of the Kneuro Microprocessor Knee System and to avoid injury, read the Instructions for Use and installation guide in their entirety before installing or operating this device. This device is to be used under the guidance of a certified prosthetist only. The patient should be provided with these Instructions for Use.

This device is required to be installed by an authorized prosthetist who has completed training and is certified by BrainRobotics. If you need help with product setup, use, or maintenance, contact the BrainRobotics support team at +1 (512) 969-2989 or support@brainrobotics.com.

Federal law (USA) restricts this device to sale by or on the order of a practitioner licensed by the laws of the state in which they practice to use or prescribe the use of the device.

Patient Information

This product can only be installed by an authorized prosthetist who has completed training and is certified by BrainRobotics. If you need help with product setup, use, or maintenance, contact your prosthetist or the BrainRobotics support team at +1 (512) 969-2989 or support@brainrobotics.com.

To ensure the device is used safely and to avoid injury, read the information for use and follow the instruction provided by your authorized prosthetist when using the Kneuro Microprocessor Knee System. This device is to be used under the guidance of a certified prosthetist only.

Any serious incident in relation to device must be reported to the manufacturer and relevant authorities.

Liability Disclaimer:

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorized modification of the device.

Warranty:

Purchase of this device includes a 36-month warranty covering all manufacturer defects. The warranty is only valid

if usage of the device is compliant with the manufacturer's recommendations detailed in these Instructions for Use. Note: Failure to comply with the Instructions for Use will result in voiding the warranty. An extended warranty and service recommendation is available.

For service or repair, contact BrainRobotics customer support at +1 (512) 969-2989 or support@brainrobotics.com.

Device Description


The Kneuro Microprocessor Knee, components, and Kneuro App make up the Kneuro Microprocessor Knee System as referred to in this Instructions for Use. The Kneuro Microprocessor Knee is a lower limb prosthetic knee designed for integration as part of an endoskeletal prosthesis system. It features an advanced sensor system which works in conjunction with a real-time control algorithm to adjust the hydraulic system, providing resistance tailored for a natural, adaptive user experience. The device is designed for everyday use. The Kneuro App is used to configure the device, adjust the settings and switch the operational modes of the device. The term 'device' is used to describe the Kneuro microprocessor-controlled knee in this document. The device is universally designed to be connected to other manufacturers components as part of the endoskeletal prosthesis system.

Each user must be evaluated by a certified prosthetist to determine suitability to operate and use device. Upon delivery, the device must be fit, aligned, and programmed by a certified prosthetist that has been authorized by BrainRobotics. The device is intended for a single user . The device has an IP67 rating.

Key components of the device include:

Figure 1

#	Part
1	Proximal pyramid adapter
2	Control button
3	Flexion stops / location of device label
4	Hydraulic actuator
5	USB-C charging port
6	Status indicator
7	Distal tube clamp adapter, 30mm

Device is shipped with the following:


Figure 2

Qty	Item
1	Kneuro knee (BR4B/BR4C)
1	Pylon, 30mm
2	10° Flexion stops
2	15° Flexion stops
2	20° Flexion stops (pre-installed)
2	25° Flexion stops
8	Flexion stop screws
1	2mm Allen key
1	1 Charging adapter
1	USB-C charging cable
1	Flexion installation guide
1	Instructions for Use
1	Warranty card
1	Certificate of quality
1	Kneuro app (via QR scannable code included in Instructions for Use)

Intended Purpose

The Kneuro Microprocessor Knee System is intended for use as a part of a lower limb endoskeleton prosthesis intended to actively assist individuals with above-knee, hip- or knee-disarticulation amputations in performing daily activities such as standing, walking, ascending and descending stairs and slopes, standing and sitting activities, and kneeling.

The Kneuro app is intended to be used with the Kneuro Microprocessor Knee via Bluetooth connection to set and manage users' access, customize parameters, and modes. Professional users are required to complete the Kneuro Certification course prior to use of the application with the device. The Kneuro app must be used to configure the Kneuro Microprocessor Knee prior to use.

The intended clinical benefit when used following required training and in accordance with Instructions for Use is to supplement the knee functionality for amputees.

Intended Users

The device is intended to be installed, calibrated and adjusted by a certified prosthetist trained by BrainRobotics.

The device is intended to be used and operated by competent individuals, trained by a certified prosthetist on safe device operation and use.

Intended Use Environment

Water – The device has an IP67 rating, indicating it is dust-tight and waterproof up to 1 meter for 30 minutes.

Do not submerge the device in saltwater or chlorinated water.

Outdoors – The device is Weatherproof. A Weatherproof device can be used in a wet or humid environment and can tolerate being submerged in fresh water for 30 minutes at a time, in a maximum depth of 1 meter.

Do not use in temperatures below 14°F (-10°C).

Indoor/Outdoor Activities – The device is intended to be used during a variety of activities in various terrains and settings.

Indications

The Kneuro Microprocessor Knee System is indicated for the following:

- Patients with mobility levels K2 - K4
- Unilateral and bilateral transfemoral amputations
- Unilateral and bilateral knee disarticulation amputations
- Unilateral hip disarticulations
- Patients weighing under 330lbs (150kg)
- Patients with the cognitive ability to follow instructions, operate and care for device

Contraindications

The Kneuro Microprocessor Knee System is contraindicated for the following:

- Users weighing more than 330lbs (150kg)
- Users lacking the cognitive ability to operate the device

Performance Characteristics

- **Adaptive Gait Control:** The Kneuro Microprocessor Knee supports adaptive control for both stance and swing phases, automatically adjusting to various walking conditions such as road surface, incline, stair navigation, and walking speed. This ensures smooth transitions and adaptability for diverse environments.
- **Customization and User-Specific Adjustments:** Using dedicated software, the control parameters of the Kneuro Microprocessor Knee can be customized to meet individual user requirements. Custom modes tailored for special activities—such as cycling, golfing, and others—are preset by a certified prosthetist via the Kneuro app and activated by the user as needed.
- **Water Ingress Protection:** The Kneuro Microprocessor Knee is IP67-rated, offering high protection against dust and temporary water submersion.
- **Power Adapter:** IP22, providing protection against limited water ingress when tilted up to 15 degrees.
- **Battery Life:** Under normal usage, battery life can last up to five days. This estimation is based on an average of 5,000 steps per day at room temperature, with idle periods factored in.
- **Safe Mode Functionality:** When the battery is depleted, the Kneuro Microprocessor Knee enters Safe Mode, allowing users to walk short distances. While adaptive control features are disabled, the knee remains secure and will not release into swing phase. However, users can still flex the knee to sit down or descend ramps and stairs comfortably. It is recommended to recharge the product promptly to restore full functionality.
- **Real-Time Adaptive Control:** The built-in sensor system continuously measures environmental and movement data, enabling real-time adjustments to the hydraulic system based on user activity and terrain. This adaptive feedback system supports a comfortable, responsive experience for daily mobility and various activities.

Device Lifetime

The expected lifetime of the device is 6 years. The device's expected lifetime may vary depending on the user's activity level, usage environment, and maintenance practices.

Do not use the device for activities involving unstable surfaces or positions requiring exceptional balance. This may result in insufficient support, risking falls or injury. Consult your certified prosthetist before attempting any such activities.

Warnings and Precautions

Warnings and precautions specific to the device should be understood by the certified prosthetist and communicated to and understood by the user prior to use.

Possible noise during movement: When using Kneuro under regular conditions, minor operational noise may arise from internal components. If you observe a significant increase in sound level during the knee's lifecycle, immediately consult your prosthetist or BrainRobotics customer support.

General warnings include:

- **Do not** ignore the safety instructions in this manual. Using this device beyond its designated scope can cause the device to malfunction and/or risk injury.
- **Do not** use the device when operating a motor vehicle.
- Immediately stop device use if there is a change or loss in device functionality or if there are signs of damage or wear affecting the structural integrity of the device. The user must contact their certified prosthetist for further instruction.
- **Do not** enter areas with strong magnetic fields or strong electrical interference (i.e., theft-prevention security systems, metal detectors, etc.) or rooms with strong magnetic fields (i.e., MRI, MRT equipment, etc.). These systems can cause irreparable damage to the device. If unavoidable, remove the device and leave it outside of the area.
- **Do not** use damaged power plugs, adapters, or charging cables. Exposed electrified components may cause electric shocks, and damaged electrical equipment can cause fires.
- **Do not** use the device while charging, as this can lead to device damage and/or personal injury.
- **Do not** submerge the device in water while charging, as this can lead to device damage and/or personal injury.
- This device contains a lithium-ion battery that is not field serviceable and requires all repairs to be done by the manufacturer. Lithium-ion batteries contain hazardous metals and should never be disposed of in residential or commercial garbage, nor should they be incinerated due to explosive materials. Refer to local regulations for proper disposal instructions.
- The device contains small parts and cords that can be choking hazards.
- **Do not** disassemble, repair, or modify this device or any related components, as it may cause malfunctions or personal injury.

General precautions include:

- Alignment, programming, and initial operation must be completed by a certified prosthetist that has been authorized by BrainRobotics.
- **Do not** use the Kneuro Microprocessor Knee when performing high-risk activities that may lead to injury, such as extreme sports, climbing, or other activities beyond typical daily use. The knee is designed to

provide stability and support for everyday mobility and walking. When engaging in activities where swing initiation is not desired, it is recommended to activate the Manual Lock.

Prosthetic Installation Instructions

Connect to the Kneuro App

Turn on the device.

Download and open the Kneuro app on an iOS or Android device.

☒ System Requirement: iOS 12 or above/Android 5.0 or above.

Clinicians must complete the Kneuro certification course before registering for the Kneuro app.

Select “Certified Prosthetist” on the welcome screen and enter login credentials OR create an account.

Clinicians will need to use their phone number or email address to register their account.

☒ There are two types of accounts: a **clinician account** and a **user account**. Assignment of account types is based on login credentials and account authorization.

In the clinician account, the certified prosthetist can access user management, parameter customization, and mode customization. The user list page is where the clinician can add new user profiles, select an existing user profile, or connect to a trial unit.

To pair the device to a user profile:

1. For **new profiles**, select the ‘+’ icon in the top right corner and select “Generate new profile”.
For **existing profiles**, select the correct user profile from the user list, OR select the ‘+’ icon in the top right corner and select “Associate existing profile” to add the profile to the user list.
For **trials**, select the ‘+’ icon in the top right corner and select “Proceed without profile”.
2. Confirm the serial number (SN) on the back of the device.
3. Select “Pair a Kneuro” to link the user’s profile for future connections.

☒ The device and your Bluetooth must be on prior to pairing.

To create a user account, the prosthetist must first create a user profile. The patient then registers using an email address or phone number, completes setup, and pairs the knee. Once connected, the user can adjust resistance parameters, switch modes, check battery status, and monitor step count.

A patient trialing a knee does not require a user account. Select “Trial Mode” from the user login screen instead.

For assistance logging in, contact BrainRobotics at +1 (512) 969-2989 or support@brainrobotics.com.

Bench Alignment - Alignment Goals

The alignment reference line should:

- Cross through the lateral bisection of the socket in estimation of the user's weight line.
- Be 0-10° anterior to the knee rotation axis.
- Cross through the posterior third of the foot.
- (5mm heel wedge under the foot without a shoe.)

Establish an alignment reference line perpendicular to the ground:

1. Place the posterior third of the foot on the alignment reference line.
(Cross reference the manufacturer's prosthetic foot recommended alignment.)
2. Place a 5mm heel wedge under the foot (without a shoe) and set the transverse plane rotation.
3. Align the knee center 0-10mm anterior to the alignment reference line. The device joint center should be in line with the contralateral joint center.

Follow alignment procedures carefully to avoid premature triggering of stance release (especially on ramp descent). Failure to do so can result in personal injury.

4. Mark a bisection line through the lateral aspect of the socket.
5. Position the socket so that the alignment reference line passes through the marked bisection in estimation of the user's weight line.
6. Conduct a Thomas Test on the user. Place the socket in 5° greater flexion than recorded hip flexion contracture.


Figure 3

7. Place the socket in 5° of adduction.
8. Set the appropriate prosthetic height.
9. Ensure at least 3mm of clearance is met when the knee is in maximum flexion. If there is not appropriate clearance, use a 2mm hex wrench to replace the pre-installed flexion stops with others included in the scope of delivery.

☒ Pylon insertion depth is 40mm, torque specifications at 10Nm.

Flexion Stop

The device comes pre-installed with the 20° flexion stops, which limits the maximum flexion angle to 125°.

The 10°, 15°, and 25° flexion stops can be installed to increase or decrease the flexion angle.

The flexion stops can be replaced with a 2mm hex wrench.

The device is not permitted for use without installing flexion stops.

☒ Perform a maximum flexion check, as the device and/or hydraulic unit may be damaged if the socket contacts the device.

☒ Ensure that there is a minimum distance of 3mm between the socket and the device.


Figure 4

If needed, install different sized flexion stops (included). [See figure 4]

See Flexion Stops Installation Guide [See figure 5]

Static Alignment

1. Power ON the device.
2. Have the user evenly load the prosthesis. Confirm the appropriate height of the prosthesis.
3. Confirm appropriate transverse plane alignment.
4. Confirm that the device is not yielding during quiet standing when weight is evenly distributed on both legs.


Figure 5

Dynamic Alignment

1. With the device powered ON, connect to the device using the Kneuro app. This is where all device parameters will be adjusted.
2. Have the user test the yielding resistance level by sitting into a chair.
3. In parallel bars, have the user take their first step with their sound limb. This will allow the user to experience stance release with the appropriate triggers on the device.
4. Evaluate both stance and swing phase of gait, ensure step-length symmetry, appropriate heel rise, and yielding.
5. Once the user is confident inside the parallel bars, move to level ground outside of the bars.
6. Instruct users on how to descend ramps and stairs, ascend stairs step-over-step, and perform walk-to-run, if appropriate.
7. The certified prosthetist must provide the user instructions for:
 - Powering ON/OFF
 - Activating and proper utilization of modes
 - Engaging manual lock
 - Charging the device
 - Feedback signals – LED, audio, and haptic
 - Kneuro user app

Operating the Device

Fully charge the device before first use.

Power ON/OFF

Press the control button for 8 seconds. The device will play 1 beep and 1 vibration.

The LED indicator light on the back of the device will turn green when powered ON.

The LED indicator light will turn off when the device is powered OFF.

Figure 6

Manual Lock

With the device in full extension, press the control button for 2 seconds to engage or disengage Manual Lock. A beep and short vibration will confirm the action.

Quick Access

This feature allows users to quickly switch between standard mode and the assigned activity mode of their choosing. If activated, switching is done by double-tapping the control button.

Before switching, the user must be seated or in a safe position for the desired mode.

The quick access feature is OFF by default. Turn ON the desired activity mode in the Kneuro app. Only one activity mode can be assigned.

Connecting to the Kneuro App

Turn ON the device.

Download and open the Kneuro app on an iOS or Android device.

System Requirement: iOS 12 or above/Android 5.0 or above

Select “I’m a Patient” on the welcome screen and enter login credentials, or create an account.

To create a user account, ensure the prosthetist has created a user profile. Patients will need to use their phone number or email address to register their account and complete account set up.

In the user account, after the device is paired, the user will be able to make small adjustments to their resistance parameters, change mode selection, check battery status, and monitor their step count.

The device and your Bluetooth must be on prior to pairing.

For assistance logging in, contact BrainRobotics at +1 (512) 969-2989 or support@brainrobotics.com.

Sitting Down

With equal weight on both legs, push weight back into the chair and initiate a yielding knee.

The user will be able to load the device and use their body weight to maintain a constant resistance in knee flexion until the desired position is reached.

Walking on the Level Ground

The device provides hydraulic resistance to support the user in stance phase. When the stance release criteria are met, the device transitions into swing phase.

The device will adjust to the user's self-selected speed.

Resistance and angular parameters can be adjusted in the Kneuro app with a clinician account.

Standing

With feet evenly placed and weight equally distributed across the sound side and prosthetic limb, instruct the user to lean forward with their upper body.

The device will provide resistance to support the user as they move into an upright position.

Sitting

When in a seated position, if the thigh is parallel to the ground for 2 seconds, the device resistance will reduce to a minimum value, allowing the user to position the prosthesis with ease.

Stance Lock

When the user moves into extension, the lock will disengage.

This function allows users to bear weight on the device comfortably.

This can be engaged or disengaged depending on the user's preference and mobility levels. Please set accordingly. Stance Lock is default ON.

Descending Ramps

Take the first step down with the prosthetic limb, allowing the device to slowly flex under the high resistance provided by the yielding support.

Please use caution when instructing the user to descend ramps.

Ensure care is taken in initial alignment procedures to ensure stance release does not trigger prematurely in ramp descent.

Users are advised not to walk on excessively steep slopes greater than 25°, as it may result in decreased yielding resistance and loss of support.

Yielding resistance can be adjusted in the Kneuro app by a certified prosthetist.

Walk-to-Run

Once the user has mastered level ground ambulation with the device, instruct the user to move into more advanced movement patterns, gradually increasing their cadence from a walk into a run.

The user should have confidence that the device will support them in activities of daily living, sports, and emergency settings.

Step-Over-Step Stair Ascent

To switch into this mode, the user must quickly extend at their hip followed by a rapid hip flexion movement to achieve the desired height of the step.

Once adequate hip flexion is obtained, the device will increase extension resistance, ensuring the foot is appropriately placed on the stair.

Upon bearing weight onto the device, it will transition to its highest flexion resistance to prevent further flexion of the device. The user will contract their hip extensors to fully extend the device and should place the contralateral limb on the next stair.

Each subsequent step will require less effort to initiate step-over-step stair ascent.

Always use a handrail when ascending stairs.

Stair Extension Resistance can be adjusted in the Kneuro app by a certified prosthetist.

Stair Descent

Take the first step with the prosthetic limb, placing the foot half off the stairs to allow for a safe and secure rollover of the prosthetic foot. Roll the foot over the edge of the stair, allowing the device to slowly flex under high resistance. Place the sound limb onto the next step. Practice this movement pattern with the user to ensure appropriate foot placement for every step.

Always use a handrail when descending stairs.

Yielding resistance can be adjusted in the Kneuro app by a certified prosthetist.

Stumble Recovery

Increased flexion dampening will occur if swing phase flexion or extension is interrupted.

This high resistance will provide the user time to prevent, or prepare, in case of a fall.

Mode Selection

The user can select modes via the Kneuro app. Ensure the user is in a seated and safe position prior to switching modes.

To switch back to standard mode, use the Kneuro app.

The user can designate one mode to the quick access feature activated by the control button. This feature must first be enabled by a clinician in the Kneuro app. See Audio, Vibration and LED Feedback Signals.

Cycle Mode

In cycle mode, there is no flexion or extension resistance. This mode has no adjustable parameters.

Ensure the user is in a seated and safe position prior to switching modes.

Golf Mode

These parameters can be adjusted based on the user preferences. The golf mode default settings are:

- Flexion resistance: 60°
- Extension resistance: 30°
- Flexion limit angle: 20°

Custom Mode A/B/C

There are three modes available for unique customization. Ensure the user is safe and specified settings are adjusted appropriately. Custom modes are disabled until programmed for the user.

- ☒ The parameters set in modes are only adjustable in the Kneuro app by a certified prosthetist. The user can select programmed modes within the app.

Battery Information

 This product is equipped with a lithium-ion battery. For your safety and to ensure proper function of the device:

- **Do not** disassemble.
- **Do not** apply pressure or puncture the device.
- **Do not** use or store in temperatures above 104°F (40° C).

The device's battery is a rechargeable lithium-ion battery.

A fully charged battery will last approximately 4-5 days, depending on usage and activity level.

- ☒ The device takes 4 hours to fully charge from a depleted state.
- ☒ The user may experience a shorter battery life when the device is used in extreme temperatures.

Charging the Battery

Charge the device daily for optimal performance.

Plug the power supply into a wall socket or car adapter.

Connect the USB-C cable to the charging port on the posterior aspect of the device.

- While charging, the LED indicator light will be yellow.
- When charging is complete, the LED indicator light will be green.

 **Do not** use the device while charging. The user must doff the prosthesis before charging to avoid electrocution and/or injury to the residual limb.

- ☒ The battery cannot be overcharged
- ☒ The exact battery level percentage is displayed in the Kneuro app

Battery Level	LED Signal	Vibration Signal	Audio Signal
Battery Level < 20%	Yellow	Medium signal x5	
*Battery Level < 5%	Blinking Yellow	Short signal x5	
Battery Level 0%	None		
*Safe mode is activated when the battery level is < 5%			


Figure 7

 The battery can be charged while the device is ON or OFF

Audio, Vibration and LED Feedback Signal

Type	Function	Trigger	Audio Signal	Vibration Signal	*LED Signal
Normal use	Power ON	Press control button for 8 seconds	Medium signal	Medium signal x1	Green
	Power OFF			-	LED off
	Switch to default mod	Double press control button or use the Kneuro app	Medium signal x3 (3 second countdown) Long signal x1 (switch completed)		Green
	Switch to custom mode (Quick access feature)		Medium signal x3 (3 second countdown) Short signal x2 (switch completed)		
	Stance Lock	Hold flexed stance for 2 seconds (20°-40°)	-	-	
	Manual lock/Unlock	Press control button for 2 seconds	Short signal x1	-	
Recharging and maintenance	Charging	-	-	-	Yellow - charging
					Green - complete
	Problem report	-	Start and stop; Short signal x2	-	White
	Firmware update				Blinking white
Idle mode	Idle	At rest for 30 min	-	-	Blinking blue
	Awake	-	-	Long signal x1	Green
Device connection	Connected	-	-	-	Green
	Disconnected	-	-	-	
Parameter Setting	Setting updated		Medium signal x 1		

Type	Function	Trigger	Audio Signal	Vibration Signal	*LED Signal
☑ For signal meanings, see Warning Signals. *If the LED indicator light is malfunctioning, contact customer support at +1 (512) 969-2989 or support@brainrobotics.com					

Warning Signals

Function	Event	Audio	Vibration	LED	User Action
Low Battery	Below 20%	Medium signal x5		Yellow	Recharge within 6 hours
	Below 5%	Short signal x5		Blinking yellow	Recharge as soon as possible
Malfunction	Critical error	Continuous for 8 seconds. Repeats until malfunction is resolved or device is powered OFF		Red	Power cycle the device. If unresolved, stop use and contact support immediately
Risk	Overheating	-	Medium signal, short signal. Repeats every 5 seconds until device cools		Reduce activity level immediately

Troubleshooting

Safe Mode

The device will enter Safe Mode if:

- Battery level is <5%
- Device is OFF

While in Safe Mode, the device will maintain a constant resistance. The pre-programmed default is set for high resistance. Programmed parameters only apply when triggered by low battery and when the device is OFF.

In Safe Mode, the Kneuro will not release into swing phase. It will still enable yielding resistance that is adjustable when programming the Kneuro ensuring that the user can safely sit down, and descend ramps and stairs.

Safe Mode will automatically disable once the battery has been recharged, or the issue is resolved.

Reducing the Safe Mode resistance is not recommended and can cause injury. If resistance is adjusted, ensure the user can stand without the knee flexing, that they are aware of the risk, and that they are prepared for the level of resistance.

The device will not switch modes while in Safe Mode.

Overheating

The device may overheat during prolonged periods of high activity. While overheating does not impact immediate performance, continued use in this state may lead to damage of the hydraulic unit.

When overheating occurs, the device will emit a medium and short vibration every 5 seconds, and the LED indicator will display a steady red light until it cools.

Reduce activity and refrain from use until the device cools down.

Do not touch the hydraulic unit when overheated.

Normal operating temperature range: 14–104°F (–10–40°C).

Hydraulic Noise or Inconsistent Resistance

Noise in the hydraulic system or an inconsistent resistance may be due to an accumulation of air in the hydraulic system.

The device has an automatic air release that will relieve the air accumulation within 10-20 steps.

To prevent this issue from occurring, please store the device in a vertical and fully extended position.

Bluetooth Not Connecting

Power the device OFF then turn back ON.

Reconnect the device to Bluetooth. (Bluetooth connectivity is always ON when the device is turned ON.)

If turning the device OFF and back ON does not work, shut down the Kneuro app.

Be sure the device is located within 15 ft of the iOS or Android device and that the Kneuro app has permission to access your location data.

Reopen the Kneuro app and attempt to reconnect again.

If Bluetooth continues not connecting to device after options above, contact BrainRobotics or your local distributor for additional troubleshooting and support.

Cleaning

The device must be cleaned to avoid corrosion or damage. To clean the device, dampen a cloth with warm tap water and wipe down the device and all components. Pat dry with a lint-free cloth. Clean the device regularly to ensure continued function.

Storage

Please store the device in a vertical position, as storing in a horizontal position may lead to accumulation of air within the hydraulic system. The device must be stored in a fully extended position.

 The device must be stored in a dry and cool environment. Storage temperature range is 23°F (-5°C) to 86°F (30°C) with maximum relative humidity at 95%.

Maintenance

It is important to have both the device and the overall prosthesis regularly inspected by a certified prosthetist. The frequency of these inspections should be determined based on the user's activity level. The manufacturer recommends a service interval of 36 months.

Regularly check the device for any signs of damage. If damage is noticed, inform your certified prosthetist immediately. If you have questions or concerns, contact BrainRobotics customer support at +1 (512) 969-2989 or support@brainrobotics.com.

Disposal

 The crossed out wheeled bin symbol with under bar shown on the device, device packaging or accompanying documents indicates the product requires appropriate treatment, collection and recycle for waste electrical and electronic equipment (WEEE) under the Directive 2012/19/EU.

This device should not be disposed of as unsorted household waste. The device contains lithium batteries and should be disposed of according to local regulations for electronic equipment.

The charger should be disposed of according to local regulations for electronic equipment or through a community electronics recycling center.

Magnetic Resonance (MR) Compatibility

The device is MR Unsafe and should not be worn near or during imaging procedures (i.e. CT scan, MRI, etc.).

Electronic Compatibility Declaration

This device requires special precautions regarding electromagnetic compatibility (EMC), and this device must be installed and used in accordance with the electromagnetic compatibility information specified in this manual.

The charging cable must be 3 feet (1 meter) in length to meet electromagnetic emission and immunity requirements.

Portable and mobile radio frequency communication equipment may affect this device.

The use of accessories and cables (transducers) other than those specified and sold as spare parts for internal components may result in increased emissions or decreased immunity of this device or system.

The use of this device adjacent to or stacked with other equipment should be avoided. If such is unavoidable, this device should be observed to ensure that it functions properly in the configuration in which it is used.

Working	During the immunity test, there is no significant change in the performance of this device during the measurement process; temporary degradation of performance may be allowed during the electrostatic discharge, but it can return to the pre-discharge operating condition within 10 seconds.
Charging	The LED indicator will show a solid yellow light while the battery is charging. The LED indicator will show a solid green light when the battery is fully charged.

Compliance

The device has been designed and tested according to the applicable international standards or in-house defined standards when no international standard exists or applies.

This device is designed for operation in the following electromagnetic environments, and the purchaser or the user shall ensure that it is used in such an electromagnetic environment.

Recommended protection distances between portable and mobile radio frequency communication equipment and this device. This device is intended for use in electromagnetic environments where radiated, radio-frequency disturbances are controlled. Depending on the maximum rated output power of the communication equipment, the purchaser or the user can prevent electromagnetic interference by maintaining a minimum distance between portable and mobile radio frequency communication equipment (transmitters) and this device as recommended below.

Electromagnetic Emissions			
Standard	Test Item	Limit	Result
EN 60601-1-2:2015 EN 60601-1-2:2015/A1:2021 IEC 60601-1-2:2014+AMD1:2020	Conducted disturbance - A C mains ports	Class B group	PASS
	Radiated disturbance	Class B group	PASS
IEC 61000-3-2:2005 +A1:2008+A2:2009	Harmonic current disturbance	Class A	N/A
IEC 61000-3-3:2013	Voltage fluctuations and flicker	-	PASS
IMMUNITY (IEC 60601-1-2:2014+AMD1:2020 and EN 60601-1-2:2015 + A1:2021)			

Electromagnetic Interference Immunity			
Basic Standard	Test Item	Test Specification	Result
IEC 61000-4-2:2008	Electrostatic discharge	Air discharge: ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV (Direct) Contact discharge: ± 8 kV (Direct/indirect)	PASS
IEC 61000-3-2:2005 +A1:2008+A2:2009	Continuous RF electromagnetic field disturbance (Professional healthcare facility environment)	3 V/m, 80%, 1 kHz, AM, 80 MHz – 2,7 GHz; and Note (2)	N/A
	Continuous RF electromagnetic field disturbance (Home healthcare environment)	10 V/m, 80%, 1 kHz, AM, 80 MHz – 2,7 GHz; and Note (2)	PASS
IEC 61000-4-4:2012	Electrical fast transient/burst (AC input ports)	2.0 kV 5/50 ns, 100 kHz-	PASS
	Electrical fast transient/burst (Input d.c. power PORT intended to be connected permanently to cables longer than 3 meters)	2.0 kV 5/50 ns, 100 kHz	N/A
	Electrical fast transient/burst (Signal input/output parts PORT- whose max cable length is less than 3m are excluded.)	1 kV 5/50 ns, 100 kHz	N/A

Electromagnetic Interference Immunity			
Basic Standard	Test Item	Test Specification	Result
IEC 61000-4-5:2014 +A1:2017	Surge (AC input ports)	0.5 kV, 1 kV (Line-to-line) 0.5 kV, 1 kV 2 kV (Line-to-ground) 1.2/50 us	PASS
	Surge (Input d.c. power PORT- intended for permanent connection to cable longer than 3m)	0.5 kV, 1 kV (Line-to-line) 0.5 kV, 1 kV 2 kV (Line-to-ground) 1.2/50 us	N/A
	Surge (Signal input/output parts PORT- intended to connect directly to outdoor cables)	2 kV (Line-to-ground)	N/A
IEC 61000-4-6:2013	Conducted disturbances induced by RF fields (AC input ports)	150 kHz - 80 MHz, 80%, 1 kHz, 3 V (unmodulated, r.m.s.) 6 V (unmodulated, r.m.s.) in ISM bands between 0,15 MHz and 80 MHz	PASS
	Conducted disturbances induced by RF fields (Input d.c. power PORT- intended for permanent connection to cable longer than 3m)	150 kHz - 80 MHz, 80%, 1 kHz, 3 V (unmodulated, r.m.s.) 6 V (unmodulated, r.m.s.) in ISM bands between 0,15 MHz and 80 MHz	N/A
	Conducted disturbances induced by RF fields (PATIENT coupling PORT-All PATIENT-COUPLED cables shall be tested, individually or bundled)	150 kHz - 80 MHz, 80%, 1 kHz, 3 V (unmodulated, r.m.s.) 6 V (unmodulated, r.m.s.) in ISM bands between 0,15 MHz and 80 MHz	N/A
	Conducted disturbances induced by RF fields (Signal input/output parts PORT- whose max cable length is less than 1m are excluded)	150 kHz - 80 MHz, 80%, 1 kHz, 3 V (unmodulated, r.m.s.) 6 V (unmodulated, r.m.s.) in ISM bands between 0,15 MHz and 80 MHz	N/A
IEC 61000-4-8:2009	Power frequency magnetic field	50 Hz, 30 A/m	PASS
IEC 61000-4-39:2017	Proximity magnetic fields	134,2 kHz, Pulse modulation 2,1 kHz, 65 A/m; 13,56 MHz, Pulse modulation 50 kHz, 7,5 A/m	N/A

Electromagnetic Interference Immunity			
Basic Standard	Test Item	Test Specification	Result
IEC 61000-4-11:2004 +A1:2017	Voltage dips	Voltage dips: 0% UT; 0,5 cycle A 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	PASS
		0% UT; 1 cycle Single phase: at 0° PASS	PASS
		70% UT; 25/30 cycles Single phase: at 0°	PASS
	Interruptions	Voltage interruptions: 0% UT; 250/300 cycle	PASS
	ISO 7637-2:2011	Electrical transient conduction along supply lines	ISO 7637-2:2004 Pulse 1
ISO 7637-2:2004 Pulse 2a			N/A
ISO 7637-2:2004 Pulse 2b			N/A
ISO 7637-2:2004 Pulse 3a			N/A
ISO 7637-2:2004 Pulse 3b			N/A

☒

(1) "N/A " denotes test is not applicable in this Test Report.

☒

(2) Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment.

RF Wireless Communications Compliance					
Test Frequency (MHz)	Band (MHz)	Service	Modulation	Max power (W)	Immunity test level (V/m)
385	380-390	TETRA 400	Pulse modulation 18 Hz	1,8	27
450	430-470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	28
710	704-787	LTE Band 13, 17	Pulse modulation 217 Hz	0,2	9
745					
780					

RF Wireless Communications Compliance					
Test Frequency (MHz)	Band (MHz)	Service	Modulation	Max power (W)	Immunity test level (V/m)
810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	28
870					
930					
1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; U MTS	Pulse modulation 217 Hz	2	28
1845					
1970					
2450	2400-2570	Bluetooth, W LAN, 802.11 b/g /n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	28
5240	5100-5800	W LAN 802.11 a/n	Pulse modulation 217 Hz	0,2	9
5500					
5785					
(3) The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz					
(4) Only applicable to the device which contains magnetically sensitive components.					
(5) The measurement result for the sample received is <PASS> according to < IEC 60601-1-2:2014+AMD1:2020, EN 60601-1-2:2015+A1:2021> when < Simple Acceptance > decision rule is applied.					

USA-Federal Communications Commission (FCC)

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy. If not installed and used in accordance with the instructions, it may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a specific installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try and correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the distance between the equipment and the receiver.
- Connect the equipment to an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Exposure to Radio Frequency Radiation – This device must not be co-located or operating in conjunction with any other antenna or transmitter.

Disclaimer

Before using this device, the prosthetist and user are expected to carefully read this "Instructions for Use" to protect their legal rights or interests. BrainRobotics reserves the right to update this document.

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorized modification of the device.

In the process of using this device, the user should strictly observe and implement the requirements in accordance with the descriptions and instructions provided in this document. For all personal injuries, accidents, property damage, legal disputes, and other adverse events caused by or during the use as mentioned in the safety instructions, the user shall bear the corresponding responsibilities and losses, and the manufacturer shall not assume any responsibility.

All names described in the documents are unconditionally protected by applicable trademark law and are the property of the respective owners. The brands, trade names, or company names described herein may be registered brands and are the property of the respective owners. Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

Bluetooth is a registered trademark of Bluetooth SIG, Inc.

BrainRobotics disclaims responsibility for any violation of laws and regulations caused directly or indirectly by the user's use of this device.

Serious Incident Reporting

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

-  'Serious incident' is defined as any incident that directly or indirectly led, may have led, or might lead to any of the following: (a) the death of a patient, user, or other person; (b) the temporary or permanent serious deterioration of a patient's, user's, or other person's state of health; or (c) a serious public health threat.

If you have any questions or concerns that were not addressed in this manual, please contact the manufacturer.

Technical Specifications

Device basic parameters	Device degree of freedom	1 DOF
	Device height	11.42in (290mm)
	Build height	9.15in (232.6mm)
	Device width	2.99 (76mm)
	Device depth	4.21in (107mm)
	Device weight	3.53lbs (1.6kg)
	Device weight limit	330 lbs (150kg)
	Device IP rating	IP67 (water resistant)
	Power adapter IP rating	IP22
	Max Bluetooth connection range	16ft (5m)
	Firmware upgrade	OTA
Data transmission	Wireless protocol	BLE5.0
	Frequency range	2.4000-2.4835GHz
	Modulation	GFSK
	Max output power	+8dBm
Battery	Battery type	Lithium-ion/ 3000mAh/ 3.7V
	Battery life (≥80% original capacity)	≥500 Cycles
	Battery level after 1h charging	~ 50%
	Battery level after 2h charging	~ 80%
	Battery level after 4h charging	100%

Power adapter	Adapter input voltage	AC 90-264V
	Adapter input power frequency	47-63Hz
	Adapter output voltage	DC 5V
	Adapter output current	0.01-2A
Reliability standard	Executive standard	ISO 10328-P7
Environment requirements	Long term storage temperature/ humidity	23-86°F (-5-30°C) Max 95% relative humidity 70-106 kPa (up to 3,000m without pressure equalization)
	Working temperature/ humidity/ air pressure	14-104°F (-10-40°C) Max 95% relative humidity 70-106 kPa (up to 3,000m without pressure equalization)
	Charging temperature/ humidity/ air pressure	50-113°F (10-45°C) Max 95% relative humidity 70-106 kPa (up to 3,000m without pressure equalization)
	Transportation temperature/ humidity	-4°F-140°F (-20-60°C) Max 95% relative humidity 70-106 kPa (up to 3,000m without pressure equalization)

Symbol Table

Symbol	Definition	Symbol	Definition
	Note		Bluetooth a registered trademark of Bluetooth SIG, Inc.
	Control button		Recycle
	USB plug-in		Consider Instructions for Use when operating device and electronic location of Instructions for Use
	Wireless connection		Single Patient Use
	Caution hot surface		United Kingdom (UK) declaration of conformity

Symbol	Definition	Symbol	Definition
	Non-ionizing radiation		Manufacturer
	Device contains electronic components and/or batteries that should not be disposed in regular waste		Type BF applied part
	Meets compliance with requirements according to "FCC Part 15" (USA)		European Authorized Representative
	Must read Instructions for Use		Importer
IP67	Dust-tight, protect against		Prescription use only (USA)
	Warning		Indicates device is a medical device
	Serial Number		Catalog number
	Direct current		Country of Origin and date of manufacture
	Fragile		Unique Device Identifier
	Protect from moisture		United Kingdom UK responsible Person
	Declaration of conformity according to the applicable European directives		Unsafe for MR environment
	This side up		

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Bluetooth ist eine eingetragene Marke von Bluetooth SIG, Inc.