ENGLISH

MD Medical Device

INTENDED USE

The device is intended to immobilize the shoulder and/or elbow

Indications for use

For shoulder conditions that may benefit from immobilizing the shoulder, such as:

- Pre or post-operative
- Soft tissue repairs/strains
- Shoulder instabilities/injuries

No known contraindications.

GENERAL SAFETY INSTRUCTIONS

Read these instructions carefully before use. Keep them for future reference.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities. Never wear the device directly over an open wound. Inspect the skin under the device's contact area before and after each use.

The patient should stop using the device and contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

USAGE

Device Application

1. Unfold device and remove all straps from the belt.
2. Wrap the belt around the waist and secure contact closure in the back (Fig. 1).
3. Fasten the shoulder strap to the waist belt in front and back, making sure the strap is over the healthy shoulder (Fig. 2, 3).
4. Adjust the shoulder strap length for proper fit (Fig. 4).
5. Fasten arm band above the elbow of the injured shoulder and secure contact patch to the waist belt (Fig. 5, 6).

6. Fasten the wrist band around the wrist and secure the contact patch to the waist belt (Fig. 7, 8, 9).

Device removal

Undo the arm and wrist band straps and take your arm out of the straps. Unfasten the belt to remove the device. Shoulder strap, arm and wrist band do not have to be removed from the belt.

Cleaning and care

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry, do not expose to direct heat.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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