1C11 Terion K2


DE	Gebrauchsanweisung	5
EN	Instructions for use	12
FR	Instructions d'utilisation	19
IT	Istruzioni per l'uso	27
ES	Instrucciones de uso	35
PT	Manual de utilização	43
NL	Gebruiksaanwijzing	51
SV	Bruksanvisning	59
DA	Brugsanvisning	66
NO	Bruksanvisning	73
FI	Käyttöohje	80
PL	Instrukcja użytkowania	87
HU	Használati utasítás	95
CS	Návod k použití	103
RO	Instrucțiuni de utilizare	110
HR	Upute za uporabu	118
SL	Navodila za uporabo	126
SK	Návod na používanie	133
BG	Инструкция за употреба	140
TR	Kullanma talimatı	149
EL	Οδηγίες χρήσης	156
RU	Руководство по применению	164
JA	取扱説明書	173
ZH	使用说明书	180
KO	사용 설명서	186

1

2

7

INFORMATION

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INFORMATION

Last update: 2015-10-29

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

1.1 Construction and Function

The 1C11 Terion K2 prosthetic foot was developed for the needs of users with mobility grades 1 and 2. The foot is suitable for walking on various surfaces.

The functional properties of the prosthetic foot are achieved through the spring element made of fibre-reinforced polymer and a functional foam base. The prosthetic foot permits perceptible plantar flexion at heel strike and a natural rollover movement.

1.2 Combination possibilities

This prosthetic component is part of the Ottobock modular system. As a result of its design, it may be combined with other prosthetic components of the modular system. You can find possible restrictions in this section. Contact the manufacturer with questions.

2 Application

2.1 Indications for use

The product is to be used solely for lower limb prosthetic fittings.

2.2 Area of application

The product is recommended for mobility grade 1 (indoor walker) and mobility grade 2 (restricted outdoor walker).

Maximum body weight according to stiffness [kg]									
Size [cm]	22	23	24	25	26	27	28	29	30
Stiffness 1	55	55	55	55	–	–	–	–	–
Stiffness 2	75	75	75	75	75	75	75	–	–
Stiffness 3	100	100	100	100	100	100	100	100	100
Stiffness 4	–	–	125	125	125	125	125	125	125
Stiffness 5	–	–	150	150	150	150	150	150	150
Stiffness 6	–	–	–	–	–	175	175	175	175

2.3 Environmental conditions

Allowable environmental conditions
Temperature range for use: -10 °C to +60 °C
Relative humidity: no restrictions
Chemicals/moisture: fresh water, salt water, chlorine water, soapsuds, sweat, urine
Solids: dust, sand, highly hygroscopic particles (e. g. talcum)
Clean the product after contact with humidity/chemicals/solids, in order to avoid increased wear and damage (see Page 18).

2.4 Service Life

This product has been tested by the manufacturer with 2 million load cycles according to ISO 22675. Depending on the patient's activity level, this corresponds to a service life of 2 to 3 years.

3 Safety

3.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

3.2 General safety instructions

 **CAUTION**

Exceeding the service life and reuse on another patient

Risk of injury due to loss of functionality as well as damage to the product

- ▶ Ensure that the approved service life is not exceeded (see section "Service life").
- ▶ Only use the product for a single patient.

 **CAUTION**

Excessive strain on the product

Risk of injury due to breakage of load-bearing components

- ▶ Use the prosthetic components according to the MOBIS classification (see section "Area of application").

 **CAUTION**

Unallowable combination of prosthetic components

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved according to the section "Combination possibilities".
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

 **NOTICE**

Mechanical damage to the product

Change in or loss of functionality due to damage

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).

- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

NOTICE

Use under unallowable environmental conditions

Damage to product due to unallowable environmental conditions

- ▶ Do not expose the product to unallowable environmental conditions (see the section "Environmental Conditions").
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e. g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

Signs of changes in or loss of functionality during use

Reduced spring effect (e.g. decreased forefoot resistance or changed rollover characteristics) is an indication of loss of functionality.

4 Scope of delivery

1C11 Terion K2		
Quantity	Designation	Reference number
1	Instructions for use	–
1	Prosthetic foot	–
1	Toe insert	–
1	Protective sock (size dependent)	22 to 28: SL=Spectra-Sock-2 29 to 30: SL=Spectra-Sock
1	Footshell	2C12=*
1	Connection cap	2C13=*

5 Preparation for use

⚠ CAUTION

Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- ▶ Observe the alignment and assembly instructions.

NOTICE

Grinding the prosthetic foot or footshell

Premature wear resulting from damage to the product

- ▶ Do not grind the prosthetic foot or footshell.

5.1 Alignment

INFORMATION

There is a plastic adapter cover on the proximal connection of the prosthetic foot. This protects the connecting section from scratches during the alignment and trial fitting.

- ▶ Remove the adapter cover before the patient leaves the workshop/fitting area.

5.1.1 Applying/removing the footshell

Inserting the foot into the footshell

- > **Recommended tools:** 2C100 changing device for footshell
- 1) Pull the Spectra sock over the prosthetic foot, avoiding any wrinkles.
- 2) Slide the prosthetic foot into the footshell. Insert the carbon spring into the slot on the toe insert (see fig. 4).
- 3) Push the prosthetic foot into the heel area of the footshell until it locks into place (see fig. 6).
- 4) Place the connection cap on the footshell (see fig. 7).

Removing the footshell

- > **Recommended tools:** 2C100 changing device for footshell
- 1) Push the heel area of the footshell to the rear and pull the prosthetic foot up.
- 2) Remove the prosthetic foot from the footshell.

5.1.2 Bench Alignment

TT bench alignment

Bench alignment process

Recommended tools:

662M4 goniometer, 743S12 heel height measuring apparatus, 743A80 50:50 gauge, alignment apparatus (e.g. 743L200 L.A.S.A.R. assembly or 743A200 PROS.A. assembly)

Perform the alignment and assembly of the prosthetic components in the alignment apparatus according to the following specifications:

Bench alignment process	
1	Heel height: Effective heel height (shoe heel height – sole thickness of fore-foot) + 5 mm
2	Sagittal alignment: Alignment reference line on a–p markings of the connection cap/footshell (see fig. 2)
3	Frontal alignment: Alignment reference line on m–l markings of the connection cap/footshell (see fig. 2)
4	Connect the prosthetic socket and prosthetic foot using the chosen adapters. Be sure to follow the adapter instructions for use during adjustment and installation.
5	Sagittal alignment: Determine the centre of the prosthetic socket with the 50:50 gauge. Align the prosthetic socket centrally to the alignment reference line. Socket flexion: Residual limb flexion of the patient + 5°
6	Observe the abduction or adduction position.

TF bench alignment

- Observe the information in the prosthetic knee joint instructions for use.

5.1.3 Static Alignment

- Ottobock recommends fitting the prosthesis using the L.A.S.A.R. Posture.
- If necessary, the alignment recommendations (TF modular leg prostheses: **646F219***, TT modular leg prostheses: **646F336***) may be requested from Ottobock.

5.1.4 Dynamic Trial Fitting

- Adapt the setup of the prosthesis in the frontal plane and the sagittal plane (e.g. by making angle or slide adjustments) so that correct heel contact and optimal roll-over are ensured.
- **TT fittings:** Make sure that physiological knee movement is achieved after the heel makes impact and the leg begins to bear weight.
- Remove the adapter cover from the prosthetic foot after completion of the dynamic trial fitting and gait training exercises.

5.2 Optional: Installing the foam cover

The connection cap serves as a detachable connection between the foam cover and footshell, and snaps precisely into the proximal brim of the footshell.

> **Required materials:** 634A58 isopropyl alcohol, 636N9 contact adhesive or 636W17 plastic adhesive

- 1) Cut the foam cover to length according to the patient's measurements, allowing for distal compression of about **1 cm** for transtibial foam covers and of about **4 cm** for transfemoral foam covers.
- 2) Pull the foam cover over the prosthesis.
- 3) Join the connection cap to the footshell and mount the prosthetic foot.
- 4) Mark the outer contours of the connection cap on the distal face of the foam cover.
- 5) Dismount the prosthetic foot and remove the connection cap from the footshell.
- 6) Clean the connection cap with isopropyl alcohol.
- 7) Glue the connection cap to the face of the foam cover as marked and allow the glued joint to set (approx. **10 minutes**).
- 8) Mount the prosthetic foot and adapt the exterior cosmetic cover shape. Take into account compression due to cosmetic stockings or SuperSkin.

6 Cleaning

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

7 Maintenance

- ▶ The prosthetic components should be inspected after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

8 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in the country of use regarding return, collection and disposal procedures.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

10 Technical data

Sizes [cm]	22	23	24	25	26	27	28	29	30
Heel height [mm]	10 ±5								
Footshell	narrow			standard					
System height with footshell [mm]	40		41		46		48		51
Structural height [mm]	58		59		65		67		70
Product weight without footshell [g]	237		283		325			423	
Product weight with footshell [g]	410	445	491	509	573	580	635	709	744
Max. body weight [kg]	100		150		175				
Mobility grade	1 + 2								

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Otto Bock HealthCare GmbH