INDLIGH

MD Medical Device

INTENDED USE

The device is intended as a connection component of a prosthetic system that replaces a missing lower limb. Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional. The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Lower limb loss, amputation, or deficiency
- No known contraindications
- The weight limit for the device is 166kg.

The device is for high impact use, e.g., walking and occasional running.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

Warning: Only use Ossur replacement parts for the device for ensured durability The device is for single patient use.

ASSEMBLY INSTRUCTIONS

Offset Adjustment:
1. Optional: Security provided from the limiting screw restricts shift movement in one direction. To provide shift movement in the other direction, remove the limiting screw and locate into the second hole
2. Loosen slide adjustment screw
3. Set the shift adapter to desired offset
4. Secure slide adjustment screws with medium strength threadlocker and torque until both have reached desired torque.

Caution: When sliding in frontal plane the slide adjustment screws must be positioned anteriorly.

Ensure the Shift Adapter is positioned with the female shift section proximally and the male shift section distally on the prosthesis. This orientation provides maximum strength from the adapter.

Warning: For secure fastening use medium strength threadlocker and apply rated torque.

Warning: When extreme alignment is necessary, use extended set screws to achieve secure fastening.

Tighten the screws with the following torque:

- Slide Adjustment Screws: 12 Nm
- Slide Limit Screws: 5 Nm

USAGE

Cleaning and care

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

Environmental Conditions

The device is Waterproof. A Waterproof device can be used in a wet or humid environment and submerged in up to 3 meters deep water for a maximum of 1 hour. It can tolerate contact with: Salt water, chlorinated water, perspiration, urine, and mild soaps. It can also tolerate occasional exposure to sand, dust, and dirt. Continuous exposure is not allowed. Dry with a cloth after contact with fresh water or humidity.

Clean with fresh water after exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

REPORT OR SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Ossur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance

This device has been tested according to ISO 10328 standard to three million load cycles.

Depending on patient activity, this may correspond to 3-5 years of use.

ISO 10328 - P7 - 166 kg *) ▲
<small>*) Body mass limit not to be exceeded! For specific conditions and limitations of use see manufacturer's written instructions on intended use!</small>
▲

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