
Instructions for Use


---

LP VARI-FLEX<sup>®</sup>





VLPExyyzs(BR)





		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	9
FR	Notice d'utilisation	13
ES	Instrucciones para el uso	17
IT	Istruzioni per l'uso	21
DA	Brugsanvisning	25
SV	Bruksanvisning	29
EL	Οδηγίες χρήσης	32
NL	Gebruiksaanwijzing	36
PT	Instruções de Utilização	40
JA	取扱説明書	44
ZH	中文说明书	47
KO	사용 설명서	50

**BENCH ALIGNMENT**

# ENGLISH

---

## LP VARI-FLEX® FOOT SYSTEM (Figure 1)

- Low profile foot module
- Build height: Pyramid adapter 66,1 – 77,2 mm (2 1/2" - 2 7/8")
- Maximum weight 166 kg (365 lbs)
- Available in 10 mm (3/8") and 19 mm (3/4") heel height
- Available in sizes 22-30 cm

**NOTE:** Select foot category according to impact level.

## BENCH ALIGNMENT (Figure 2 & 3)

- Fit foot with selected cover.
- Adjust to appropriate heel height (using the shoe).
- Introduce appropriate socket angles flexion/ extension and abduction/adduction.
- Divide the foot cover length into 3 equal portions.
- The load line should fall at the junction of the posterior and the middle third portions as shown in the diagram.

## DYNAMIC ALIGNMENT

The heel stores energy after initial contact, slowly releasing it at midstance. The forward momentum generated by this action results in the toe being loaded for optimum energy release at terminal stance.

The heel to toe action can be influenced by

- A-P positioning of foot
- Dorsi-Plantarflexion
- Heel Stiffness

## CATEGORY SELECTION CHART

Please refer to the selection charts below to determine the appropriate stiffness required according to Össur recommendations.

WEIGHT KG	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147	148-166
WEIGHT LBS	99-115	116-130	131-150	151-170	171-194	195-221	222-256	257-287	288-324	325-365
Low Impact Level	1	1	2	3	4	5	6	7	8	9
Moderate Impact Level	1	2	3	4	5	6	7	8	9	N/A
High Impact Level	2	3	4	5	6	7	8	9	N/A	N/A

## HEEL STIFFNESS

### Heel too soft

#### Symptoms

- Foot will come to flat position too early (amputee feels he/she is sinking into a hole).
- Extra energy is required to climb up over the toe.
- Toe will feel too stiff.
- Knee may hyperextend.

#### Action

- Add polyurethane heel wedge.
- Shift socket anterior (or foot posterior).

## Heel too hard

### Symptoms

- Rapid heel to toe movement
- At initial contact the amputee has poor control of his prosthesis
- Minimal energy return feeling
- Knee may become unstable

### Action

- Replace heel wedge with a smaller one.
- Shift socket posterior or foot anterior.

## HEEL WEDGES CM150011 (Figure 4)

The wedge can influence the heel to toe function. Start altering the heel resistance to improve foot response by adding a heel wedge. The 2°, 4° & 6° wedges can be interchanged to customize the stiffness and achieve the desired functional characteristics. A combination of wedges may be used

1. 2° – yellow
2. 4° – red
3. 6° – blue

## Temporary Wedge Placement

- Cut the wedge to the width of the foot module.
- Roughen the upper and lower surface of the wedge with abrasive paper.
- Place the wedge in the angle of heel and foot module.
- Secure in position bond to upper surface (**Figure 4**) with tape wrapped around the foot module (**Figure 5**).

## Permanent Wedge Placement

- Apply adhesive on the upper side of the wedge only.
- Locate in the foot/heel junction and position before adhesive sets.
- For split toe feet install the heel wedge then remove a thin slice in the middle by cutting with a sharp knife through the split in the carbon foot module (**Figure 6**).

Instant adhesive is necessary to bond the urethane heel wedge. The adhesive cures in 15-20 seconds. For removal the adhesive may be softened by soaking in acetone or cyanoacrylate adhesive remover.

## HEEL ASSEMBLY

If it is necessary to replace the heel component or heel bolts use Loctite 410 and torque to the following specifications:

Category	Bolt Size	Heel Bolt Torque
1-6	1/4"	14 Nm (10ft-lbs)
7-9	5/16"	28 Nm (21ft-lbs)

**NOTE:** After service the bolts should be tightend, allowed to set for a few minutes and then tightened again. A biannual inspection of attachment bolts is recommended:

Check for corrosion and sign of wear. Replace bolts if needed.

## HEEL DIVIDER

The heel divider should be placed in the posterior of the splitted heel part. To secure its position a drop of instant adhesive can be used on one side.

## FLEX-FOOT SOCK

The toe end of the sock is divided by a thread and fits within the toe split. Upon fitting the sock, hold the top of the toe thread and slide the sock into the toe split. The thread of the sock should rest on top of foot module and not underneath the toe split. The sock should be loose from the heel area and secured proximally with the cable tie provided.

## SHOE HORN

When removing and installing the cover use shoehorn to avoid damage to foot or cover (**Figure 7 & 8**).


**NOTE:** Tampering with the pyramid attachment bolts will void the warranty of the product

## LIABILITY

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

## COMPLIANCE

This component has been tested according to ISO 10328 standard to two million load cycles. Depending on the amputee's activity this corresponds to a duration of use of two to three years. We recommend carrying out regular yearly safety checks.

In the standard mentioned, test levels (P) are assigned to a certain maximal body masses (m in kg). In some cases, which are marked with, no test level is assigned to the product related maximal body mass. In these cases, the test loads have been adapted adequately on the basis of the specified load level.

Category Össur high activity						
Category	Weight (Kg)	Lable text				
1	52	P3	ISO 10328 -	P3	52	kg
2	59	P3	ISO 10328 -	P3	59	kg
3	68	P4	ISO 10328 -	P4	68	kg
4	77	P4	ISO 10328 -	P4	77	kg

Category Össur high activity						
Category	Weight (Kg)	Lable text				
5	88	P4	ISO 10328 -	P5	88	kg
6	100	P5	ISO 10328 -	P5	100	kg
7	116	P5	ISO 10328 -	P6	116	kg
8	130	P6	ISO 10328 -	P6	130	kg
9	147	P6	ISO 10328 -	P6	147	kg

- This manual is intended for use by a certified prosthetist.

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**EN – Caution:** Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Össur products are used with other recommended Össur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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