



Instructions for Use





PRO-FLEX[®] LP

Product number: PLP0xyyz

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	10
FR	Notice d'utilisation	15
ES	Instrucciones para el uso	20
IT	Istruzioni per l'uso	25
NO	Bruksanvisning	30
DA	Brugsanvisning	35
SV	Bruksanvisning	40
EL	Οδηγίες χρήσης	45
FI	Käyttöohjeet	50
NL	Gebruiksaanwijzing	55
PT	Instruções de Utilização	60
JA	取扱説明書	65
KO	사용 설명서	70

ENGLISH

This manual is intended for use by a certified Prosthetist.

PRODUCT DESCRIPTION

Pro-Flex LP is a prosthetic foot device with integrated male pyramid and three composite blades. It features a sandal toe and a full length composite blade. The bottom blade is secured together with heel bolts.

The device consists of the following components (**Figure 1**):

1. Male Pyramid Adapter
2. Pyramid cover
3. Top blade
4. Middle blade
5. Bottom blade
6. Bottom blade bolts
7. Heel wedge kit
8. Spectra Sock
9. Footcover

INDICATIONS FOR USE

The foot is a non-invasive, single patient, reusable prosthetic device. It is made for walking and everyday use of low to high active users.

The Pro-Flex LP is designed for single patient use, everyday activities and not for extreme sports. More details see “technical specifications” and “use environment”.

MEDICAL INDICATION

Lower extremity amputations, for example due to:

- Trauma
- Vascular Disease
- Cancer
- Congenital Defects

The Pro-Flex LP is tested and validated for exoskeletal prosthetic fittings of lower extremities amputations.

TECHNICAL SPECIFICATIONS

Following specs are for a size 27 cat 5 foot:

- Product weight with footcover : 647 g / 1,42lbs
- Foot cover weight: 165 g / 0,36 lbs
- Build height: 73 mm / 2 7/8"
- Maximum user weight: 166 kg / 365 lbs

CATEGORY SELECTION CHART

Please refer to the selection charts below to determine the appropriate stiffness required according to Össur recommendations.

Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147	148-166
Weight lbs	99-115	116-130	131-150	151-170	171-194	195-221	222-256	257-275	288-324	325-365
Low Impact Level	1	1	2	3	4	5	6	7	8	9
Moderate Impact Level	1	2	3	4	5	6	7	8	9	N/A
High Impact Level	2	3	4	5	6	7	8	9	N/A	N/A

- Foot size 22–30

SAFETY PRECAUTIONS

Caution: Statement that alerts the user to the possibility of a problem with the device associated with its use or misuse.

Precaution: Should be taken to avoid the hazard

Note: Possibility of a technical damage or malfunction of the device.

Caution: Make sure that the user knows if he/she feels or notices any change of function or alignment of the device, he/she should contact his healthcare practitioner. This noticeable function loss may be the result of reduced stiffness or loss of support.

Note: The foot is designed to be used with a spectra sock and a foot cover.

BENCH ALIGNMENT (Figure 4)

- Fit foot with selected cover.
- Adjust to appropriate heel height (using the shoe).
- Introduce appropriate socket angles flexion/ extension and abduction/adduction.
- Divide the foot cover into 3 equal portions.
- The load line should fall at the junction of the posterior and the middle third portions as shown in figure 3.
- The 1/3 position is marked on the inside wall of the foot cover.

DYNAMIC ALIGNMENT

The heel stores energy after initial contact, slowly releasing it at midstance. The forward momentum generated by this action results in the toe being loaded for optimum energy release at terminal stance. The heel to toe action can be influenced by:

- A-P positioning of foot
- Dorsi-Plantarflexion
- Heel Stiffness

HEEL STIFFNESS

Heel too soft

Symptoms

- Foot will come to flat position too early (amputee feels he/she is sinking into a hole).
- Extra energy is required to climb up over the toe.
- Toe will feel too stiff
- Knee may hyperextend.

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Heel too hard

Symptoms

- Rapid heel to toe movement.
- At initial contact the amputee has poor control of his prosthesis.
- Minimal energy return feeling.
- Knee may become unstable.

Actions

- Replace heel wedge with a smaller one.
- Shift socket posterior (or foot anterior).

HEEL WEDGES

A wedge can influence the heel to toe function. Start altering the heel resistance to improve foot response by adding a heel wedge. The small, medium and large wedges can be interchanged to customize the stiffness and achieve the desired functional characteristics.

Temporary Wedge Placement (Figure 5)

- Cut the wedge to the width of the foot module
- Roughen the upper and lower surface of the wedge with abrasive paper.
- Place the wedge in the angle of heel and foot module
- Secure in position with tape wrapped around the foot module

Permanent Wedge Placement (Figure 6)

- Apply adhesive on the lower side of the wedge only.
- Locate in the foot/heel junction and position before adhesive sets.
- For split toe feet install the heel wedge then remove a thin slice in the middle by cutting with a sharp knife through the split in the carbon foot module.

Instant adhesive is necessary to bond the urethane heel wedge. The adhesive cures in 15-20 seconds.

PROTECTIVE FILM ON PYRAMID

The protective film protects the pyramid from scratching. Remove the film when the fitting is complete.

FOOT COVERS AND SPECTRA SOCKS (Figure 3)

The foot is designed to be used with foot covers and spectra socks. When installing and removing the foot covers use shoehorn to avoid damage to the foot or the foot cover.

Please view the instructions for use for the foot covers before fitting them and ensure that the user knows how to use them.

It is important that the middle section of the spectra sock covers the toplade to make sure it doesn't intervene with the movement of the blades.

The thread of the sock should rest on top of foot module and not underneath the toe split (**Figure 7**).

Use a shoehorn to avoid damaging to the foot or the foot cover.

1. Installing – Make sure that the foot module is placed completely into the foot cover.

2. Removing the foot from the foot cover.

Note: Use spectra socks as shown in figure 1.

MAINTAINANCE & CLEANING

Note: The device should always be used with spectra socks and foot covers to minimize dust and dirt from entering between carbon blades. It is important to fix the spectra sock around pyramid location in order to prevent noise issues coming from sand entering between the plates. The device is weather-proof but not resistant to corrosion. Therefore, the device should not come into contact with salt water or chlorinated water. The device is not designed to be used under extreme conditions like diving or jumping into water. If for some reason the foot components get wet parts should be dried with a lint free cloth. Should the device come into contact with salt water, chlorinated water, dust or dirt it must be rinsed with fresh water immediately.

DEBRIS BETWEEN BLADES

Noise could occur if sand or debris are present in between the top and middle blades, same could happen between bottom blade and middle blade.

Use compressed air to remove debris in between the carbon parts.

USE ENVIRONMENT

- Environmental operating temperature: -15 °C (+5 °F) to 50 °C (122 °F).
- Operating humidity: 10-95% relative humidity
- Shipping and storage temperature: -20 °C (-4 °F) to 70 °C (158 °F).

LIABILITY


The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

CE CONFORMITY

This device meets the requirements of the 93/42/EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guideline.

COMPLIANCE

This component has been tested according to ISO 10328 standard to two million load cycles. Depending on the amputee's activity this corresponds to a duration of use of two to three years. We recommend carrying out regular yearly safety checks.

In the standard mentioned, test levels (P) are assigned to a certain maximal body masses (m in kg). In some cases, which are marked with, no test level is assigned to the product related maximal body mass. In these cases, the test loads have been

adapted adequately on the basis of the specified load level.

ISO 10328 - label		
Category	Weight (Kg)	Label text
1	52	ISO 10328-P3-52-kg
2	59	ISO 10328-P3-59-kg
3	68	ISO 10328-P3-68-kg
4	77	ISO 10328-P4-77-kg
5	88	ISO 10328-P4-88-kg
6	100	ISO 10328-P5-100-kg
7	116	ISO 10328-P6-116-kg
8	130	ISO 10328-P6-130-kg
9	147	ISO 10328-P6-147-kg

WARRANTY

Össur provides a warranty of 36 months for this foot module and 6 months for the foot cover.

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Draft

EN – Caution: Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Össur products are used with other recommended Össur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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