Össur Americas
27051 Towne Centre Drive
Foothill Ranch, CA 92610, USA
Tel: +1 (949) 382 3883
Tel: +1 800 233 6263
ossurusa@ossur.com

Össur Canada
2150 – 6900 Graybar Road
Richmond, BC
V6W OA5 , Canada
Tel: +1 604 241 8152

Össur Europe BV
Ekkersrijt 4106-4114
P.O. Box 120
5690 AC Son en Breugel
The Netherlands
Tel: +800 3539 3668
Tel: +31 499 462840
info-europe@ossur.com

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










INDICATION FOR USE

- Anterior Cruciate ligament instability, posterior cruciate ligament instability, medial collateral ligament injury, lateral collateral ligament injury.
- Conditions requiring unicompartmental load control.

WARNING

This knee brace is a supportive brace only and is not intended or guaranteed to prevent knee injury.

Purchaser agrees to the following:

1. Össur accepts no responsibility for any injury which might occur while using Össur knee bracing products.
2. With use of this type of device, it is possible to receive various types of skin irritation, including rashes and pressure points. If you experience any unusual reactions while using this product, consult your health care provider immediately.
3. If you have any questions as to how this knee brace applies to your own specific knee condition, consult your health care provider.

If you do not agree to the above conditions, contact your brace provider before using this product.

BRACE APPLICATION

1. Unfasten all straps of the brace.
2. With the knee bent in approximately 90 degrees, place the brace over your knee such that the patella is centered between the condyle pads. **(Figure 2)**
Note: This places the hinge slightly above the joint line so the brace can properly seat itself when you stand. **(Figure 3)**
Note: All straps can be trimmed or adjusted to length.
3. Fasten the lower cruciate strap just below the knee. Fasten the strap snugly above the calf muscle as this will aid in the brace suspension. **(Figure 4)**
4. Next, snugly fasten the bottom strap by inserting it through the D-ring.
5. Next, snugly fasten the upper cruciate strap above the knee.
6. Next, snugly fasten the top strap. **Note:** Over-tightening the two top straps can cause the brace to migrate.
7. Last, snugly fasten the anterior tibial strap **(Figure 5)**. If the brace is too anterior, loosen the anterior tibial strap and tighten the lower

EN – Caution: This product has been designed and tested based on single patient usage and is not recommended for multiple patient use. If any problems occur with the use of this product, immediately contact your medical professional.

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cruciate and bottom straps. If the frame is too posterior, loosen the lower cruciate and bottom straps AND tighten the anterior tibial strap. Proper anterior/posterior placement will help prevent migration and ensure proper brace suspension.

Warning: If brace causes discomfort or irritation, discontinue use immediately and contact your health care provider.

BRACE MAINTENANCE

The Rebound DUAL brace is designed for minimum maintenance and care. Follow these easy maintenance steps to keep your brace in good condition.

CLEANING

- Liners, straps and pads may be removed for cleaning
- Do NOT machine wash
- Do NOT place in dryer
- Hand-wash using mild non-detergent soap
- Rinse thoroughly
- Air-dry

If migration occurs, clean lower liners and pads with an alcohol wipe or damp cloth. If brace is used in salt water or chlorinated water, rinse well and air dry.

HINGE MAINTENANCE

- Do NOT disassemble the hinge
- Check regularly for foreign materials in hinge (i. e. dirt or grass)
- Carefully clean dirt from hinge using water or compressed air
- Lubricate with silicone spray or dry graphite if needed

ACCESSORIES

The following accessories are available for your Rebound DUAL brace at an additional charge.

To order call your Össur Customer Service Department.

- SmartDosing kit
- Flexion stop kit
- Liner kit
- Strap kit
- D-ring kit
- Extension Stop Kit

SETTING UP THE D-RING

The field-serviceable d-rings can be changed from internal to external d-rings or replaced if a d-ring is broken.

1. While holding the d-ring firmly against the frame, remove the screw (**Figure 6**). The following instructions show replacement of internal d-ring; the d-rings can also be replaced or reconfigured to be external.
2. Remove the d-ring system from the frame.
3. Snap the female insert into the d-ring and place the d-ring/female insert assembly into the frame (internal or external). (**Figure 7**)
4. Insert new male insert on the opposite side of the frame into the d-ring assembly. (**Figure 8**)

Warning: Once the D-ring insert has been used it cannot be used

again. Do NOT re-use the same insert to replace a broken D-ring or to change from internal to external d-rings. **Note:** Replace with the appropriate male insert. The outer surface of the insert has either no ring or one ring. Use the same type as the originally installed version.

5. Tighten down the screw.

REFURBISHING

From time to time, your Rebound DUAL should be sent to Össur for refurbishing. The frequency of this will depend on how often the brace is worn and for what activities. Your brace supplier should help you to determine when your brace needs to be sent in for service, but once a year is a good rule of thumb.

As part of the refurbishment process, Össur will check your brace for signs of wear and rebuild it with all new components, including rivets, buckles, straps and padding. There is a usual turnaround time of three to five working days. If you have any questions, contact your local Customer Service Department.

FIT ADJUSTMENTS

Contouring the Brace:

The Rebound DUAL brace can be adjusted for individual anatomies. The aircraft grade aluminum frame can be cold formed to accommodate different thigh and tibia shapes using bending irons as needed.

Adjusting the Varus / Valgus angle of the Brace:

The Rebound DUAL brace can be adjusted for individual varus or valgus angles and to create an offloading force to either the medial or lateral compartment. To adjust the angle, loosen the screws on the desired side of the brace and slide the upper frame to the desired position with a steady force. Adjust the medial upper frame to offload the medial compartment or adjust the lateral upper frame to offload the lateral compartment. The slide bar has a lettered scale and may be moved to increase or decrease the amount of varus or valgus force on the leg. **(Figure 9)** The lettered scale does not relate to degrees of offloading and should be used as a guide. Retighten the screws until snug (recommended torque value of 3 in-lb). All adjustments should be determined by a healthcare provider, in accordance with the patient's comfort level. Each time the brace is adjusted, the healthcare provider should note the letter for a record of progression.

Adjusting the Height of the Brace:

The Rebound DUAL brace can be adjusted for individual height requirements. The brace height can be adjusted by loosening the screws on both sides of the upper frame and pulling the frame up with a steady force. Ensure the brace is properly aligned by checking that both hinge arms are located on the same position on the medial and lateral hinge arm. Retighten the screws until snug.

Widening or narrowing fit at knee joint:

The Rebound DUAL brace can be adjusted for individual knee width. To widen or narrow the width at the knee, the medium condyle pads can be removed and replaced with the thick or thin version shipped with the Rebound DUAL. Any combination of pads can be used to achieve a more comfortable fit.

Adjusting fit at the tibia:

The tibial pad can be installed to fill in any gaps around the tibia and achieve the desired fit. To add the tibial pad, pull back the lower liner and place the tibial pad on top of the hook-and-loop fastener to which the lower liner was originally fastened. Press the lower liner back into place on top of the pad. Adjust the padding until the Rebound DUAL brace fits snugly against the tibia.

Changing Extension Stops

Extension stops are available for varying degrees of extension control.

1. Remove the pad from the condylar area.
2. Use a short Philips head screwdriver and remove the retaining screw.
3. Slide the 0° stop out of the hinge.
4. Insert desired extension stop.
5. Secure with retaining screw (DO NOT OVER TIGHTEN)
6. Repeat on the other hinge.
7. To add flexion stop: follow steps 1, insert desired flexion stop, and then follow steps 5–6.

LIMITED WARRANTY

This warranty covers all defects in materials and workmanship of your new Rebound DUAL Knee Brace. This warranty covers the frame and hinges for 2 years, and custom warranty is lifetime. The warranty covers all other materials, including straps, sleeves, liners and pads for six months from the date of original purchase.

Össur will repair or replace (at its option) any defective part covered under warranty at no charge. Repairs necessitated by accident, negligence, misuse or abuse by the purchaser are not covered. Also, consequential and incidental damages are not recoverable under this warranty. These limitations may not apply as some states do not allow the exclusion or limitation of incidental or consequential damages. Only the original purchaser is entitled to claim service under this warranty.

To obtain service under this warranty call our Customer Service Department.

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