A CONTENIDO (Dib.1)

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


A PACK CONTENTS (Fig.1)

- a **1K150** double female pyramid connector including :
 - a M36 threaded extension piece with 1.5 pitch (1) and its screw fittings
 - a M36 tapped sleeve with 1.5 pitch (2) and its screw fittings
- 4 headless hexagonal socket screws M8 x 20 (not shown) [See Use §E => Fig. 5b]
- this leaflet of instructions for use.

B USE AND CONNECTIONS (Fig.2)

- The **1K150** adjustable double connector is convenient for patients up to 100 kg, carried load included.
- The **1K150** connector should be installed between two male pyramid components (3 and 4).

C WHICH LENGTH SHOULD WE CUT OUT (Fig.3)

- Determine the build height of the product according to its location in the prosthesis (in the following operations this dimension will be called Y).
- Screw the threaded extension piece (1) into the tapped sleeve (2) until you reach the thread root. Measure dimension X.
- The length to be cut will be calculated as follows : **Length to be cut = X-Y.**

D CUTTING TO LENGTH THE THREADED EXTENSION PIECE (Fig.4)

Using the 5C413 cutting tool is highly recommended to cut the threaded extension piece to length.

Caution : The cutting tool is 12 mm thick. This thickness indicates the minimal length of the re-cut threaded piece. NEVER re-cut the threaded piece so that it becomes shorter than 12 mm because this would lead the component to get fragile.

- Screw the threaded extension piece (1) into the **5C413** cutting tool so that the extension piece protrudes out side **F** of the tool by the length that should be cut.
- Install the assembly «threaded extension piece and cutting tool» into a vice and firmly tighten the jaws to lock the threaded extension piece.
- Use a hacksaw to cut the threaded extension piece (1) along side **F** of the **5C413** cutting tool.
- Loosen the jaws of the vice and remove the assembly «threaded extension piece and cutting tool» from the vice.
- Unscrew the threaded extension piece and remove it from the cutting tool.
- Carefully burr the threaded extension piece at the section level.

E ASSEMBLING THE 1K150 CONNECTOR ONTO THE PROTHESIS (Fig.5a / 5b)

- Re-screw the threaded extension piece (1) into the tapped sleeve (2) until you reach the thread root.
- Install the **1K150** connector onto the prosthesis.
- Check that the (2a) screws are in the sagittal plane and the (2b) screws in the frontal plane. If there are not so, you can slightly unscrew the tapped sleeve (2) so that the (2a) and (2b) screws be well positioned.




Caution : Never unscrew the tapped sleeve by more than a quarter turn. Otherwise the component would get fragile.

- Tighten the (2c) screw of the sleeve up to 7.4 Nm.
- Proceed to the alignment adjustments.
- If one of the screws (1a or 2a or 2b) is too short, replace it by one of the 20 mm long screws included in the pack (fig. 5b).
- Cover the pyramid locking screws (1a / 2a / 2b) with XC047 adhesive and tighten them up to a 15 Nm tightening torque.

F RECYCLING

The (1) and (2) components are made of aluminum alloy and the (1a), (2a), (2b), (2c) screws are made of steel. Each of these components should be recycled according to the laws in force.

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