













Instructions for Use

REBOUND[®] DUAL ST

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	9
FR	Notice d'utilisation	13
ES	Instrucciones para el uso	17
IT	Istruzioni per l'uso	21
NO	Bruksanvisning	25
DA	Brugsanvisning	29
SV	Bruksanvisning	33
EL	Οδηγίες Χρήσης	37
FI	Käyttöohjeet	41
NL	Gebruiksaanwijzing	45
PT	Instruções de Utilização	49
PL	Instrukcja użytkownika	53
CS	Návod k použití	57
TR	Kullanım Talimatları	61
RU	Инструкция по использованию	65
JA	取扱説明書	69
ZH	中文说明书	73
KO	사용 설명서	77

PRODUCT DESCRIPTION

List of Parts (Figure 1):

1. Upper Thigh Strap #2
2. Lower Thigh Strap #3
3. Calf Strap #5
4. Upper Shin Strap #4
5. Lower Shin Strap #1
6. Upper Frame
7. Lower Frame
8. D ring
9. Condyle pad
10. PCL Slot

INDICATIONS FOR USE

The brace is indicated for combined and non-combined instabilities of the:

- Anterior cruciate ligament (ACL)
- Medial collateral ligament (MCL)
- Lateral collateral ligament (LCL)
- Posterior cruciate ligament (PCL)

CONTRAINDICATIONS

None known.

GENERAL SAFETY INSTRUCTIONS

This knee brace is a supportive knee brace only. It is not intended or guaranteed to prevent knee injury.

If brace causes discomfort, irritation, or is damaged, discontinue use immediately and contact your health care provider.

The purchaser agrees to the following:

1. Össur accepts no responsibility for any injury which might occur while using Össur knee bracing products.
2. With use of this type of device, it is possible to receive various types of skin irritation, including rashes and pressure points. If you experience any unusual reactions while using this product, consult your health care provider immediately.
3. If you have any questions as to how this knee brace applies to your own specific knee condition, consult your health care provider.
4. If brace causes discomfort, irritation, or is damaged, discontinue use immediately and contact your health care provider.
5. This brace is single patient use only. If you do not agree to the above conditions, contact your brace provider before using this product.

FITTING INSTRUCTIONS

Brace Application

1. Unfasten straps 1#,2#,3# and 4#. Loosen strap #5 (Figure 2).
2. Bend your knee to between 30–45 degrees, step through the brace, and then slide it up your leg. Align middle of the hinge with top of the knee cap. Make sure the condyle pads are in full contact with the knee (Figure 3).

3. Insert lower shin strap #1 through the D ring and fasten tightly (**Figure 4**).
4. Insert upper thigh strap #2 through the D ring and fasten (**Figure 5**).
5. Then fasten lower thigh strap #3. Tightening the lower thigh strap #3 will pull the hinges towards the back of the knee.

NOTE: Do not overtighten the two top straps as this can make the brace move. (**Figure 6**).

6. Insert upper shin strap #4 through the D ring and tighten until the hinges are aligned with the midline of the knee. Readjust strap #3 if necessary (**Figure 7**).
7. Tightly fasten calf strap #5, while making sure it fits comfortably. After fitting, the calf strap will only need to be readjusted to improve brace comfort/ suspension (**Figure 8**).
8. When standing, the brace should fit the leg with no larger gaps. The hinge should be aligned with the middle of the knee cap. Seen from the side, the hinges are aligned with the knee midline (**Figure 9**).

Troubleshooting:

If the brace is too far in the front, loosen the upper shin strap #4, and then tighten the calf strap #5 and lower thigh strap #3.

If the frame is too posterior, loosen the calf #5 and lower thigh #3 straps and tighten the upper shin strap #4.

BRACE ADJUSTMENTS

Adjusting Pad and Strap Length:

- For optimal fit and comfort, ensure the strap pads are centered on the leg when tightening straps.
- If a strap pad is too long, remove the pad from the strap and use scissors to shorten the pad.
- To adjust the length of a strap, detach the alligator clip at the end of the strap.
- Adjust the length as necessary and reattach the alligator tab to the strap.

Adjusting the Brace Contouring:

- The Rebound DUAL ST brace can be adjusted to individual anatomies.
- Use bending irons to cold form the aircraft grade aluminum frame to fit different thigh and calf shapes.

Widening or narrowing fit at knee joint:

- The Rebound DUAL ST brace can be adjusted to individual knee widths.
- Adjust to knee width by removing the medium condyle pads and replacing them with the thick or thin version
- Any combination of pads can be used to achieve a more comfortable fit.

Replacing and Reconfiguring D rings:

D rings are field-serviceable. D rings can be moved from the interior to the exterior or replaced if broken.

To move a D ring from the interior to the exterior:

1. Remove strap from D ring.
2. Bend D ring mounted on the brace interior and push it inwards to release from the keyhole lock. Remove from the rivet (**Figure 10**).
3. Attach D ring to the brace exterior by pushing down onto rivet through the keyhole lock. Slide and click D ring into place. Place keyhole lock with indented side facing up (**Figure 11**).
4. Replace strap. If needed, repeat steps for other D rings.

Range of Motion adjustment (ROM):

The brace has 0° extension stops attached to the hinges. If needed, range of motion can be adjusted using the extension and flexion stop kits. Extension and flexion stop kits are available separately for order.

To install a kit:

1. Remove the condyle pads from the medial and lateral hinges. Remove retaining screws.
2. Using a short Philips head screwdriver, slide the 0° stop out of the hinge on both sides.
3. Re-apply the desired extension stop and gently secure it with the retaining screw.
4. Repeat steps to change the flexion stop.

NOTE: Do not overtighten retaining screw.

PCL Strap Attachment:

The PCL kit is designed to counteract posterior subluxation of the tibia. Trim and adjust the PCL strap for proper fit.

To attach the PCL strap:

1. Take the end of the PCL strap with the alligator clip and loop it through the PCL slot on the brace frame. Alligator clip will close on the inside of the strap. Trim strap if necessary (**Figure 12**).
2. Loop the other end of the PCL strap through the frame. Fasten strap and secure it onto the outside of the strap (**Figure 13**).
3. Place the adjustable pad in the middle of the strap. Trim the strap pad if needed.
4. Fully extend the leg. Tighten the strap, securing the pad.

PRODUCT MATERIAL LIST

Aluminum, stainless steel, polycarbonate, polyamide, silicone, thermoplastic elastomer, polyester, acetal, glass-filled nylon, low density polyethylene, ethylene vinyl acetate, spandex, polyurethane, polyethylene.

CLEANING & MAINTENANCE

Follow the steps below to keep your brace in a good condition.

CLEANING

Liners, straps and pads can be removed for cleaning. Hand-wash using mild non-detergent soap.

Rinse thoroughly

Air-dry

If the brace moves, clean lower liners and pads with an alcohol wipe or damp cloth.

NOTES:

- Do not machine wash
- Do not tumble dry

Hinge Maintenance

Regularly check the hinge regularly for dirt, debris or similar contamination. Remove dirt, debris or similar contamination carefully using water or compress air.

NOTE: Do not disassemble the hinge.

DISPOSAL

All components of the product and packaging should be disposed of in accordance with respective national environmental regulations. Users should contact their local governmental office for information on how these items can be disposed of in an environmentally sound manner.

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