
3R60=P

DE	Gebrauchsanweisung	12
EN	Instructions for use	22
FR	Instructions d'utilisation	33
IT	Istruzioni per l'uso	44
ES	Instrucciones de uso	55
PT	Instruções de utilização	66
NL	Gebruiksaanwijzing	77
SV	Bruksanvisning	88
DA	Brugsanvisning	98
NO	Bruksanvisning	108
FI	Käyttöohje	119
PL	Instrukcja użytkowania	129
HU	Használati útmutató	140
CS	Návod k použití	150
RO	Instrucţiuni de utilizare	161
HR	Upute za uporabu	171
TR	Kullanım kılavuzu	182
EL	Οδηγίες χρήσης	192
RU	Руководство по применению	204
UK	Інструкція щодо використання	216
JA	取扱説明書	228
ZH	使用说明书	238

1

1 3R60=P

2 2R38


3 2Z11=KIT


4 647G1868

2


3

743Y580=*
(3R80, 3R106=HD, =KD, =ST)

743A220
PROS.A. Assembly


4

704Y14=*


719R5



2


718R1

3


7


8


9

3D-Mode



10


2D-Mode


11



12



13


14








15

17

2

< 5 mm

3

≥ 5 mm

Product combinations	
Reference number	Combinable
3S27 (for hip disarticulation)	✓
3S107	✓
7E7	✓
7E9	✓
7E10	✗

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INFORMATION

Date of last update: 2026-07-28

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.

- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

The modular EBS knee joint 3R60=P has the following key features:

- Polycentric prosthetic knee joint
- Stance phase control
 - Polycentric kinematics
 - Adjustable ergonomically balanced stride for limited knee flexion from 0 – 10° at heel strike
- Swing phase control
 - Adjustable swing phase flexion resistance (hydraulic damping)
 - Adjustable swing phase extension resistance (hydraulic damping)
 - Hydraulic unit with extension assist function (spring force)
- Large flexion angle
 - Up to a maximum of 165°

2 Safety-related information

CAUTION!

Risk of injury and risk of product damage

- ▶ To prevent mechanical damage, treat the product with care.
- ▶ Observe the information on product combinations in the instructions for use for the products.
- ▶ Observe the lifetime of the product (see page 24).
- ▶ Observe the product's operating conditions (see page 23).
- ▶ Do not expose the product to restricted environmental conditions (see page 24).
- ▶ Check the product for damage and readiness for use prior to each use.
- ▶ Do not reach into the joint mechanism. There is a risk of pinching.
- ▶ Parts of the product may heat up during the activity. Be careful when touching the product. There is a risk of burns.
- ▶ If there are changes in functionality, reduce all activities so the product can cool down.
- ▶ Do not continue using the product if oil comes out. Have the product repaired.
- ▶ Follow the maintenance instructions (see page 31).
- ▶ Observe the instructions for cleaning and care (see page 30).
- ▶ Do not use the product if it is damaged or if you are not sure whether it is working properly. Have the product inspected, repaired or replaced by the manufacturer or a specialist workshop.
- ▶ The product is intended for use by one person only and may not be reused by other persons.

3 Intended use

3.1 Indications for use

The product is to be used **solely** for lower limb exoprosthetic fittings.

3.2 Operating conditions

- Maximum permissible body weight: **125 kg**

The product is recommended for:

- Mobility grade 2: restricted outdoor walker
- Mobility grade 3: unrestricted outdoor walker

3.3 Environmental conditions

Storage and transport

Storage temperature: -20 °C to +60 °C, relative humidity: 20 % to 90 %, no vibrations or impacts

INFORMATION

Clean the product after contact with humidity/chemicals/solids, in order to avoid increased wear and tear and damage (see page 30).

Allowable environmental conditions

- Temperature range: -10 °C to +45 °C
- Relative humidity: 20 % to 90 %
- Contact with soap solution (e.g. when taking a shower)
- Contact with perspiration
- Submersion in fresh water/chlorinated water (e.g. swimming pool):
 - Maximum immersion depth: 2 m
- Submersion in salt water with a low salt content (14 days/year, 1 hour/day):
 - Permissible salt content: max. 3.5 % (e.g. ocean)
 - Maximum immersion depth swimming pool: 1 m
 - Maximum immersion depth sea: 0.5 m
- Contact with dust, drifting sand (e.g. walking on the beach)
- Contact with salty, condensing air
- UV resistant

Unallowable environmental conditions

- Contact with liquid-binding particles (e.g. talcum)
- Contact with large amounts of sand and dust (e.g. digging, kneeling in sand, construction site work)
- Contact with acids
- Contact with urine
- Submersion in salt water with excessive salt content
- Salt content: > 3.5 % (e.g. brine bath, Dead Sea)
- Cleaning agents and disinfectants with solvents, chlorine and phosphate
- High water pressure (e.g. diving, jumping into water)

3.4 Product combinations

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

Permissible and impermissible product combinations with the 3R60=P prosthetic component are listed in the tables (see page 10).

3.5 Lifetime

This product was tested by the manufacturer with 3 million load cycles. Depending on the user's activity level, this corresponds to a maximum lifetime of 5 years.

4 Scope of delivery

The section is intended for qualified personnel.

3R60=P (see fig. 1)		
Item	Reference number	Designation
①	3R60=P	Prosthetic knee joint
②	2R38	Tube adapter
③	2Z11=KIT	Adapter protector
④	647G1868	Instructions for use

5 Accessories

The section is intended for qualified personnel.

Reference number	Designation
4Z59=6.3X6	Stop bumper
4G2014=1	Flexion stops to limit maximum flexion to 135°

6 Preparing the product for use

The section is intended for qualified personnel.

CAUTION

Incorrect alignment, assembly or adjustment

Injury due to incorrectly installed or adjusted as well as damaged prosthetic components

- Observe the alignment, assembly and adjustment instructions.

6.1 Optional: Using the adapter protector

The adapter protector protects the connecting section of the prosthetic joint from scratches during alignment and fitting.

- > The instructions refer to the following illustration: see fig. 2
- 1) Follow the instructions in the accompanying document for the adapter protector.
- 2) Remove the adapter protector from the product before handing it over to the user.
- 3) Readjust the prosthetic alignment after removing the adapter protector.

6.2 Optional: Using the 743A220 PROS.A.* Assembly alignment apparatus

An alignment apparatus such as the 743A220 PROS.A. Assembly can be used for the alignment of modular leg prostheses.

- > The instructions refer to the following illustration: see fig. 3
- 1) Follow the instructions in the accompanying document for the alignment device.
- 2) Mount the 743Y580=* adapter inserts with the labels, **3R80**, **3R106=HD**, **=KD**, **=SD**.
- 3) Carry out the alignment of the prosthesis.

6.3 Shortening the tube adapter

CAUTION

Incorrect processing of the tube

Risk of injury due to damage to the tube

- Do not clamp the tube in a vice.
- To shorten the tube, use only a tube cutter or a cutting device.

- > The instructions refer to the following illustration: see fig. 4
- > **Required materials and tools:** 704Y14=* cutting device or 719R5 tube cutter, 718R1 tube deburrer
- 1) Shorten the tube adapter with a cutting device or a tube cutter ①.

- 2) Deburr ❸ the inside and outside of the cut edge with a tube deburrer ❷.

6.4 Installing the tube adapter

- > The instructions refer to the following illustration: see fig. 5
 - > **Required materials and tools:** 634A58 isopropyl alcohol, 710D20 torque wrench, hexagon socket screw adapter 4 mm
- 1) Clean the contact surfaces of the tube and of the tube receiver using a degreasing cleaning agent (such as isopropyl alcohol).
 - 2) **CAUTION! Insert the tube adapter into the tube receiver up to the stop. Otherwise, components may become overloaded and break.**
Insert the tube adapter into the tube receiver up to the stop ❶.
 - 3) Tighten the cap screw **without** additional Loctite to **10 Nm** ❷.

6.5 Bench Alignment

- Align the prosthesis using the alignment images (see fig. 7; see fig. 8) as well as the instructions for use of all prosthetic components used.
 - The prosthesis should be aligned so as to enable the user to stand safely.
 - > **Required materials and tools:** alignment apparatus (e.g. 743A220 PROS.A. Assembly), 743A80 50:50 gauge, 743S12 heel height measuring device
- 1) Position the prosthetic foot in the alignment apparatus.
 - 2) **Heel height**
→ The effective heel height of everyday shoe (x) **(+ 5 mm)**.
 - 3) **Frontal plane (see fig. 7)**
→ Align the prosthetic foot with the alignment reference line: **Position the middle of the prosthetic foot on the alignment reference line (external rotation approx. 5 – 7°).**
 - 4) **Sagittal plane (see fig. 8)**
→ Align the prosthetic foot with the alignment reference line: **Place the alignment reference line behind the middle of the prosthetic foot (see table in figure).**
 - 5) Position the prosthetic knee joint and prosthetic socket in the alignment apparatus.
 - 6) **Frontal plane (see fig. 7)**
→ Align the prosthetic socket with the alignment reference line.
→ Take the individual adduction position into account.
 - 7) **Sagittal plane (see fig. 8)**
→ Align the prosthetic socket with the alignment reference line.
→ Determine the centre of the prosthetic socket proximally and distally with the 50:50 gauge in the direction of movement and mark the centre line.
→ Mark the socket reference point on the centre line proximal to the tuberosity reference point **(30 mm)**.
→ Align the prosthetic socket so the alignment reference line runs through the socket reference point.
→ Adjust the socket flexion by rotating around the socket reference point: **individual residual limb flexion/flexion contracture (x) (+ 5°).**

6.6 Static Alignment

Ottobock recommends optimising the alignment with the 743L500 3D L.A.S.A.R. Posture measuring device. The 743L100 L.A.S.A.R. Posture can be used as an option.

- The following instructions refer to the illustrations: see fig. 9, see fig. 10

Optimising the prosthetic alignment:

Make sure that the user is standing comfortably and that the reference points in the illustrations are observed.

Noting the lines:

Dotted lines = auxiliary lines

Solid lines = load lines

Maintaining distances:

Use the distances between the auxiliary line and load line shown in the illustrations for static prosthetic alignment.

3D mode (see fig. 9)

- ❶ **Prerequisite:** The prosthetic knee joint is fully extended.
- ❷ The auxiliary line lies on the sagittal reference point of the prosthetic knee joint (front bottom pivot point).
- ❸ The auxiliary line lies on the frontal reference point (anterior superior iliac spine).
- ❹ The auxiliary line lies on the frontal reference point of the prosthetic knee joint (pyramid).

2D mode (see fig. 10)

- ❶ **Prerequisite:** The prosthetic knee joint is fully extended.
- ❷ The auxiliary line lies on the sagittal reference point of the prosthetic knee joint (front bottom pivot point).
- ❸ The auxiliary line lies on the frontal reference point (anterior superior iliac spine).
- ❹ The auxiliary line lies on the frontal reference point of the prosthetic knee joint (pyramid).
- ❺ The auxiliary line lies on the frontal reference point of the prosthetic foot (pyramid).

6.7 Dynamic Trial Fitting

The optimal gait pattern is developed during dynamic fitting. The prosthesis is aligned optimally in the frontal plane and the sagittal plane for this purpose.

6.7.1 Checking the settings

- 1) Manually check the settings of the product prior to the dynamic fitting.
- 2) In case of deviations, reset the product to the factory settings.

6.7.2 Factory settings

Function	Adjustment position	Factory settings	Meaning
Ergonomically balanced stride	Adjustment wheel for EBS unit (see fig. 11)	Clockwise, to the stop, then pre-tension the bumper anticlockwise by 2 turns	Minimum resistance
Swing phase flexion resistance	Adjustment screw F (see fig. 6)	Middle setting (12 o'clock)	Medium resistance
Swing phase extension resistance	Adjustment screw E (see fig. 13)	Anticlockwise, to the stop	Minimum resistance

6.7.3 Adjusting the ergonomically balanced stride (EBS)

INFORMATION

The user should do exercises on the parallel bars to learn how the ergonomically balanced stride unit works. Additional exercises on level ground and on the ramp help to optimise the adjustment of the prosthesis and find a good compromise between both situations.

The ergonomically balanced stride unit produces controlled stance phase movement in the range between 0° and 10°.

> The instructions refer to the following illustration: see fig. 11

> **Required materials and tools:** hexagon socket wrench 4 mm

► Adjust the EBS pre-tension so that the prosthetic knee joint goes into stance phase flexion at heel strike, not too fast and not too slow.

→ **For more stance phase flexion**, turn the adjustment screw clockwise ❶.

→ **For less stance phase flexion**, turn the adjustment screw anticlockwise ②.

Adjustment		Meaning	
Turning clockwise	=	Pre-tension is reduced	More stance phase flexion
Turning anticlockwise		Pre-tension is increased	Less stance phase flexion

6.7.4 Swing phase flexion resistance

> The instructions refer to the following illustration: see fig. 12; see fig. 14

> **Required materials and tools:** hexagon socket wrench 2 mm

► Adjust the swing phase flexion resistance so that the prosthetic foot swings through sufficiently at a slower walking speed and does not swing through too far at a faster walking speed.

→ **For more resistance**, turn adjustment screw **F** clockwise.

→ **For less resistance**, turn adjustment screw **F** anticlockwise.

6.7.5 Swing phase extension resistance

CAUTION! Set the swing phase extension resistance so that the prosthetic knee joint reaches full extension, even when taking slow and small steps

> The instructions refer to the following illustrations: see fig. 13; see fig. 14

> **Required materials and tools:** hexagon socket wrench 2 mm

► Adjust the swing phase extension resistance so the prosthetic knee joint reaches full extension but does not swing too hard against the extension stop.

→ **For more resistance**, turn adjustment screw **E** clockwise.

→ **For less resistance**, turn adjustment screw **E** anticlockwise.

6.8 Using flexion stops

- Use flexion stops with protruding prosthetic components (e.g. Quickchange adapter...) if they collide with the prosthetic knee joint at maximum flexion (see fig. 6).
- If there is no risk of collision, no flexion stops are needed.

Only for hip disarticulation:

- Using flexion stops is not permitted.
- The 4R57 adapter may be used even if it touches the prosthetic knee joint

Installing flexion stops

> The instructions refer to the following illustration: see fig. 16

> **Required materials and tools:** flexion stop set (4G2014=1), flat-head screwdriver, 2 mm hexagon key, Loctite 241 (636K13)

1) Remove the cover caps on both sides ① of the product:

→ **NOTICE! When prying off the cover caps, take care not to damage the protective stickers beneath them with the screwdriver.**

→ Carefully insert the flat-head screwdriver into both openings in the cover cap, one after the other ②.

→ Carefully pry off the cover cap ③, only levering the screwdriver towards the knee joint.

→ Repeat these steps on the opposite cover cap.

2) Position the flexion stops in the product ④.

3) Apply Loctite 241 ⑤ to the threads of the screws supplied.

4) Use the hexagon key ⑥ to hand-tighten the coated screws.

5) Attach the cover caps on both sides ⑦ of the product.

→ Installation of the flexion stops is complete.

→ After assembly of the flexion stops, the prosthetic knee joint has a reduced **flexion angle of 135°**.

6.9 Optional: Installing the foam cover

CAUTION

Use of liquid-binding substances (e.g. talcum)

Risk of injury, damage to the product due to lack of lubrication

- ▶ Prevent contact of the product with liquid-binding substances.

- > **Prerequisite:** Use the **3S107** foam cover for the product.
- > **Required materials and tools:** 519L5 silicone spray
- 1) Spray the silicone spray directly onto the contact surfaces in the foam cover to improve the gliding characteristics and reduce noise.
- 2) Attach the foam cover to the prosthesis.
 - The foam cover must not impede full knee extension at terminal swing.
 - After attaching the foam cover, perform another **dynamic fitting** (see page 27) to check the settings.
- 3) Finish the prosthesis.

6.10 Finishing the prosthesis

Finish working on the prosthesis using the reference image (see fig. 15) and the following instructions.

- 1) Check that the joint collides with the socket or prosthetic components only at the tube clamp. If not, install flexion stops (see page 28).
 - 2) Remove the adapter protector from the product before handing it over to the user.
 - 3) Readjust the prosthetic alignment after removing the adapter protector.
 - 4) Replace set screws that are too short with longer ones to prevent components from coming loose accidentally.
 - 5) Replace set screws that are too long with shorter ones to prevent collisions and product damage.
 - 6) Clean the threads of the screw connections on the prosthetic socket and foot using a degreasing cleaner.
 - 7) Only use Loctite 241 (636K13) on the screw connections of the prosthetic socket and foot.
 - 8) Finish all screw connections. Observe the tightening torques and thread locks for all prosthetic components.
- The prosthesis is completed.

7 Use

The section is intended for users.

7.1 Information for use

INFORMATION

Use in salt water

Using the prosthetic knee joint in salt water constitutes an increased load. Despite thorough cleaning and care, signs of corrosion may occur (e.g. on set screws) because invisible salt residues are present. As long as the function of the product is not impaired, it does not have to be sent in for service.

- ▶ Clean the product thoroughly after use in salt water (see page 30).
- ▶ Contact qualified personnel if the functionality of the product is impaired.

INFORMATION

Temporary change in walking sensation after kneeling down

The knee joint may feel somewhat softer for a short while when walking after kneeling down and standing up. This is normal and disappears within 10 minutes.

- ▶ **Prior to use**, check that **full flexion**, **extension** and **damping resistance** are working.
- In case of insufficient functionality, have the prosthesis adjusted by qualified personnel.

7.2 Practising sitting down

- > Practise sitting down with the user so they can become familiar with the behaviour of the joint.
- > Have the user complete the following steps:
 - 1) Hold on to a suitable object such as the chair arm supports with both hands.
 - 2) Shift your weight to the sound side.
 - 3) Place weight on the **forefoot** of the prosthesis to let the prosthetic knee joint flex.
 - 4) Sit down carefully.

7.3 Walking on stairs or ramps

⚠ CAUTION

Walking step-over-step on stairs or ramps not possible

Risk of falling due to unintentional buckling of the prosthetic knee joint

- ▶ Make sure that the prosthetic knee joint is fully extended prior to heel strike.

Walking step-over-step on stairs or ramps is not possible because the prosthetic knee joint flexes in an uncontrolled manner when it is flexed and the heel makes contact at the same time.

The prosthetic knee joint provides optimum stability only if fully extended prior to loading response.

Walking down stairs or ramps

- > **Prerequisite:** Practise walking correctly on stairs or ramps.
 - 1) Extend the prosthetic knee joint in front of the stairs or ramp.
 - 2) Hold on to the railing for safety.
 - 3) First set the prosthetic leg down in the extended state.
 - 4) Repeat with the sound leg.
 - 5) Maintain the correct sequence until the end.
- The user should be able to walk safely on stairs and ramps.

Walking up stairs or ramps

- > **Prerequisite:** Practise walking correctly on stairs or ramps.
 - 1) Extend the prosthetic knee joint in front of the stairs or ramp.
 - 2) Hold on to the railing for safety.
 - 3) Take the first step with the sound leg.
 - 4) Then continue with the prosthetic leg in the extended state.
 - 5) Maintain the correct sequence until the end.
- The user should be able to walk safely on stairs and ramps.

8 Cleaning

The section is intended for users.

Requirements for cleaning agents and disinfectants

- Solvent-free
- Chlorine-free
- Phosphate-free

Cleaned using soap and water

- > **Required materials and tools:** 2 microfibre cloths, pH-neutral soap (e.g. 453H10 Derma Clean), tap water
- 1) **In case of light soiling (e.g. small deposits):** Moisten a microfibre cloth with water.
- 2) **In case of heavy soiling (e.g. stubborn deposits):** Moisten a microfibre cloth with pH-neutral soap and water.
- 3) Clean the product with the moistened cloth.
- 4) Dry the product with a dry microfibre cloth.
- 5) Allow the product to air dry.

Cleaning after contact with chlorine water

- > **Required materials and tools:** microfibre cloth, tap water
- 1) Rinse the joint thoroughly under running tap water.
- 2) Allow the joint to drain and dry it thoroughly with a microfibre cloth.
- 3) Check the functionality of the joint (flexing, extending).

Cleaning after contact with salt water

NOTICE! Remove salt crystals from the prosthetic knee joint to prevent corrosion.

- > **Required materials and tools:** 2 microfibre cloths, warm tap water
 - 1) Fully immerse the joint in warm tap water and flex and extend it several times.
 - If salt crystals stick, carefully wipe them off with a microfibre cloth. Immerse the joint in water again and rinse it off.
 - 2) Remove the joint from the water and allow it to drain.
 - 3) Dry it thoroughly with a microfibre cloth.
 - 4) Check the functionality of the joint (flexing, extending).
- After maximum use in salt water (**1 hour/day for a maximum of 14 days/year**), have it checked by qualified personnel.

Disinfecting the product

- > **Required materials and tools:** 2 microfibre cloths, water-based disinfectant
- 1) Apply disinfectant to the microfibre cloth.
- 2) Disinfect the product with the moistened cloth.
- 3) Dry the product with a dry microfibre cloth.
- 4) Allow the product to air dry.

9 Maintenance

The section is intended for qualified personnel.

- ▶ **Every 12 months:** Have maintenance performed by qualified personnel. Reduce the maintenance interval based on the use of the prosthesis.
 - In case of insufficient functionality, adjust the prosthesis.
- ▶ Check the prosthesis for **wear and tear, functionality and unusual noises**.
- ▶ Do not lubricate or grease the prosthetic joint.
- Repair work must be performed exclusively by manufacturer service.
- ▶ Using the prosthetic knee joint in salt water constitutes an increased load. Despite thorough cleaning and care, signs of corrosion may occur because invisible salt residues are present. As long as the function of the product is not impaired, it does not have to be sent in for service.
- ▶ Check all product functions for safety. Check flexion and extension and the ergonomically balanced stride (EBS) in particular.
- ▶ In case of functional limitations and damage, send the product to the manufacturer's service department.

9.1 Replacing the stop bumper

Replace the stop bumpers in case of visible, pronounced wear and tear or loss of comfort.

- > The following instructions refer to the illustration: see fig. 17
- > **Required materials and tools:** flat-head screwdriver, 2x 4Z59=6.3X6 replacement buffer
- 1) Move the joint to maximum flexion ❶.
- 2) Remove the cover caps on both sides of the product:
 - **NOTICE! When prying off the cover caps, take care not to damage the protective stickers beneath them with the screwdriver.**
 - Carefully insert the flat-head screwdriver into both openings in the cover cap, one after the other ❷.
 - Carefully pry off the cover cap ❸, only levering the screwdriver towards the knee joint.
 - Repeat these steps on the opposite cover cap.
- 3) Carefully pry out the stop bumpers using the flat-head screwdriver ❹.
- 4) Insert the new stop bumpers ❺.
- 5) Move the joint to full extension ❻.
- 6) Install the cover caps on both sides of the joint ❼.
- At full extension, the joint has a damped stop and the stop bumpers are positioned correctly.

9.2 Checking the sacrificial anode

The sacrificial anode ❶ protects the product against electrochemical corrosion by corroding itself first, thereby protecting the product from damage (see fig. 18).

► Check the sacrificial anode **every 12 months**:

- 1) Remove dirt and crystals from the area around the sacrificial anode.
- 2) Check the wear and tear of the sacrificial anode using callipers:
 - Result **< 5 mm** ❷: The product can continue to be used. Replacement by the manufacturer's service department is not required.
 - Result **≥ 5 mm** ❸: The sacrificial anode is too worn. Have a new sacrificial anode installed by the manufacturer's service department.

10 Disposal

Do not dispose of the product in household waste. Improper disposal is harmful to health and the environment. Observe the instructions of the local authority regarding return, collection and disposal.

11 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

11.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

11.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices.

The EU declaration of conformity is available at the following internet address:
<http://www.ottobock.com>

11.3 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

12 Technical data

Proximal connection	Pyramid
---------------------	---------

Distal connection [ø mm]	30
Knee flexion angle [°]	165
EBS [°]	10
Knee flexion angle with flexion stops [°]	135
Total weight [g]	1065
Overall height [mm]	232
Build height [mm]	217
Distal build height [mm]	202
Proximal build height [mm]	15
Mobility grade	2, 3
Permissible body weight, maximum [kg]	125

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