

## 4R72 Double adapter



Article number	4R72=32	4R72=45	4R72=60	4R72=75
<b>Material</b>	Titanium			
<b>System height</b>	69 mm	82 mm	97 mm	112 mm
<b>Weight</b>	85 g	95 g	110 g	125 g
<b>Max. body weight</b>	150 kg			



647H34

- For the use in water and in humid environments the included set screws have to be changed to anti corrosive Titanium set screws 506G5=\* (e.g. 506G5=M8X12). The Titanium set screws are available in the dimensions M8X12, M8X14 and M8X16.

## 4R104 Double adapter, movable

The adapter connects 2 prosthetic components with pyramid adapter and also allows a translation adjustment medially or laterally in the frontal plane, or dorsally or ventrally in the sagittal plane.



Article number	4R104=60	4R104=75
<b>Material</b>	Titanium	
<b>System height</b>	97 mm	112 mm
<b>Weight</b>	215 g	225 g
<b>Displacement</b>	+/- 11 mm	
<b>Max. body weight</b>	100 kg	



647H137



## Double adapter



Article number	4R76	4R78
<b>Material</b>	Stainless steel	
<b>System height</b>	-32 mm	-30 mm
<b>Weight</b>	95 g	115 g
<b>Max. body weight</b>	150 kg	



647G300

About this catalogue

Modular lower limb prostheses

Lower limb prostheses for children

Initial/interim prostheses

Waterproof walking devices

Sport prostheses

Prosthetic feet

Adapter

Knee joints

Hip joints

Socket Technologies

Cosmetic covers

Exoskeletal

Index