

How to Size the Push care Wrist Brace

Measure around the wrist as shown.

NOTE: This measurement is only a guideline for determining size. If your measurement is near the size break point, you should try on both sizes to determine the best fit.

Select: RIGHT or LEFT



SIZE	INCHES	CENTIMETERS
1	5 ¹ / ₈ - 6	13 - 15
2	6 - 6 ⁵ / ₈	15 - 17
3	6 ⁵ / ₈ - 7 ¹ / ₂	17 - 19
4	7 ¹ / ₂ - 8 ¹ / ₄	19 - 21