



STERLING™ GEL LINER

The College Park Sterling is a seamless TPE gel liner designed to conform to the patient's limb to provide optimal comfort. The Thermoplastic Elastomer (TPE) gel can accommodate limb irregularities. This durable liner can be easy to don and doff with an optimal TPE surface. It features a soft, extra thick distal pad which may improve user comfort, protect sensitive distal ends, and conform to irregular distal shapes.



STERLING LINER PART NUMBER

	THICKNESS*	TYPE**	SIZE*	QUANTITY
SLK G	3,6 OR 9	L OR C	U	S-XL
				1 OR 2

*See charts above. Pin (L)ock or (C)ushion.

STERLING™ GEL SLEEVE

The new Sterling sleeve features Thermoplastic Elastomer (TPE) Gel, designed to use body heat and mold to a user's limb for a near-custom fit. It features a seamless textile to help prevent snagging, with a thicker middle section for durability.



PART NUMBER	DESCRIPTION
SGS 3541	Sterling Suspension Sleeve, Size 35-41 cm
SGS 4149	Sterling Suspension Sleeve, Size 41-49 cm
SGS 4955	Sterling Suspension Sleeve, Size 49-55 cm
SGS 5562	Sterling Suspension Sleeve, Size 55-62 cm
SGS 6276	Sterling Suspension Sleeve, Size 62-76 cm

TO DETERMINE THE CORRECT SIZE:

Measure the circumference of the residual limb 15 cm above mid-patella and choose the corresponding device size. See the sizing chart below.

SIZE	MEASUREMENT (CM)
35	35-41 cm
41	41-49 cm
49	49-55 cm
55	55-62 cm
62	62-76 cm