

Newport[®] ULTRA[™]

The new Newport ULTRA hip orthosis from Orthomerica is designed to provide protection and support for patients at risk for or with hip instability.

INDICATIONS

- Post-operative or injury repair hip patients who will wear the orthosis for protection during healing
- After reposition of hip luxation
- Muscular imbalance
- Patients who have been, or are likely to be non-compliant with rigid, bulky hip orthosis
- Patients whose cognitive function impairs their ability to follow hip precautions
- Arthroscopic hip procedures
- Labral repair, with or without gluteus medius repair
- Post-operative proximal hamstring repair
- Patients at risk to dislocate

CONTRAINDICATIONS

- Global instability
- Patients at risk for an anterior dislocation
- Chronic dislocation

FEATURES/BENEFITS

- No tools required to adjust brace
- Locking screw provided to prevent patient from adjusting flexion/extension range
- Protects the hip joint capsule by providing a limited and safe range of motion following surgical procedures
- Well-tolerated by patients
- Virtual joint is easily set to control hip flexion, extension, adduction and abduction
- Proprioceptive reminder
- Universal for right or left side
- Easily customized to provide an exact fit
- One size fits most
- Laterally placed abduction paddle helps to maintain the head of the femur in the acetabulum
- Compression belt design with pull tabs anchors the hip joint and thigh cuff to promote optimal alignment
- Accommodates hip development through two abduction/adduction joint mechanisms



Newport[®] ULTRA[™] hip orthosis

one size fits most waists 25"-53" (64cm-135cm)

Item	Model #
Complete Orthosis	3663

ORTHOMERICA[®]

(800) 446-6770 | www.orthomerica.com

6333 North Orange Blossom Trail • Orlando FL 32810 • USA



Emergo Europe
Prinsessegracht 20
2514 AP The Hague
The Netherlands