



DOUBLE CONNECTOR



• Double Pyramid Receiver

Avoid using a tube in tight spaces.

Nr	Max. Patient's Weight	Material	Weight	TH	BH
1K207-032	100 kg	Aluminium	85 g	32 mm	32 mm
1K207-045	100 kg	Aluminium	100 g	45 mm	45 mm
1K207-060	100 kg	Aluminium	115 g	60 mm	60 mm
1K207-075	100 kg	Aluminium	130 g	75 mm	75 mm



• Double Pyramid Receiver - 125 kg

Avoid using a tube in tight spaces.

Nr	Max. Patient's Weight	Material	Weight	TH	BH
1K207-P6-032	125 kg	Aluminium	100 g	32 mm	32 mm
1K207-P6-045	125 kg	Aluminium	130 g	45 mm	45 mm
1K207-P6-060	125 kg	Aluminium	145 g	60 mm	60 mm
1K207-P6-075	125 kg	Aluminium	160 g	75 mm	75 mm

• Telescopic Double Pyramid Connector

Double connector including two pyramid receivers. Can be adjusted in height by shortening, it allows the CPO to accurately adjust the prosthesis height. It can also provide rotation adjustment.



Nr	Max. Patient's Weight	Material	Weight	TH	BH
1K150	100 kg	Aluminium	109 to 170 g	40 to 100 mm	40 to 100 mm


Double connector including a receiver and a pyramid. It allows for individual height adjustment of the prosthesis as its height can be shortened.



Nr	Max. Patient's Weight	Material	Weight	TH	BH
1K151	100 kg	Aluminium	141 to 201 g	44 to 89 mm	29 to 89 mm

5C413 cutting tool : recommended to cut the 1K150 connector to the desired height.

• Double Male Pyramid Connector



Nr	Max. Patient's Weight	Material	Weight	TH	BH
1K190	125 kg	Stainless Steel	120 g	36 mm	7 mm