

Push[®] med Patella Brace

Comfortable, precise patella tendon pressure

- Distributes comfortable pressure to the patellar tendon via skin compatible, anatomically shaped, viscoelastic polymer pad
- Prevents downward slippage with the non-slip pad
- Ensures even distribution of pressure to the patellar tendon with the semi-rigid strap reinforcement
- Can be machine washed on gentle cycle at low temperature
- Creates soft, comfortable feel and wicks moisture through the Sympress[™] microfibers



Indications

- Patellofemoral pain
- Tendinitis of patella tendon (Jumper's knee)
- Osgood-Schlatter disease
- Subluxation/proximal hypermobility of tibiofibular joint

Measuring / Sizing

Bend your knee 90° and measure around the narrowest part of your knee.

One Size: Fits measurements of 11-17 ³/₄ inches (28-45 cm)

Easy Application! Pull brace up the leg until the polymer pad is directly over the patellar tendon (just below the patella). Adjust until snug but comfortable.



Insurance Information

- The innovative characteristics for some of the Push braces do not fit current L-code identifier categories; therefore, there is no L-code identifier for the Push med Patella Brace.
- BraceLab is not a provider. We do not file insurance nor can we provide information on insurance coverage or reimbursement rates. Please contact your insurance for filing information.