

Tamarack[®] Prosthetic Sleeve Protector

with **GlideWear[™]**

Targeted protection from skin irritation and sealing sleeve damage
in below-knee prosthetic vacuum systems.

Thoughtfully designed in collaboration with Jack Richmond, CPOA.

Part No.	Size	Knee Circumference
GW-PSP-S	Small (S/P)	9 - 13 in 23 - 33 cm
GW-PSP-M	Medium (M)	12 - 16 in 30 - 40 cm
GW-PSP-L	Large (L/G)	14 - 18 in 36 - 45 cm
GW-PSP-XL	XLarge (XL/XG)	17 - 22 in 43 - 56 cm

**Sleeve Protector
Installation Video**



CARE INSTRUCTIONS

Wash with mild detergent.



MATERIALS

85% Nylon, 15% Spandex,
Polyurethane Coating
DOES NOT CONTAIN LATEX



Tamarack Habilitation Technologies, Inc.

1670-94th Lane NE, Blaine, MN 55449 | Tel: 866.795.0057 | E-mail: info@tamarackhti.com

© 2026 Tamarack Habilitation Technologies, Inc. | ISO 13485:2016 Registered QMS | REV 08/2026

GlideWear is a registered trademark of MIPS, manufactured in the USA by Tamarack HTI

Patented US8646459 - US9241853 - GB2478206 - OA2738162 | Other Patents Pending | All Rights Reserved

