


phlebology measuring charts

MEASURING INSTRUCTIONS

1. Select product/style from measuring chart. 2. Measure at indicated locations on diagram. 3. Locate proper product size based on measurements.

mediven® comfort 15–20, 20–30, 30–40 mmHg

Circumference	I	II	III	IV	V	VI	VII
2 EXTRA-WIDE CALF (excludes sheer & soft)	N/A	N/A	16½"–19½" 42–49.5 cm	17¾"–20¾" 45–53 cm	19"–22" 48–56 cm	20¼"–23¼" 51.5–59 cm	21½"–24½" 54.5–62 cm
2 CALF	11"–14" 28–35.5 cm	11¾"–15¼" 30–39 cm	12¾"–16½" 32.5–42 cm	13¾"–17¾" 35–45 cm	14½"–19" 37–48.5 cm	15½"–20¼" 39.5–51.5 cm	16¼"–21½" 41.5–54.5 cm
1 ANKLE	7"–7¾" 18–20 cm	8"–8¾" 20.5–22.5 cm	9"–9¾" 23–25 cm	10"–10¾" 25.5–27.5 cm	11"–11¾" 28–30 cm	12"–13" 30.5–33 cm	13¼"–14¼" 33.5–36 cm
LEG LENGTH	B CALF ▶ STANDARD: 15"–19" (38–47.5 cm) PETITE: 12"–15" (31–38 cm)						



For All Calf Styles:

Measure **1** **2** and **B**