


phlebology measuring charts

MEASURING INSTRUCTIONS

1. Select product/style from measuring chart.
2. Measure at indicated locations on diagram.
3. Locate proper product size based on measurements.

mediven® comfort 15–20, 20–30, 30–40 mmHg

Circumference	I	II	III	IV	V	VI	VII
3 PANTY & THIGH	16½"–22½" 42–57 cm	17¾"–24¼" 45–61.5 cm	19¼"–26¼" 49–66.5 cm	20½"–28¾" 52–73 cm	22"–31" 56–79 cm	23¾"–33¼" 60–84.5 cm	25½"–34½" 64–87.5 cm
2 CALF	11"–14" 28–35.5 cm	11¾"–15¼" 30–39 cm	12¾"–16½" 32.5–42 cm	13¾"–17¾" 35–45 cm	14½"–19" 37–48.5 cm	15½"–20¼" 39.5–51.5 cm	16¼"–21½" 41.5–54.5 cm
1 ANKLE	7"–7¾" 18–20 cm	8"–8¾" 20.5–22.5 cm	9"–9¾" 23–25 cm	10"–10¾" 25.5–27.5 cm	11"–11¾" 28–30 cm	12"–13" 30.5–33 cm	13¼"–14¼" 33.5–36 cm
LEG LENGTH	C PANTY & THIGH ▶ STANDARD: 28"–34" (71–86.5 cm) PETITE: 23"–28" (58.5–71 cm)						



For All Panty & Thigh Styles:

Measure **1** **2** **3** and **C**