


Sizing Instructions: measure across the proximal phalanges of digits 2-5 just proximal to the PIP joints. Please make sure the fingers are in the adducted position. Do not wrap the tape down the lateral or medial side of the hand.

Consider the length of the fingers when deciding between 2 different sizes.

(i.e.. shorter fingers = smaller size).

NOTE: Product must be fitted by a therapist or orthotist



Item No.	Side	Width Across Digits 2-5	Size
SaeboStretch			
0815-55-093	Left	2 1/4" to 2 3/4"	S
0815-55-101	Left	2 3/4" to 3 1/4"	M
0815-55-119	Left	3 1/4" to 4"	L
0815-55-127	Right	2 1/4" to 2 3/4"	S
0815-55-135	Right	2 3/4" to 3 1/4"	M
0815-55-143	Right	3 1/4" to 4"	L