

Alignment Unit

Transtibial and Transfemoral Endoskeletal Alignment Unit

FILLAUER

- Coupling ring with disconnect and full transverse rotation of 360 degrees
- Adjustments automatically compensate for wear
- Euro-pattern compatibility
- Metric threads

125321 Transtibial Alignment Unit Complete Assembly

Includes:

- A. **124105** Universal Coupling Ring, with Top Plate
- B. **124140** A-P and M-L Slide Unit, with Rotation and Pyramid
- C. **124170** Pylon Tube Clamp, with Pyramid Receiver
- D. **125032** Pylon Tube, 30mm Diameter x 10" Long
- E. **125101** Lower Pylon Tube Clamp, Internal Locking (3/8"-16x4" Bolt, Part # K20104 and Flat Washer Included)

