

Specify Side of Body – Left or Right

| 39SRH – W/ FINGER SEPARATORS                                                                                                           |      |               |               |                |         |
|----------------------------------------------------------------------------------------------------------------------------------------|------|---------------|---------------|----------------|---------|
| Model #                                                                                                                                | SIZE | WIDTH OF PALM |               | PRODUCT LENGTH |         |
|                                                                                                                                        |      | IN            | CM            | IN             | CM      |
| 39SRH-S                                                                                                                                | S    | 2.5 – 3"      | 6.4 – 7.6 cm  | 11"            | 27.9 cm |
| 39SRH-M                                                                                                                                | M    | 3 – 3.5"      | 7.6 – 8.9 cm  | 11.5"          | 29.2 cm |
| 39SRH-L                                                                                                                                | L    | 3.5 – 3.75"   | 8.9 – 9.5 cm  | 13"            | 33 cm   |
| 39SRH-XL                                                                                                                               | XL   | 3.75 – 4.25"  | 9.5 – 10.8 cm | 14.5"          | 36.8 cm |
| ASSIGNED PDAC CODE L3807 OR OTS L3809   MSPC 2.0  |      |               |               |                |         |

Measurement Instructions

Measure width of patients palm across knuckles. (Ends of hand)