

NEW SEMI-FLEXIBLE LIGAMENT BRACE

SecuTec® Genu Flex

INNOVATIVE PROTECTION AND STABILIZATION OF THE KNEE




FLEX-LINK TECHNOLOGY

OPTIMUM FIT THROUGHOUT THERAPY AND REHABILITATION.

NATURALLY ADJUSTS TO
CHANGING LEG SHAPES.

GIVES THE MUSCLES
MAXIMUM FREEDOM
OF MOVEMENT.



► Flexible in the direction of movement and
► tension-stable in the direction of therapy

Flex-Links perfectly adapt to movement.
The SecuTec Genu flex stays securely in position
during activity.



SecuTec® Genu Flex

A PARADIGM SHIFT IN LIGAMENT BRACING.

INDICATIONS

- Rupture of the anterior and / or posterior cruciate ligament ACL / PCL)
- After ligament surgery / ligament reconstruction
- Severe and / or complex instability (traumatic, degenerative, e.g. "unhappy triad")
- Collateral ligament injuries
- To restrict the knee joint's range of motion
- After meniscus surgery
- For non-surgical treatment or long-term rehabilitation
- Functional instabilities

TOOL-FREE RANGE OF MOTION CONTROL:

- Flexion: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°
- Extension: 0°, 10°, 20°, 30°, 45°
- The joint can be immobilized in steps of 0°, 10°, 20°, 30°, 45°.

L1845/L1852
PDAC verified



NEW

Telescoping
Hinge

Adjustable hinge
for customized
ROM limitation
and immobilization

INNOVATIVE SEMI-FLEXIBLE FRAME

Adapts to individual leg
circumference
through Flex-Links



SUCCESSFUL THERAPY

STABILIZATION. PROTECTION.

WEARING COMFORT.



STABILIZATION AND RELIEF

- Stabilization according to 4-point principle by combining a semi-flexible frame construction and non-elastic straps
- Innovative Flex-Modules easily adjust to thigh and calf circumference
- Adjustable over the course of treatment



ALL DAY WEARING COMFORT

- Lightweight, low profile design can be worn under most pants
- Semi-flexible frame adapts to the leg anatomy
- Soft microfiber padding
- Easy donning and doffing



- Proven therapy success
- Limited inventory required due to adjustable Flex-Link-Design that easily adjusts to leg circumference
- Cost effective reimbursement

BAUERFEIND KNEE SOLUTIONS FOR ALL WALKS OF LIFE.

Bauerfeind is known globally for world class quality, innovation and manufacturing of orthopedic braces that provide proven therapeutic benefits and are extremely comfortable to wear. With the new SecuTec GenuFlex, Bauerfeind now offers a complete range of sports-medicine solutions for the management of ligament injuries.



GenuTrain

**GenuTrain S/S Pro
L1832/L1833**

**SecuTec Genu Flex
L1845/L1852**

**SecuTec Genu
L1845/L1852**

**SofTecGenu
L1832/L1833**

Post-injury
stabilization



Post-operative
immobilization



Protection during
rehabilitation



Long-term
functional bracing




Long-term joint control
(Improved
proprioception)



SIZE AND ORDER INFORMATION

SecuTec® Genu Flex



Size	1	2	3	4	5	Order number
Circumfer. 6 3/4" above knee	14 1/2 - 17 3/4	17 3/4 - 20 3/4	20 3/4 - 24	24 - 27 1/4	27 1/4 - 32	12146781050+0°
Circumfer. 17 cm above knee	37 - 45	45 - 53	53 - 61	61 - 69	69 - 81	
Circumfer. 6" below knee	11 - 14 1/4	12 1/4 - 15 3/4	14 1/4 - 17 1/4	15 3/4 - 19	17 1/4 - 20 1/2	
Circumfer. 15 cm below knee	28 - 36	32 - 40	36 - 44	40 - 48	44 - 52	

+6 = right, +7 = left, ° = size

BAUERFEIND USA, INC.

75 14th St NE

Suite 2350

Atlanta, GA 30309

P 1 800-423-3405

P 404-201-7800

F 404-201-7839

info@bauerfeindusa.com



Bauerfeind supports are used and endorsed by the National Basketball Athletic Trainers Association.