

Specifications

| Article number | Side | Size | Operating voltage | Opening width | Proportional gripping force | Proportional speed | Weight with inner hand | For |
|----------------|-------|-------|-------------------|---------------|-----------------------------|--------------------|------------------------|--------------------|
| 8E44=6-L7 1/4 | Left | 7 1/4 | 6/7.2 V | 100 mm | 0 - 90 N | 15 - 130 mm/sec | ~308 g | Women, Adolescents |
| 8E44=6-L7 3/4 | Left | 7 3/4 | 6/7.2 V | 100 mm | 0 - 90 N | 15 - 130 mm/sec | ~308 g | Men |
| 8E44=6-L8 1/4 | Left | 8 1/4 | 6/7.2 V | 100 mm | 0 - 90 N | 15 - 130 mm/sec | ~308 g | Men |
| 8E44=6-R7 1/4 | Right | 7 1/4 | 6/7.2 V | 100 mm | 0 - 90 N | 15 - 130 mm/sec | ~308 g | Women, Adolescents |
| 8E44=6-R7 3/4 | Right | 7 3/4 | 6/7.2 V | 100 mm | 0 - 90 N | 15 - 130 mm/sec | ~308 g | Men |
| 8E44=6-R8 1/4 | Right | 8 1/4 | 6/7.2 V | 100 mm | 0 - 90 N | 15 - 130 mm/sec | ~308 g | Men |

Myo Glove Sizing for Hands

| Hand Size | For | Standard Glove Article Number | Skin Natural Glove Article Number |
|-----------|-----------------------|-------------------------------|-----------------------------------|
| 5 | Children 1-3 yrs old | 8S20=136x41 (L/R) color | 8S20N=136x41 (L/R) color |
| 5 1/2 | Children 3-6 yrs old | 8S20=147x45 (L/R) color | 8S20N=147x45 (L/R) color |
| 6 | Children 5-10 yrs old | 8S20=162x56 (L/R) color | 8S20N=162x56 (L/R) color |
| 6 1/2 | Children 8-13 yrs old | 8S20=177x64 (L/R) color | 8S20N=177x64 (L/R) color |
| 7 | Women, Adolescents | 8S13=7 (L/R) color | 8S13=7 (L/R) color |
| 7 1/4 | Men, Adolescents | 8S11=190x76 (L/R) color | 8S11N=190x76 (L/R) color |
| 7 1/4 | Women | 8S12=190x78 (L/R) color | 8S12N=190x78 (L/R) color |
| 7 3/4 | Men | 8S11=210x78 (L/R) color | 8S11N=210x78 (L/R) color |
| 8 1/4 | Men | 8S11=225x80 (L/R) color | 8S11N=225x80 (L/R) color |